

EXHIBIT G-2
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY
STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000
STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST Daniel McGinnis AGENCY VL50
DATE 04-08-09 TIME 0003
INTOXILYZER SERIAL # 80-001019 LOCATION Prescott Jail

- DM* 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution .10 AC.
OR
Pour a ~~standard alcohol concentration solution~~ _____ AC, into a clean dry simulator and assemble the simulator. Ensure that a ~~tight seal has been made. Turn on the simulator and allow temperature to reach 34°C ± 0.2°C~~
- DM* 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
DM 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
DM 4. Air blank completed.
DM 5. Calibration check completed. Test results 0. 098 AC.
DM 6. Air blank completed.
DM 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE *Daniel McGinnis* 543/QAS #38426

INTOXILYZER 8000

Location: YCSO PRESCOTT JAIL

Serial Number: 80-001019

04/08/2009

00:03:55

Standard Lot#: 818201G

Last Changed By: D. MCGINNIS #8086

QAS: D. MCGINNIS #8086

YAVAPAI CSO

Test	g/210L	Time
Air Blank	0.000	00:04:23
Cal Check	0.098	00:04:41
Air Blank	0.000	00:05:11

 S43/QAS #33426
Operator's Signature

EXHIBIT G-4
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY
STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000

STANDARD QUALITY ASSURANCE PROCEDURE

QA SPECIALIST Daniel McGinnis AGENCY YCSO

DATE 04-08-09 TIME 0006

INTOXILYZER SERIAL # 80-001019 LOCATION Prescott Jail

1. Display Reads "PUSH BUTTON TO START"

DIAGNOSTIC TESTS

1. Clock time check.
2. Date check.

OPERATIONAL TESTS

1. Alcohol-free subject test result 0. 000 AC.
2. Error recognition logic system functioning.
Not a Successfully Completed Test Sequence printed
3. Proper sample recognition system.
Not a Successfully Completed Test Sequence printed
Deficient sample printed.
4. Standard Calibration Check standard 0. 10 AC. Result 0. 017 AC.

Instrument is operating properly and accurately. YES DM NO

COMMENTS

All OK

SIGNATURE Daniel McGinnis 243/QAS #33426

INTOXILYZER 8000

Location: YCSO PRESCOTT JAIL

Serial Number : 80-001019

04/08/2009

00:06:02

Standard Lot#: 818201G

Last Changed By: D. MCGINNIS #8086

QAS: D. MCGINNIS #8086

YAVAPAI CSO

Operator: D. MCGINNIS #8086

YAVAPAI CSO

Subject: SQAP,TEST

DOB: 09/09/1990

Sex: M

Weight: 190

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	0.000	00:06:56
Diagnostic Test	Pass	00:07:24
Air Blank	0.000	00:07:51
0.100 Cal Check	0.098	00:08:10
Air Blank	0.000	00:08:39
Subject Test	0.000	00:09:09
Air Blank	ABT*	00:09:20
Air Blank	RFI**	00:09:27

*Sequence Aborted

**RFI Detect

Not a Successfully Completed Test Sequence

INTOXILYZER 8000

Location: YCSO PRESCOTT JAIL

Serial Number : 80-001019

04/08/2009

00:10:17

Standard Lot#: 818201G

Last Changed By: D. MCGINNIS #8086

QAS: D. MCGINNIS #8086

YAVAPAI CSO

Operator: D. MCGINNIS #8086

YAVAPAI CSO

Subject: SQAP,TEST

DOB: 09/09/1990

Sex: M

Weight: 190

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	0.000	00:11:18
Diagnostic Test	Pass	00:11:46
Air Blank	0.000	00:12:13
0.100 Cal Check	0.098	00:12:32
Air Blank	0.000	00:13:02
Subject Test	0.042*	00:16:10
Air Blank	0.000	00:16:42
Air Blank	0.000	00:17:10
Subject Test	XXX**	00:17:49
Air Blank	IPS***	00:17:57

*Deficient Sample

**Invalid Sample

Begin new deprivation period

***Improper Sample

Not a Successfully Completed Test Sequence

EXHIBIT G-2
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY
STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000
STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST Daniel McGinnis AGENCY YCSO
DATE 04-08-09 TIME 0019
INTOXILYZER SERIAL # 80-00109 LOCATION Prescott Jail

- DM* 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution 10 AC.
OR
Pour a standard alcohol concentration solution _____ AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$
- DM* 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
DM 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
DM 4. Air blank completed.
DM 5. Calibration check completed. Test results 0. 017 AC.
DM 6. Air blank completed.
DM 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE

Daniel McGinnis 543/QAS#33426

INTOXILYZER 8000

Location: YCSO PRESCOTT JAIL

Serial Number: 80-001019

04/08/2009

00:19:01

Standard Lot#: 818201G


Last Changed By: D. MCGINNIS #8086

QAS: D. MCGINNIS #8086

YAVAPAI CSO

Test	g/210L	Time
Air Blank	AMB*	00:19:29
Air Blank	0.000	00:19:58

*Ambient Fail


Operator's Signature 543/QAS #33426

INTOXILYZER 8000

Location: YCSO PRESCOTT JAIL

Serial Number: 80-001019

04/08/2009

00:20:08

Standard Lot#: 818201G

Last Changed By: D. MCGINNIS #8086

QAS: D. MCGINNIS #8086

YAVAPAI CSO

Test	g/210L	Time
Air Blank	0.000	00:20:35
Cal Check	0.097	00:20:54
Air Blank	0.000	00:21:24


Operator's Signature