

EXHIBIT G-2
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY
STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000
STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST D. McGinnis AGENCY YCSO
DATE 06-03-09 TIME 2128
INTOXILYZER SERIAL # 80-001019 LOCATION Prescott

- [Signature]* 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution .10 AC.
OR
Pour a standard alcohol concentration solution AC into a clean dry simulator and assemble the simulator. ~~Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach 34°C ± 0.2°C~~
- [Signature]* 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
- [Signature]* 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
- [Signature]* 4. Air blank completed.
- [Signature]* 5. Calibration check completed. Test results 0. 097 AC.
- [Signature]* 6. Air blank completed.
- [Signature]* 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE *[Signature]* S43/OAS# 33426

INTOXILYZER 8000

Location: YCSO PRESCOTT JAIL

Serial Number: 80-001019

06/03/2009

21:28:56

Standard Lot#: 818201G

Last Changed By: D. MCGINNIS #8086

QAS: D. MCGINNIS #8086

YAVAPAI CSO

Test	g/210L	Time
Air Blank	0.000	21:29:23
Cal Check	0.097	21:29:42
Air Blank	0.000	21:30:12


 843/QAS33426
Operator's Signature

EXHIBIT G-4
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY
STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000

STANDARD QUALITY ASSURANCE PROCEDURE

QA SPECIALIST Daniel McGinnis AGENCY YCSO
DATE 06-03-09 TIME 2131
INTOXILYZER SERIAL # 80-001019 LOCATION Prescott

Jim 1. Display Reads "PUSH BUTTON TO START"

DIAGNOSTIC TESTS
Jim 1. Clock time check.
Jim 2. Date check.

OPERATIONAL TESTS
Jim 1. Alcohol-free subject test result 0. 000 AC. ~~_____~~
Jim 2. Error recognition logic system functioning. ~~_____~~
Not a Successfully Completed Test Sequence printed
Jim 3. Proper sample recognition system. ~~_____~~
Not a Successfully Completed Test Sequence printed
Deficient sample printed.
Jim 4. Standard Calibration Check standard 0. 10 AC. Result 0. 097 AC.

Instrument is operating properly and accurately. YES Jim NO _____

COMMENTS All OK

SIGNATURE Daniel McGinnis S93/QAS # 33426

INTOXILYZER 8000

Location: YCSO PRESCOTT JAIL

Serial Number : 80-001019

06/03/2009

21:31:15

Standard Lot#: 818201G

Last Changed By: D. MCGINNIS

#8086

QAS: D. MCGINNIS

#8086

YAVAPAI CSO

Operator: D. MCGINNIS

#8086

YAVAPAI CSO

Subject: SOAP,TEST

DOB: 09/09/1990

Sex: M

Weight: 200

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	0.000	21:32:25
Diagnostic Test	Pass	21:32:52
Air Blank	0.000	21:33:20
0.100 Cal Check	0.097	21:33:39
Air Blank	0.000	21:34:08
Subject Test	0.000	21:34:31
Air Blank	ABT*	21:34:39
Air Blank	RFI**	21:34:48

*Sequence Aborted

**RFI Detect

Not a Successfully Completed Test Sequence

INTOXILYZER 8000

Location: YCSO PRESCOTT JAIL

Serial Number : 80-001019

06/03/2009

21:35:06

Standard Lot#: 818201G

Last Changed By: D. MCGINNIS #8086

QAS: D. MCGINNIS #8086

YAVAPAI CSO

Operator: D. MCGINNIS #8086

YAVAPAI CSO

Subject: SOAP,TEST

DOB: 09/09/1990

Sex: M Weight: 200

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	0.000	21:36:00
Diagnostic Test	Pass	21:36:27
Air Blank	0.000	21:36:54
0.100 Cal Check	0.097	21:37:13
Air Blank	0.000	21:37:43
Subject Test	0.086*	21:40:51
Air Blank	0.000	21:41:22
Air Blank	0.000	21:41:50
Subject Test	XXX**	21:42:38
Air Blank	IPS***	21:42:46

*Deficient Sample

**Invalid Sample

Begin new deprivation period

***Imple

Not a Successfully Completed Test Sequence

EXHIBIT G-2
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY
STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000

STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST Daniel McGinnis AGENCY YCSO
DATE 06-03-09 TIME 2148
INTOXILYZER SERIAL # 80-001019 LOCATION Prescott

1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution .10 AC.
OR

~~Pour a standard alcohol concentration solution AC into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$~~

2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
4. Air blank completed.
5. Calibration check completed. Test results 0. 097 AC.
6. Air blank completed.
7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE Daniel McGinnis 543/QAS # 33426

INTOXILYZER 8000

Location: YCSO PRESCOTT JAIL

Serial Number: 80-001019

06/03/2009

21:44:23

Standard Lot#: 818201G

Last Changed By: D. MCGINNIS #8086

QAS: D. MCGINNIS #8086

YAVAPAI CSO

Test	g/210L	Time
Air Blank	AMB*	21:44:51
Air Blank	0.000	21:45:19

*Ambient Fail


Operator's Signature

INTOXILYZER 8000

Location: YCSO PRESCOTT JAIL

Serial Number: 80-001019

06/03/2009

21:48:24

Standard Lot#: 818201G

Last Changed By: D. MCGINNIS #8086

QAS: D. MCGINNIS #8086

YAVAPAI CSO

Test	g/210L	Time
Air Blank	0.000	21:48:52
Cal Check	0.097	21:49:10
Air Blank	0.000	21:49:40

 543/QAS # 33426
Operator's Signature