

EXHIBIT G-2
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY
STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000
STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST Daniel McGinnis AGENCY YCSO

DATE 08-22-09 TIME 0217

INTOXILYZER SERIAL # 80-001019 LOCATION Prescott Jail

- [Handwritten initials]* 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution 10 AC.
OR
Pour a standard alcohol concentration solution _____ AC, into a clean dry simulator and assemble the simulator. ~~Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach 34°C ± 0.2°C~~
- [Handwritten initials]* 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
- [Handwritten initials]* 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
- [Handwritten initials]* 4. Air blank completed.
- [Handwritten initials]* 5. Calibration check completed. Test results 0. 097 AC.
- [Handwritten initials]* 6. Air blank completed.
- [Handwritten initials]* 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE *[Handwritten Signature]* 543/QAS 33426

INTOXILYZER 8000

Location: YCSO PRESCOTT JAIL

Serial Number: 80-001019

08/22/2009

02:17:55

Standard Lot#: 818201G

Last Changed By: D. MCGINNIS #8086

QAS: D. MCGINNIS #8086

YAVAPAI CSO

Test	g/210L	Time
Air Blank	0.000	02:18:22
Cal Check	0.097	02:18:41
Air Blank	0.000	02:19:10


Operator's Signature