

**EXHIBIT G-2**  
**THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)**  
**ARIZONA DEPARTMENT OF PUBLIC SAFETY**  
**STANDARD QUALITY ASSURANCE PROCEDURES**  
**INTOXILYZER MODEL 8000**  
**STANDARD CALIBRATION CHECK PROCEDURE**

QA SPECIALIST DAN RAISS AGENCY YAVAPAI Co. S.O.  
 DATE 2-6-08 TIME 1110  
 INTOXILYZER SERIAL # 80-001095 LOCATION YCSO-MAYER

- OR 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution 100 AC.  
 OR  
 Pour a standard alcohol concentration solution \_\_\_\_\_ AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach 34°C ± 0.2°C
- OR 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
- OR 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
- OR 4. Air blank completed.
- OR 5. Calibration check completed. Test results 0. 098 AC.
- OR 6. Air blank completed.
- OR 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE Daniel Q. Raiss #5-31/QAS #33026

INTOXILYZER 8000

Location: YCSO-MAYER

Serial Number: 80-001095

02/06/2008

11:15:24

-----  
Standard Lot#: 610201I

Last Changed By: D. RAISS #S-31  
-----

QAS: D. RAISS #3693

YAVAPAI CSO  
-----

Test	g/210L	Time
Air Blank	0.000	11:15:51
Cal Check	0.098	11:16:10
Air Blank	0.000	11:16:39

*Donald Q. Raiss #S-31/QAS #33026*  
Operator's Signature