

EXHIBIT G-2
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY
STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000
STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST Daniel McGinnis AGENCY YCSO

DATE 07-26-09 TIME 0826

INTOXILYZER SERIAL # 80-01095 LOCATION Meyer

1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution 100 AC.
OR

Pour a standard alcohol concentration solution _____ AC, into a clean dry simulator and assemble the simulator. ~~Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach 34°C ± 0.2°C~~

2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"

3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.

4. Air blank completed.

5. Calibration check completed. Test results 0. 098 AC.

6. Air blank completed.

7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE Daniel McGinnis S43/QAS # 33426

INTOXILYZER 8000

Location: YCSO-MAYER

Serial Number: 80-001095

07/26/2009

02:26:52

Standard Lot#: 733201G

Last Changed By: D. RAISS #3693

QAS: D. MCGINNIS #8086

YAVAPAI CSO

Test	g/210L	Time
Air Blank	0.000	02:27:20
Cal Check	0.098	02:27:39
Air Blank	0.000	02:28:08


Operator's Signature

EXHIBIT G-4
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY
STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000
STANDARD QUALITY ASSURANCE PROCEDURE

QA SPECIALIST Daniel McGinnis AGENCY YCSO

DATE 07-26-09 TIME 0229

INTOXILYZER SERIAL # 80-001095 LOCATION Mayer

DM 1. Display Reads "PUSH BUTTON TO START"

DIAGNOSTIC TESTS

- DM 1. Clock time check.
DM 2. Date check.

OPERATIONAL TESTS

- DM 1. Alcohol-free subject test result 0. 000 AC.
DM 2. Error recognition logic system functioning.
Not a Successfully Completed Test Sequence printed
DM 3. Proper sample recognition system.
Not a Successfully Completed Test Sequence printed
Deficient sample printed.
DM 4. Standard Calibration Check standard 0. 100 AC. Result 0. 098 AC.

Instrument is operating properly and accurately. YES DM NO _____

COMMENTS All OK

SIGNATURE Daniel McGinnis S43/QA# 33426

INTOXILYZER 8000
Location: YCSO-MAYER
Serial Number : 80-001095
07/26/2009 02:29:14

Standard Lot#: 733201G
Last Changed By: D. RAISS #3693

QAS: D. MCGINNIS #8086
YAVAPAI CSO

Operator: D. MCGINNIS #8086
YAVAPAI CSO

Subject: SQAP,TEST
DOB: 09/09/1990
Sex: M Weight: 200

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	0.000	02:30:07
Diagnostic Test	Pass	02:30:34
Air Blank	0.000	02:31:01
0.100 Cal Check	0.097	02:31:20
Air Blank	0.000	02:31:49
Subject Test	0.000	02:32:22
Air Blank	ABT*	02:32:30
Air Blank	RFI**	02:32:39

*Sequence Aborted
**RFI Detect

Not a Successfully Completed Test Sequence

INTOXILYZER 8000
Location: YCSO-MAYER
Serial Number : 80-001095
07/26/2009

02:32:51

Standard Lot#: 733201G
Last Changed By: D. RAISS #3693

QAS: D. MCGINNIS #8086
YAVAPAI CSO

Operator: D. MCGINNIS #8086
YAVAPAI CSO

Subject: SQAP,TEST
DOB: 09/09/1990
Sex: M Weight: 200

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	0.000	02:33:46
Diagnostic Test	Pass	02:34:14
Air Blank	0.000	02:34:41
0.100 Cal Check	0.098	02:35:00
Air Blank	0.000	02:35:29
Subject Test	0.028*	02:38:38
Air Blank	0.000	02:39:09
Air Blank	0.000	02:39:36
Subject Test	XXX**	02:40:03
Air Blank	IPS***	02:40:13

*Deficient Sample

**Invalid Sample

Begin new deprivation period

***Imple

Not a Successfully Completed Test Sequence

EXHIBIT G-2
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY
STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000

STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST Daniel McGinnis AGENCY YCSO

DATE 07-26-09 TIME 0242

INTOXILYZER SERIAL # 80-001095 LOCATION Mayer

DM 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution, 100 AC.
OR

~~Pour a standard alcohol concentration solution _____ AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$~~

- DM 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
- DM 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
- DM 4. Air blank completed.
- DM 5. Calibration check completed. Test results 0. 0.08 AC.
- DM 6. Air blank completed.
- DM 7. Remove printed record. Attach the record to the completed checklist.


SIGNATURE Daniel McGinnis S43/QAS # 33426

INTOXILYZER 8000
Location: YCSO-MAYER
Serial Number: 80-001095
07/26/2009 02:42:05

Standard Lot#: 733201G
Last Changed By: D. RAISS #3693

QAS: D. MCGINNIS #8086
YAVAPAI CSO

Test	g/210L	Time
Air Blank	0.000	02:42:32
Cal Check	0.098	02:42:51
Air Blank	0.000	02:43:21


Operator's Signature