

EXHIBIT G-2  
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY  
STANDARD QUALITY ASSURANCE PROCEDURES  
INTOXILYZER MODEL 8000  
STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST Daniel McGinnis AGENCY YCSO

DATE 10-17-09 TIME 2135

INTOXILYZER SERIAL # 80-001095 LOCATION Mayer

- DM* 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution .10 AC.  
OR  
Pour a standard alcohol concentration solution \_\_\_\_\_ AC, into a clean dry simulator and assemble the simulator. ~~Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach 34°C ± 0.2°C~~
- DM* 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
- DM* 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
- DM* 4. Air blank completed.
- DM* 5. Calibration check completed. Test results 0. 097 AC.
- DM* 6. Air blank completed.
- DM* 7. Remove printed record. Attach the record to the completed checklist.


SIGNATURE *Daniel McGinnis* 543/QAS # 33426

INTOXILYZER 8000  
Location: YCSO-MAYER  
Serial Number: 80-001095  
10/17/2009 21:35:57

-----  
Standard Lot#: 901401I  
Last Changed By: D. MCGINNIS #8086  
-----

QAS: D. MCGINNIS #8086  
YAVAPAI CSO  
-----

| Test      | g/210L | Time     |
|-----------|--------|----------|
| Air Blank | 0.000  | 21:36:24 |
| Cal Check | 0.097  | 21:36:43 |
| Air Blank | 0.000  | 21:37:13 |

  
Operator's Signature