

EXHIBIT G-2  
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY  
STANDARD QUALITY ASSURANCE PROCEDURES  
INTOXILYZER MODEL 8000  
STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST RICHARD J. BYRNE III AGENCY AZ DPS  
DATE 10/20/10 TIME 2108  
INTOXILYZER SERIAL # 80-000529 LOCATION WICKENBURG P.D.

- RS* 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution .100 AC.  
OR  
Pour a standard alcohol concentration solution N/A AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach  $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$
- RS* 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
- RS* 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
- RS* 4. Air blank completed.
- RS* 5. Calibration check completed. Test results 0.899 AC.
- RS* 6. Air blank completed.
- RS* 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE Richard J. Byrne III #6386

INTOXILYZER 8000

Location: WICKENBURG PD

Serial Number: 80-000929

Core Version: 8105.48

10/20/2010

21:08:43

-----  
Standard Lot#: 590854

Last Changed By: R. BYRNE, III #6386  
-----

QAS: R. BYRNE, III #6386

AZDPS  
-----

Operator: R. BYRNE, III #6386

AZDPS  
-----

Test	g/210L	Time
Air Blank	0.000	21:09:15
0.100 Cal Check	0.099	21:09:36
Air Blank	0.000	21:10:05

Preventative Maintenance  
Performed -- Timer Reset