

EXHIBIT G-2  
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY  
STANDARD QUALITY ASSURANCE PROCEDURES  
INTOXILYZER MODEL 8000

STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST Daniel McInnis AGENCY YCSO

DATE 11-14-10 TIME 0458

INTOXILYZER SERIAL # 80-001095 LOCATION Mayer

- 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution, 10 AC OR  
~~Pour a standard alcohol concentration solution \_\_\_\_\_ AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach 34°C ± 0.2°C~~
- 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
- 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
- 4. Air blank completed.
- 5. Calibration check completed. Test results 0. 097 AC.
- 6. Air blank completed.
- 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE Daniel McInnis 543/QAS # 33426

INTOXILYZER 8000  
Location: YCSO-MAYER  
Serial Number: 80-001095  
Core Version: 8105.48  
11/14/2010 04:58:23

-----  
Standard Lot#: 901401I  
Last Changed By: D. MCGINNIS #8086

-----  
QAS: D. MCGINNIS #8086  
YAVAPAI CSO

-----  
Operator: D. MCGINNIS #8086  
YAVAPAI CSO

-----  
Test g/210L Time  
-----  
Air Blank 0.000 04:58:54  
0.100 Cal Check 0.097 04:59:15  
Air Blank 0.000 04:59:44

Preventative Maintenance  
Performed -- Timer Reset

EXHIBIT G-4  
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY  
STANDARD QUALITY ASSURANCE PROCEDURES  
INTOXILYZER MODEL 8000  
STANDARD QUALITY ASSURANCE PROCEDURE

QA SPECIALIST Daniel McGinnis AGENCY YCSO  
DATE 11-14-10 TIME 0500  
INTOXILYZER SERIAL # 80-001095 LOCATION Mayer

DM 1. Display Reads "PUSH BUTTON TO START"

DIAGNOSTIC TESTS

- DM 1. Clock time check.
- DM 2. Date check.

OPERATIONAL TESTS

- DM 1. Alcohol-free subject test result 0. 000 AC.
- DM 2. Error recognition logic system functioning.  
Not a Successfully Completed Test Sequence printed
- DM 3. Proper sample recognition system.  
Not a Successfully Completed Test Sequence printed  
Deficient sample printed.
- DM 4. Standard Calibration Check standard 0. 10 AC. Result 0. 097 AC.

Instrument is operating properly and accurately. YES  NO

COMMENTS All OK

SIGNATURE Daniel McGinnis 543/QA#33427

INTOXILYZER 8000  
Location: YCSO-MAYER  
Serial Number: 80-001095  
Core Version: 8105.48  
11/14/2010 05:00:49

---

Standard Lot#: 901401I  
Last Changed By: D. MCGINNIS #8086

---

QAS: D. MCGINNIS #8086  
YAVAPAI CSO

---

Operator: D. MCGINNIS #8086  
YAVAPAI CSO

---

Subject: SQAP,TEST  
DOB: 09/09/1990  
Sex: M Weight: 200

---

15 Minute Deprivation Period? Yes

---

Test	g/210L	Time
Air Blank	0.000	05:01:58
Diagnostic Test	Pass	05:02:32
Air Blank	0.000	05:02:59
0.100 Cal Check	0.097	05:03:20
Air Blank	0.000	05:03:49
Subject Test	0.000	05:04:23
Air Blank	ABT*	05:04:31
Air Blank	RFI**	05:04:52

\*Sequence Aborted  
\*\*RFI Detect

**Not a Successfully Completed Test Sequence**

INTOXILYZER 8000  
Location: YCSO-MAYER  
Serial Number: 80-001095  
Core Version: 8105.48  
11/14/2010 05:05:53

---

Standard Lot#: 901401I  
Last Changed By: D. MCGINNIS #8086

---

QAS: D. MCGINNIS #8086  
YAVAPAI CSO

---

Operator: D. MCGINNIS #8086  
YAVAPAI CSO

---

Subject: SQAP,TEST  
DOB: 09/09/1990  
Sex: M Weight: 200

---

15 Minute Deprivation Period? Yes

---

Test	g/210L	Time
Air Blank	0.000	05:06:52
Diagnostic Test	Pass	05:07:26
Air Blank	0.000	05:07:53
0.100 Cal Check	0.097	05:08:14
Air Blank	0.000	05:08:43
Subject Test	0.128*	05:11:51
Air Blank	0.000	05:12:25
Air Blank	0.000	05:12:52
Subject Test	XXX**	05:13:43
Air Blank	0.000	05:14:25

\*Deficient Sample

\*\*Invalid Sample

Begin new deprivation period

**Not a Successfully Completed Test Sequence**

INTOXILYZER 8000  
Location: YCSO-MAYER  
Serial Number: 80-001095  
Core Version: 8105.48  
11/14/2010

05:15:59

---

Standard Lot#: 901401I  
Last Changed By: D. MCGINNIS #8086

---

QAS: D. MCGINNIS #8086  
YAVAPAI CSO

---

Operator: D. MCGINNIS #8086  
YAVAPAI CSO

---

Subject: SQAP,TEST  
DOB: 09/09/1990  
Sex: M Weight: 200

---

15 Minute Deprivation Period? Yes

---

Test	g/210L	Time
Air Blank	IPS*	05:16:49
Air Blank	0.000	05:17:17

\*Improper Sample

**Not a Successfully Completed Test Sequence**

EXHIBIT G-2  
FIELD REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (1)

ARIZONA DEPARTMENT OF PUBLIC SAFETY  
STANDARD QUALITY ASSURANCE PROCEDURES  
INTOXILYZER MODEL 8000  
STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST Daniel McGinnis AGENCY YCSO

DATE 11-14-10 TIME 0518

INTOXILYZER SERIAL # 80-001095 LOCATION Mayer

- DM* 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution .10 AC  
OR  
Pour a standard alcohol concentration solution ~~\_\_\_\_\_~~ AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach 34°C + 0.2°C
- DM* 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
- DM* 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
- DM* 4. Air blank completed.
- DM* 5. Calibration check completed. Test results 0. 097 AC.
- DM* 6. Air blank completed.
- DM* 7. Remove printed record? Attach the record to the completed checklist.

SIGNATURE *Daniel McGinnis* 543/QAS# 33426

INTOXILYZER 8000  
Location: YCSO-MAYER  
Serial Number: 80-001095  
Core Version: 8105.48  
11/14/2010 05:18:46

-----  
Standard Lot#: 9014011  
Last Changed By: D. MCGINNIS #8086  
-----

QAS: D. MCGINNIS #8086  
YAVAPAI CSO  
-----

Operator: D. MCGINNIS #8086  
YAVAPAI CSO  
-----

Test	g/210L	Time
Air Blank	0.000	05:19:16
0.100 Cal Check	0.097	05:19:37
Air Blank	0.000	05:20:06

Preventative Maintenance  
Performed -- Timer Reset