

EXHIBIT G-2  
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (1)

ARIZONA DEPARTMENT OF PUBLIC SAFETY  
STANDARD QUALITY ASSURANCE PROCEDURES  
INTOXILYZER MODEL 8000  
STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST Daniel McGinnis AGENCY YCSO

DATE 01-12-11 TIME 0906

INTOXILYZER SERIAL # 80-001075 LOCATION Meyer

1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution 10 AC.  
OR  
Pour a ~~standard alcohol concentration solution~~ AC, into a clean dry simulator and assemble the simulator. ~~Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach 34°C + 0.2°C~~
2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
4. Air blank completed.
5. Calibration check completed. Test results 0. 097 AC.
6. Air blank completed.
7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE Daniel McGinnis 543/WAS# 33926

INTOXILYZER 8000  
Location: YCSO-MAYER  
Serial Number: 80-001095  
Core Version: 8105.48  
01/12/2011 04:06:14

-----  
Standard Lot#: 9014011  
Last Changed By: D. MCGINNIS #8086

-----  
QAS: D. MCGINNIS #8086  
YAVAPAI CSO

-----  
Operator: D. MCGINNIS #8086  
YAVAPAI CSO

-----  
Test g/210L Time  
-----  
Air Blank 0.000 04:06:47  
0.100 Cal Check 0.097 04:07:08  
Air Blank 0.000 04:07:37

Preventative Maintenance  
Performed -- Timer Reset