

INTOXILYZER 8000

Location: YCSO-MAYER

Serial Number: 80-001095

Core Version: 8105.48

02/08/2011

20:21:50

Standard Lot#: 20910100A1

Last Changed By: D. MCGINNIS #8086

QAS: D. MCGINNIS #8086

YAVAPAI CSO

Operator: D. MCGINNIS #8086

YAVAPAI CSO

Test	g/210L	Time
Air Blank	0.000	20:23:14
0.100 Cal Check	0.098	20:23:35
Air Blank	0.000	20:24:04

Preventative Maintenance

Performed -- Timer Reset

EXHIBIT G-4
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY
STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000
STANDARD QUALITY ASSURANCE PROCEDURE

QA SPECIALIST Daniel McInnis AGENCY YCSO

DATE 02-08-11 TIME 2039

INTOXILYZER SERIAL # 80-001095 LOCATION Mayer

DM 1. Display Reads "PUSH BUTTON TO START"

DIAGNOSTIC TESTS

- DM 1. Clock time check.
- DM 2. Date check.

OPERATIONAL TESTS

- DM 1. Alcohol-free subject test result 0. 00 AC.
- DM 2. Error recognition logic system functioning.
Not a Successfully Completed Test Sequence printed
- DM 3. Proper sample recognition system.
Not a Successfully Completed Test Sequence printed
Deficient sample printed.
- DM 4. Standard Calibration Check standard 0. 10 AC. Result 0. 097 AC.

Instrument is operating properly and accurately. YES DM NO _____

COMMENTS All OK

SIGNATURE Daniel McInnis 543/QAS# 33426

INTOXILYZER 8000

Location: YCSO-MAYER

Serial Number: 80-001095

Core Version: 8105.48

02/08/2011

20:39:33

Standard Lot#: 20910100A1

Last Changed By: D. MCGINNIS #8086

QAS: D. MCGINNIS #8086

YAVAPAI CSO

Operator: D. MCGINNIS #8086

YAVAPAI CSO

Subject: SQAP,TEST

DOB: 09/09/1990

Sex: M

Weight: 200

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	0.000	20:40:27
Diagnostic Test	Pass	20:41:01
Air Blank	0.000	20:41:29
0.100 Cal Check	0.097	20:41:49
Air Blank	0.000	20:42:19
Subject Test	0.072*	20:45:26
Air Blank	0.000	20:45:59
Air Blank	0.000	20:46:26
Subject Test	XXX**	20:46:57
Air Blank	IPS***	20:47:05

*Deficient Sample

**Invalid Sample

Begin new deprivation period

***Improper Sample

Not a Successfully Completed Test Sequence

INTOXILYZER 8000
Location: YCSO-MAYER
Serial Number: 80-001095
Core Version: 8105.48
02/08/2011

20:33:04

Standard Lot#: 20910100A1
Last Changed By: D. MCGINNIS #8086

QAS: D. MCGINNIS #8086
YAVAPAI CSO

Operator: D. MCGINNIS #8086
YAVAPAI CSO

Subject: SQAP,TEST
DOB: 09/09/1990
Sex: M Weight: 200

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	0.000	20:34:04
Diagnostic Test	Pass	20:34:38
Air Blank	0.000	20:35:06
0.100 Cal Check	0.097	20:35:27
Air Blank	0.000	20:35:56
Subject Test	0.000	20:36:20
Air Blank	0.000	20:36:47
Wait	ABT*	
Air Blank	RFI**	20:37:25

*Sequence Aborted

**RFI Detect

Not a Successfully Completed Test Sequence

EXHIBIT G-2
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY
STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000
STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST Daniel McQuinnis AGENCY YCSO
DATE 02-08-11 TIME 2050
INTOXILYZER SERIAL # 80-001095 LOCATION Mayer

1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution, 10 AC.
OR
Pour a standard alcohol concentration solution _____ AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$
2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
4. Air blank completed.
5. Calibration check completed. Test results 0. 097 AC.
6. Air blank completed.
7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE

Daniel McQuinnis S4B/QAS # 83426

