

EXHIBIT G-2  
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY  
STANDARD QUALITY ASSURANCE PROCEDURES  
INTOXILYZER MODEL 8000

STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST LOGAN MOODY AGENCY CAVAPAI CO. S.O.  
DATE 1/10/13 TIME 1555  
INTOXILYZER SERIAL # 80-001095 LOCATION YUCO - MAGER

1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution .100 AC.  
OR  
~~Pour a standard alcohol concentration solution \_\_\_\_\_ AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach  $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$~~
2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"  
 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.  
 4. Air blank completed.  
 5. Calibration check completed. Test results 0. 093 AC.  
 6. Air blank completed.  
 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE LMoody 7486

INTOXILYZER 8000

Location: YCSO-MAYER

Serial Number: 80-001095

Core Version: 8105.48

01/10/2013

15:55:21

-----  
Standard Lot#: 01011100A1

Last Changed By: L. MOODY #7486  
-----

QAS: L. MOODY #7486

YAVAPAI CSO  
-----

Operator: L. MOODY #7486

YAVAPAI CSO  
-----

Test	g/210L	Time
Air Blank	0.000	15:57:26
0.100 Cal Check	0.093	15:57:46
Air Blank	0.000	15:58:16

Preventative Maintenance  
Performed -- Timer Reset