

EXHIBIT G-2
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY
STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000
STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST LOGAN MOODY #7486 AGENCY YAVAPAI CO. S.O.

DATE 4/6/2014 TIME 1452

INTOXILYZER SERIAL # 80-001095 LOCATION YCSO - MUGER

- (67) 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution .100 AC.
OR
~~Pour a standard alcohol concentration solution~~ AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach 34°C ± 0.2°C
- (67) 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
- (67) 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
- (67) 4. Air blank completed.
- (67) 5. Calibration check completed. Test results 0.097 AC.
- (67) 6. Air blank completed.
- (67) 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE Logan Moody #7486

INTOXILYZER 8000
Location: YCSD-MAYER
Serial Number: 80-001095
Core version: 8105.48
04/06/2014 14:52:16

Standard Lot#: 32012100A1
Last Changed By: L. MOODY #7486

OAS: L. MOODY #7486
YAUAPA: CSO

Operator: L. MOODY #7486
YAUAPA: CSO

Test	g/210L	Time
Air Blank	0.000	14:53:49
0.100 Cal Check	0.397	14:54:10
Air Blank	0.000	14:54:40

Preventative Maintenance
Performed -- Timer Reset

EXHIBIT G-4
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY

STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000

STANDARD QUALITY ASSURANCE PROCEDURE

QA SPECIALIST LOGAN MOODY #7486 AGENCY CAVAPAI CO. S.D.

DATE 4/6/2014 TIME 1456

INTOXILYZER SERIAL # 80-001095 LOCATION YCSO - MAYER

CS 1. Display Reads "PUSH BUTTON TO START"

DIAGNOSTIC TESTS

- CS 1. Clock time check.
- CS 2. Date check.

OPERATIONAL TESTS

- CS 1. Alcohol-free subject test result 0. 000 AC.
- CS 2. Error recognition logic system functioning.
Not a Successfully Completed Test Sequence printed
- CS 3. Proper sample recognition system.
Not a Successfully Completed Test Sequence printed
Deficient sample printed.
- CS 4. Standard Calibration Check standard 0. 100 AC. Result 0. 097 AC.

Instrument is operating properly and accurately. YES NO

COMMENTS

SIGNATURE Logan Moody 7486

INTOXILYZER 8000
Location: YCSO-MAYER
Serial Number: 80-001095
Core Version: 8105.46
04/06/2014 14:56:38

Standard Lot#: 32012100A1
Last Changed By: L. MOODY #7486

DAS: L. MOODY #7486
YAUAPAI CSO

Operator: L. MOODY #7486
YAUAPAI CSO

Subject: SOAP, TEST
DOB: 01/01/2000
Sex: M weight: 260

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	0.000	14:57:21
Diagnostic Test	Pass	14:57:59
Air Blank	0.000	14:58:26
0.100 Cal Check	0.097	14:58:47
Air Blank	0.000	14:59:16
Subject Test	0.000	14:59:32
Air Blank	ABT*	14:59:42
Air Blank	RFI**	14:59:49

*Sequence Aborted

**RFI Detect

Not a Successfully
Completed
Test Sequence

INTOXILYZER 8000
Location: YCSO-MAYER
Serial Number: 80-001095
Core Version: 8105.48
04/06/2014 15:00:03

Standard Lot#: 3201210091
Last Changed By: L. MOODY #7486

DAS: L. MOODY #7486
YAUAPAI CSO

Operator: L. MOODY #7486
YAUAPAI CSO

Subject: SGAP,TEST
DOB: 01/01/2000
Sex: M Weight: 200

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	0.000	15:00:48
Diagnostic Test	Pass	15:01:25
Air Blank	0.000	15:01:53
0.100 Ca ²⁺ Check	0.097	15:02:13
Air Blank	0.000	15:02:43
Subject Test	0.035*	15:05:50
Air Blank	0.000	15:06:21
Air Blank	0.000	15:06:48
Subject Test	XXX**	15:07:04
Air Blank	IPS***	15:07:11

*Deficient Sample

**Invalid Sample

Begin new deprivation period

***Improper Sample

Not a Successfully
Completed
Test Sequence

EXHIBIT G-2
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY

STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000

STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST LOGAN MORGAN 7446 AGENCY YAVAPAI CO. S.D.

DATE 4/16/2014 TIME 1507

INTOXILYZER SERIAL # 80-001095 LOCATION YCSO MAGER

(✓) 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution 100 AC.
OR

~~Pour a standard alcohol concentration solution _____ AC. into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$~~

(✓) 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"

(✓) 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.

(✓) 4. Air blank completed.

(✓) 5. Calibration check completed. Test results 0.097 AC.

(✓) 6. Air blank completed.

(✓) 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE: LMORGAN 7446

INTOXILYZER 8000
Location: YCSO-MAYER
Serial Number: 80-001095
Core Version: 9105.48
04/06/2014 15:07:28

Standard Lot#: 32012100A1
Last Changed By: L. MOODY #7486

QAS: L. MOODY #7486
YAUAPAI CSO

Operator: L. MOODY #7486
YAUAPAI CSO

Test g/210L Time

Air Blank 0.000 15:07:56
0.100 Cal Check 0.097 15:08:17
Air Blank 0.000 15:08:47

Preventative Maintenance
Performed -- Timer Reset