

EXHIBIT G-2  
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY  
STANDARD QUALITY ASSURANCE PROCEDURES  
INTOXILYZER MODEL 8000

STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST LOGAN MOORE AGENCY YAVAPAI CO. S.O.  
DATE 7/7/2014 TIME 1432  
INTOXILYZER SERIAL # 80-001095 LOCATION YCSO - MAYER

- (67) 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution 100 AC.  
OR  
~~Pour a standard alcohol concentration solution AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach  $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$~~
- (67) 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
- (67) 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
- (67) 4. Air blank completed.
- (67) 5. Calibration check completed. Test results 0. 097 AC.
- (67) 6. Air blank completed.
- (67) 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE L Moore 7486

11:23:19

INTOXICAPZER 8000  
Location: YCSO-MAYER  
Serial Number: 80-001095  
Core Version: 8105.48  
07/07/2014 14:32:36

-----  
Standard Lot#: 32012100A1  
Last Changed By: L. MOODY #7486  
-----

QAS: L. MOODY #7486  
YAUAPAI CSO  
-----

Operator: L. MOODY #7486  
YAUAPAI CSO  
-----

Test	g/210L	Time
Air Blank	0.000	14:33:47
0.100 Cal Check	0.097	14:34:08
Air Blank	0.000	14:34:38

Preventative Maintenance  
Performed -- Timer Reset

EXHIBIT G-4  
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY  
STANDARD QUALITY ASSURANCE PROCEDURES  
INTOXILYZER MODEL 8000

STANDARD QUALITY ASSURANCE PROCEDURE

QA SPECIALIST LOGAN MOODY AGENCY YAVAPAI CO. S.O.  
DATE 7/07/2014 TIME 1435  
INTOXILYZER SERIAL # 80-001095 LOCATION YCSO-MAYER

(M) 1. Display Reads "PUSH BUTTON TO START"

DIAGNOSTIC TESTS

- (M) 1. Clock time check.  
(M) 2. Date check.

OPERATIONAL TESTS

- (M) 1. Alcohol-free subject test result 0, 000 AC.  
(M) 2. Error recognition logic system functioning.  
Not a Successfully Completed Test Sequence printed  
(M) 3. Proper sample recognition system.  
Not a Successfully Completed Test Sequence printed  
Deficient sample printed.  
(M) 4. Standard Calibration Check standard 0, 100 AC. Result 0, 096 AC.

Instrument is operating properly and accurately. YES M NO \_\_\_\_\_

COMMENTS 600 psi  
\_\_\_\_\_  
\_\_\_\_\_  
\_\_\_\_\_

SIGNATURE Logan Moody 7486

INTOXILYZER 8000  
 Location: YCSO-MAYER  
 Serial Number: 80-001095  
 Core Version: 8105.48  
 07/07/2014 14:35:22

-----  
 Standard Lot#: 32012100A1  
 Last Changed By: L. MOODY #7486  
 -----

OAS: L. MOODY #7486  
 YAUAPAI CSO  
 -----

Operator: L. MOODY #7486  
 YAUAPAI CSO  
 -----

Subject: SOAP,TEST  
 DOB: 01/01/2000  
 Sex: M Weight: 200  
 -----

15 Minute Deprivation Period? Yes  
 -----

Test	g/210L	Time
Air Blank	0.000	14:36:17
Diagnostic Test	Pass	14:36:54
Air Blank	0.000	14:37:22
0.100 Cal Check	0.096	14:37:43
Air Blank	0.000	14:38:12
Subject Test	0.000	14:38:34
Air Blank	ABT*	14:38:45
Air Blank	RFI**	14:38:52

\*Sequence Aborted  
 \*\*RFI Detect

Not a Successfully  
 Completed  
 Test Sequence

INTOXILYZER 8000  
 Location: YCSO-MAYER  
 Serial Number: 80-001095  
 Core Version: 8105.48  
 07/07/2014 14:38:59

-----  
 Standard Lot#: 32012100A1  
 Last Changed By: L. MOODY #7486  
 -----

OAS: L. MOODY #7486  
 YAUAPAI CSO  
 -----

Operator: L. MOODY #7486  
 YAUAPAI CSO  
 -----

Subject: SOAP,TEST  
 DOB: 01/01/2000  
 Sex: M Weight: 200  
 -----

15 Minute Deprivation Period? Yes  
 -----

Test	g/210L	Time
Air Blank	0.000	14:40:00
Diagnostic Test	Pass	14:40:37
Air Blank	0.000	14:41:05
0.100 Cal Check	0.096	14:41:25
Air Blank	0.000	14:41:55
Subject Test	0.048*	14:45:02
Air Blank	0.000	14:45:33
Air Blank	0.000	14:46:00
Subject Test	XXX**	14:47:40
Air Blank	IPS***	14:47:47

\*Deficient Sample  
 \*\*Invalid Sample  
 Begin new deprivation period  
 \*\*\*Improper Sample

Not a Successfully  
 Completed  
 Test Sequence

EXHIBIT G-2  
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY  
STANDARD QUALITY ASSURANCE PROCEDURES  
INTOXILYZER MODEL 8000

STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST LOGAN MOORE AGENCY YAVAPAI CO. S.D.  
DATE 7/7/2014 TIME 1448  
INTOXILYZER SERIAL # 80-001095 LOCATION YCSO - MAZER

- (M) 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution 100 AC.  
OR  
~~Pour a standard alcohol concentration solution \_\_\_\_\_ AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach  $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$~~
- (M) 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"  
(M) 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.  
(M) 4. Air blank completed.  
(M) 5. Calibration check completed. Test results 0. 096 AC.  
(M) 6. Air blank completed.  
(M) 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE Logan Moore 7486

INTOXILYZER 8000  
Location: YCSD-MAYER  
Serial Number: 80-001095  
Core Version: 8105.48  
07/07/2014 14:46:19

-----  
Standard Lot#: 32012106A1  
Last Changed By: L. MOODY #7486  
-----

QAS: L. MOODY #7486  
YAUAPAI CSO  
-----

Operator: L. MOODY #7486  
YAUAPAI CSO  
-----

Test	g/210L	Time
Air Blank	0.000	14:46:48
0.100 Cal Check	0.096	14:49:09
Air Blank	0.000	14:49:33

Preventative Maintenance  
Performed -- Timer Reset