

EXHIBIT G-2
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY
STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000

STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST LOGAN MOODY 7486 AGENCY YAVAPAI CO SO.

DATE 10/6/2014 TIME 1639

INTOXILYZER SERIAL # 80-001095 LOCATION YCSO - MAYER

- (M) 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution 100 AC.
OR
~~Pour a standard alcohol concentration solution AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$~~
- (M) 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
- (M) 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
- (M) 4. Air blank completed.
- (M) 5. Calibration check completed. Test results 0.074 AC.
- (M) 6. Air blank completed.
- (M) 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE LMoody7486

DPS Form Exh G-2 (Rev 05-01)

* CAL. CHECK OUT OF TOLERANCE *

EXHIBIT G-2
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY
STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000
STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST LOGAN MOODY 7486 AGENCY YAVAPAI CO. S.O.
DATE 10/6/14 TIME 1642
INTOXILYZER SERIAL # 80-001095 LOCATION YCSO - MAGER

- (M) 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution: 100 AC.
OR
~~Pour a standard alcohol concentration solution _____ AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach 34°C ± 0.2°C~~
- (M) 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
- (M) 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
- (M) 4. Air blank completed.
- (M) 5. Calibration check completed. Test results 0. 0.81 AC.
- (M) 6. Air blank completed.
- (M) 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE Logan Moody

DPS Form Exh G-2 (Rev 05-01)

X CAL CHECK OUT OF TOLERANCE X

EXHIBIT G-2
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY
STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000

STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST LOGAN MOODY 7486 AGENCY CAVAPAE CO. S.O.
DATE 10/6/14 TIME 1646
INTOXILYZER SERIAL # 90-001095 LOCATION CYCSO - MAGER

- (M) 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution 100 AC.
OR
Pour a standard alcohol concentration solution ~~AC~~ into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$
- (M) 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
(M) 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
(M) 4. Air blank completed.
(M) 5. Calibration check completed. Test results 0. 082 AC.
(M) 6. Air blank completed.
(M) 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE LMoody 7486

DPS Form Exh G-2 (Rev 05-01)

* CAL CHECK OUT OF TOLERANCE *
INSTRUMENT REMOVED FROM SERVICE

INTOXILYZER 8000
Location: YCSD-MAYER
Serial Number: 80-001095
Core Version: 8105.48
10/06/2014 16:39:55

Standard Lot#: 32012100A1
Last Changed By: L. MOODY #7486

OAS: L. MOODY #7486
YAUAPAI CSO

Operator: L. MOODY #7486
YAUAPAI CSO

Test	g/210L	Time
Air Blank	0.000	16:40:24
0.100 Cal Check	0.074*	16:40:45
Air Blank	0.000	16:41:13

*Calibration Check Out of Tolerance

INTOXILYZER 8000
Location: YCSD-MAYER
Serial Number: 80-001095
Core Version: 8105.48
10/06/2014 16:42:17

Standard Lot#: 32012100A1
Last Changed By: L. MOODY #7486

OAS: L. MOODY #7486
YAUAPAI CSO

Operator: L. MOODY #7486
YAUAPAI CSO

Test	g/210L	Time
Air Blank	0.000	16:42:47
0.100 Cal Check	0.081*	16:43:07
Air Blank	0.000	16:43:42

*Calibration Check Out of Tolerance

INTOXILYZER 8000
Location: YCSD-MAYER
Serial Number: 80-001095
Core Version: 8105.48
10/06/2014 16:46:07

Standard Lot#: 32012100A1
Last Changed By: L. MOODY #7486

OAS: L. MOODY #7486
YAUAPAI CSO

Operator: L. MOODY #7486
YAUAPAI CSO

Test	g/210L	Time
Air Blank	0.000	16:46:30
0.100 Cal Check	0.082*	16:46:57
Air Blank	0.000	16:47:31

*Calibration Check Out of Tolerance

INTOXILYZER 8000
Location: YC50-MAYER
Serial Number: 80-001095
Core Version: 8185.48
10/07/2014 21:54:52

Standard Lot#: 32012100A1
Last Changed By: L. MOODY #7486

QAS: L. MOODY #7486
YAUAPAI CSO

Operator: L. MOODY #7486
YAUAPAI CSO

Test	g/210L	Time
Air Blank	0.000	21:55:21
0.100 Cal Check	0.097	21:55:42
Air Blank	0.000	21:56:11

Preventative Maintenance
Performed -- Timer Reset

EXHIBIT G-4
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY
STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000

STANDARD QUALITY ASSURANCE PROCEDURE

QA SPECIALIST LOGAN MOORE #7486 AGENCY YAVAPAI CO. S.O.

DATE 10/07/2014 TIME 2138

INTOXILYZER SERIAL # 80-001095 LOCATION YCSO - MAYER

1. Display Reads "PUSH BUTTON TO START"

DIAGNOSTIC TESTS

1. Clock time check.

2. Date check.

OPERATIONAL TESTS

1. Alcohol-free subject test result 0. 000 AC.

2. Error recognition logic system functioning.
Not a Successfully Completed Test Sequence printed

3. Proper sample recognition system.
Not a Successfully Completed Test Sequence printed
Deficient sample printed.

4. Standard Calibration Check standard 0. 100 AC. Result 0. 097 AC.

Instrument is operating properly and accurately. YES NO

COMMENTS

550 psi

INSTRUMENT PLACED BACK IN SERVICE

SIGNATURE LMoore #7486

INTOXILYZER 6000
Location: YCSO-MAYER
Serial Number: 80-001095
Core Version: 6105.48
10/07/2014 21:39:23

Standard Lot#: 32012100A1
Last Changed By: L. MOODY #7486

QAS: L. MOODY #7486
YAUAPAI CSO

Operator: L. MOODY #7486
YAUAPAI CSO

Subject: SOAP,TEST
DOB: 01-01-2000
Sex: M Weight: 150

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	0.000	21:39:18
Diagnostic Test	Pass	21:39:55
Air Blank	0.000	21:40:23
0.100 Cal Check	0.095	21:40:43
Air Blank	0.000	21:41:13
Subject Test	0.000	21:41:31
Air Blank	ABT*	21:41:41
Air Blank	RF1**	21:41:45

*Sequence aborted
**RF1 detect

Not a Successfully
Completed
Test Sequence

INTOXILYZER 6000
Location: YCSO-MAYER
Serial Number: 80-001095
Core Version: 6105.48
10/07/2014 21:42:43

Standard Lot#: 32012100A1
Last Changed By: L. MOODY #7486

QAS: L. MOODY #7486
YAUAPAI CSO

Operator: L. MOODY #7486
YAUAPAI CSO

Subject: SOAP,TEST
DOB: 01-01-2000
Sex: M Weight: 150

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	0.000	21:43:30
Diagnostic Test	Pass	21:44:07
Air Blank	0.000	21:44:35
0.100 Cal Check	0.097	21:44:55
Air Blank	0.000	21:45:25
Subject Test	INT*	21:45:54
Air Blank	0.000	21:46:27

*Interferent Detect

Not a Successfully
Completed
Test Sequence

INTOXILYZER 6000
Location: YCSO-MAYER
Serial Number: 80-001095
Core Version: 6105.48
10/07/2014 21:42:48

Standard Lot#: 32012100A1
Last Changed By: L. MOODY #7486

QAS: L. MOODY #7486
YAUAPAI CSO

Operator: L. MOODY #7486
YAUAPAI CSO

Subject: SOAP,TEST
DOB: 01-01-2000
Sex: M Weight: 150

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	0.000	21:47:34
Diagnostic Test	Pass	21:48:11
Air Blank	0.000	21:48:36
0.100 Cal Check	0.096	21:48:59
Air Blank	0.000	21:49:29
Subject Test	0.074*	21:52:35
Air Blank	0.000	21:53:16
Air Blank	0.000	21:53:35
Subject Test	XXX**	21:53:55
Air Blank	IFS***	21:54:07

*Deficient Sample
**Invalid Sample
Begin new deprivation period
***Improper Sample

Not a Successfully
Completed
Test Sequence

EXHIBIT G-2
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY
STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000
STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST LOGAN MOODY #7486 AGENCY YAVAPAI CO. S.O.
DATE 10/07/2014 TIME 2136
INTOXILYZER SERIAL # 80-001095 LOCATION YCSO - MAYER

- CM 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution .100 AC.
OR
~~Pour a standard alcohol concentration solution _____ AC into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$~~
- CM 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
- CM 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
- CM 4. Air blank completed.
- CM 5. Calibration check completed. Test results 0.096 AC.
- CM 6. Air blank completed.
- CM 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE Logan Moody #7486

INTOXILYZER 8000
Location: YCSD-DAYER
Serial Number: 00-001095
Core Version: 8105.48
10/07/2014 21:36:20

Standard Lot#: 32012100A1
Last Changed By: L. MOODY #7486

QAS: L. MOODY #7486
YAUAPAI CSO

Operator: L. MOODY #7486
YAUAPAI CSO

Test	g/210L	Time
Air Blank	0.000	21:36:51
0.100 Cal Check	0.096	21:37:12
Air Blank	0.000	21:37:41

Preventative Maintenance
Performed -- Timer Reset