

EXHIBIT G-4
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY
STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000

STANDARD QUALITY ASSURANCE PROCEDURE

QA SPECIALIST S. MARTIN #20 AGENCY CUMD
DATE 2-2-00 TIME 1406
INTOXILYZER SERIAL # 030-001096 LOCATION CUMD

() 1. Display Reads "PUSH BUTTON TO START"

DIAGNOSTIC TESTS

- () 1. Clock time check.
() 2. Date check.

OPERATIONAL TESTS

- () 1. Alcohol-free subject test result 0. 000 AC.
() 2. Error recognition logic system functioning.
Not a Successfully Completed Test Sequence printed
() 3. Proper sample recognition system.
Not a Successfully Completed Test Sequence printed
Deficient sample printed.
() 4. Standard Calibration Check standard 0. 100 AC. Result 0. 090 AC.

Instrument is operating properly and accurately. YES NO

COMMENTS

SIGNATURE S. MARTIN #20

LOT #706502 I
EXP 3-9-01

INTOXILYZER 8000

Location: CVMO

Serial Number: 80-001096

02/02/2008

14:06:25

Standard Lot#: 706502I

Last Changed By: S. MARTINKO #25

QAS: S. MARTINKO #25

CAMP VERDE MO

Test	g/210L	Time
Air Blank	0.000	14:06:52
Cal Check	0.098	14:07:11
Air Blank	0.000	14:07:42

S. MARTINKO #25

Operator's Signature

INTOXILYZER 8000
Location: CVMO
Serial Number : 80-001096
02/02/2008

14:07:46

Standard Lot#: 706502I
Last Changed By: S. MARTINKO #25

QAS: S. MARTINKO #25
CAMP VERDE MO

Operator: S. MARTINKO #25
CAMP VERDE MO

Subject: SQAP,TEST,
DOB: 10/30/1974
Sex: M Weight: 240

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	0.000	14:08:40
Diagnostic Test	Pass	14:09:07
Air Blank	0.000	14:09:35
0.100 Cal Check	0.098	14:09:53
Air Blank	0.000	14:10:25
Subject Test	0.000	14:10:48
Air Blank	0.000	14:11:16
Five Minute Wait	ABT*	
Air Blank	0.000	14:11:50

*Sequence Aborted

Not a Successfully Completed Test Sequence

INTOXILYZER 8000

Location: CVMO

Serial Number : 80-001096

02/02/2008

14:15:21

Standard Lot#: 706502I

Last Changed By: S. MARTINKO #25

OAS: S. MARTINKO #25

CAMP VERDE MO

Operator: S. MARTINKO #25

CAMP VERDE MO

Subject: SQAP,TEST,

DOB: 10/30/1974

Sex: M

Weight: 240

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	IPS*	14:15:58
Air Blank	RFI**	14:16:05

*Improper Sample

**RFI Detect

Not a Successfully Completed Test Sequence

INTOXILYZER 8000

Location: CVMO

Serial Number : 80-001096

02/02/2008

14:16:58

Standard Lot#: 706502I

Last Changed By: S. MARTINKO #25

QAS: S. MARTINKO #25

CAMP VERDE MO

Operator: S. MARTINKO #25

CAMP VERDE MO

Subject: SOAP, TEST,

DOB: 10/30/1974

Sex: M

Weight: 240

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	0.000	14:18:11
Diagnostic Test	Pass	14:18:38
Air Blank	0.000	14:19:05
0.100 Cal Check	0.098	14:19:24
Air Blank	0.000	14:19:56
Subject Test	0.202*	14:23:04
Air Blank	0.000	14:23:36
Air Blank	0.000	14:24:04
Subject Test	XXX**	14:24:28
Air Blank	0.000	14:25:02

*Deficient Sample

**Invalid Sample

Begin new deprivation period

Not a Successfully Completed Test Sequence

EXHIBIT G-2
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY
STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000
STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST S. MARTINEZ AGENCY CUMD
DATE 2-2-09 TIME 1406
INTOXILYZER SERIAL # 800-001096 LOCATION CUMD

1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution 100% AC.
OR
Pour a standard alcohol concentration solution _____ AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$
2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
4. Air blank completed.
5. Calibration check completed. Test results 0. 090 AC.
6. Air blank completed.
7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE S. MARTINEZ #21

DPS Form Exh G-2 (Rev 05-01)

LOT # 706502 I
EXP 3-9-09

INTOXILYZER 8000
Location: CVMO
Serial Number: 80-001096
02/02/2008 14:25:26

Standard Lot#: 7065021
Last Changed By: S. MARTINKO #25

QAS: S. MARTINKO #25
CAMP VERDE MO

Test	g/210L	Time
Air Blank	0.000	14:25:53
Cal Check	0.096	14:26:12
Air Blank	0.000	14:26:43

S. Martinko

Operator's Signature

EXHIBIT G-2
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY
STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000
STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST S. MARTINEK AGENCY CUMD
DATE 2-2-04 TIME 1359
INTOXILYZER SERIAL # 90-001736 LOCATION CUMD

1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution 100 AC.
OR
Pour a standard alcohol concentration solution _____ AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$
2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
4. Air blank completed.
5. Calibration check completed. Test results 0. 090 AC.
6. Air blank completed.
7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE S. MARTINEK #21

DPS Form Exh G-2 (Rev 05-01)

LOT A706562 I
EXP. 3-9-04

INTOXILYZER 8000
Location: DPS
Serial Number : 80-001736
02/02/2008

14:01:19

Standard Lot#: 706502I
Last Changed By: S. MARTINKO #25

QAS: S. MARTINKO #25
CAMP VERDE MO

Operator: S. MARTINKO #25
CAMP VERDE MO

Subject: SQAP,TEST,
DOB: 10/30/1974
Sex: M Weight: 240

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	0.000	14:02:11
Diagnostic Test	Pass	14:02:38
Air Blank	0.000	14:03:06
0.100 Cal Check	0.097	14:03:24
Air Blank	0.000	14:03:54
Subject Test	0.000	14:04:08
Air Blank	0.000	14:04:36
Five Minute Wait	ABT*	
Air Blank	0.000	14:05:10

*Sequence Aborted

Not a Successfully Completed Test Sequence

INTOXILYZER 8000
Location: DPS
Serial Number : 80-001736
02/02/2008 14:09:09

Standard Lot#: 706502I
Last Changed By: S. MARTINKO #25

QAS: S. MARTINKO #25
CAMP VERDE MO

Operator: S. MARTINKO #25
CAMP VERDE MO

Subject: SQAP, TEST,
DOB: 10/30/1974
Sex: M Weight: 240

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	IPS*	14:09:49
Air Blank	RFI**	14:09:57

*Improper Sample

**RFI Detect

Not a Successfully Completed Test Sequence

INTOXILYZER 8000
Location: DPS
Serial Number : 80-001736
02/02/2008

14:13:25

Standard Lot#: 706502I
Last Changed By: S. MARTINKO #25

QAS: S. MARTINKO #25
CAMP VERDE MO

Operator: S. MARTINKO #25
CAMP VERDE MO

Subject: SQAP,TEST,
DOB: 10/30/1974
Sex: M Weight: 240

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	0.000	14:14:38
Diagnostic Test	Pass	14:15:06
Air Blank	0.000	14:15:33
0.100 Cal Check	0.098	14:15:52
Air Blank	0.000	14:16:21
Subject Test	0.125*	14:19:30
Air Blank	0.000	14:20:01
Air Blank	0.000	14:20:29
Subject Test	XXX**	14:20:44
Air Blank	0.000	14:21:17

*Deficient Sample

**Invalid Sample

Begin new deprivation period

Not a Successfully Completed Test Sequence

EXHIBIT G-2
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY
STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000
STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST S. MARTINKO AGENCY CUMD
DATE 2-2-09 TIME 1359
INTOXILYZER SERIAL # 000-001736 LOCATION CUMD

- (/) 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution 100 AC.
OR
Pour a standard alcohol concentration solution _____ AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$
- (/) 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
- (/) 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
- (/) 4. Air blank completed.
- (/) 5. Calibration check completed. Test results 0. 0.06 AC.
- (/) 6. Air blank completed.
- (/) 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE S. MARTINKO

DPS Form Exh G-2 (Rev 05-01)

LOT # 906501 I
Exp. 2-9-09

INTOXILYZER 8000

Location: DPS

Serial Number: 80-001736

02/02/2008

14:21:36

Standard Lot#: 706502I

Last Changed By: S. MARTINKO #25

QAS: S. MARTINKO #25

CAMP VERDE MO

Test	g/210L	Time
Air Blank	0.000	14:22:04
Cal Check	0.096	14:22:22
Air Blank	0.000	14:22:52

S. MARTINKO #25
Operator's Signature