

EXHIBIT G-2
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY

STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000

STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST Russell Grover AGENCY CVMO
DATE 7-29-2008 TIME 1133
INTOXILYZER SERIAL # 80-001096 DPS LOCATION CVMO

- (✓) 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution .100 AC.
OR
Pour a standard alcohol concentration solution _____ AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach 34°C ± 0.2°C
- (✓) 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
- (✓) 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
- (✓) 4. Air blank completed.
- (✓) 5. Calibration check completed. Test results 0. 096 AC.
- (✓) 6. Air blank completed.
- (✓) 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE Russell Grover #22

DPS Form Exh G-2 (Rev 05-01)

Lot # 706502 I

EXP DATE 3-9-09

[Faint, illegible text at the bottom of the page]


INTOXILYZER 8000
Location: CVMO
Serial Number: 80-001096
07/29/2008

11:32:18

Standard Lot#: 706502I
Last Changed By: S. MARTINKO #25

QAS: R. GROVER #22
CAMP VERDE MO

Test	g/210L	Time
Air Blank	0.000	11:32:46
Cal Check	0.096	11:33:05
Air Blank	0.000	11:33:37



Operator's Signature

EXHIBIT G-2
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY
STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000
STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST Russell Grover AGENCY C.V.M.O.
DATE 7-29-08 TIME 18:16
INTOXILYZER SERIAL # 80-00173608 LOCATION C.V.M.O.

- (✓) 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution -100 AC.
OR
Pour a standard alcohol concentration solution _____ AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach 34°C ± 0.2°C.
- (✓) 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
- (✓) 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
- (✓) 4. Air blank completed.
- (✓) 5. Calibration check completed. Test results 0. 095 AC.
- (✓) 6. Air blank completed.
- () 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE R. Grover #22

DPS Form Exh G-2 (Rev 05-01)

Lot # 706502 I

EXP DATE 3-9-09

[Faint, illegible text at the bottom of the page]


INTOXILYZER 8000
Location: DPS
Serial Number: 80-001736
07/29/2008

18:15:36

Standard Lot#: 706502I
Last Changed By: S. MARTINKO #25

QAS: R. GROVER #22
CAMP VERDE MO

Test	g/210L	Time
Air Blank	0.000	18:16:04
Cal Check	0.095	18:16:22
Air Blank	0.000	18:16:52



Operator's Signature