

EXHIBIT G-2  
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY  
STANDARD QUALITY ASSURANCE PROCEDURES  
INTOXILYZER MODEL 8000  
STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST Russell Grover AGENCY C.V.M.O.  
DATE 8-29-09 TIME 10:35  
INTOXILYZER SERIAL # 80-001096 DPS LOCATION C.V.M.O.

- (✓) 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution .100 AC.  
OR  
Pour a standard alcohol concentration solution \_\_\_\_\_ AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach  $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$
- (✓) 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
- (✓) 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
- (✓) 4. Air blank completed.
- (✓) 5. Calibration check completed. Test results 0. 096 AC.
- (✓) 6. Air blank completed.
- (✓) 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE Russell Grover #22

DPS Form Exh G-2 (Rev 05-01)

Lot # 706502 I

EXP DATE 3-9-09

EXHIBIT G-4  
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY  
STANDARD QUALITY ASSURANCE PROCEDURES  
INTOXILYZER MODEL 8000

STANDARD QUALITY ASSURANCE PROCEDURE

QA SPECIALIST Russell Grover AGENCY CVMO  
DATE 8-29-08 TIME 10:40  
INTOXILYZER SERIAL # 80-001096 OPS LOCATION C.V.M.O.

() 1. Display Reads "PUSH BUTTON TO START"

DIAGNOSTIC TESTS

- () 1. Clock time check.  
() 2. Date check.

OPERATIONAL TESTS

- () 1. Alcohol-free subject test result 0. 000 AC.  
() 2. Error recognition logic system functioning.  
Not a Successfully Completed Test Sequence printed  
() 3. Proper sample recognition system.  
Not a Successfully Completed Test Sequence printed  
Deficient sample printed.  
() 4. Standard Calibration Check standard 0. 100 AC. Result 0. 094 AC.

Instrument is operating properly and accurately. YES  NO

COMMENTS

NONE

SIGNATURE Russell Grover



INTOXILYZER 8000  
Location: CVMO  
Serial Number : 80-001096  
08/29/2008

10:37:52

Standard Lot#: 706502I  
Last Changed By: S. MARTINKO #25

QAS: R. GROVER #22  
CAMP VERDE MO

Operator: R. GROVER #22  
CAMP VERDE MO

Subject: SQAP,TEST

DOB: 08/28/2008

Sex: M

Weight: 225

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	0.000	10:39:06
Diagnostic Test	Pass	10:39:33
Air Blank	0.000	10:40:00
0.100 Cal Check	0.095	10:40:19
Air Blank	0.000	10:40:51
Subject Test	0.000	10:41:07
Air Blank	0.000	10:41:35
Five Minute Wait	ABT*	
Air Blank	0.000	10:42:12

\*Sequence Aborted

Not a Successfully Completed Test Sequence

*R. Grover #22 8/29/08*

INTOXILYZER 8000

Location: CVMO

Serial Number : 80-001096

08/29/2008

10:43:23

Standard Lot#: 706502I

Last Changed By: S. MARTINKO #25

OAS: R. GROVER #22

CAMP VERDE MO

Operator: R. GROVER #22

CAMP VERDE MO

Subject: SQAP,TEST

DOB: 08/28/2008

Sex: M

Weight: 225

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	IPS*	10:44:11
Air Blank	0.000	10:44:39

\*Improper Sample

Not a Successfully Completed Test Sequence

*R. Grover #22 8/29/08*

INTOXILYZER 8000  
Location: CVMO  
Serial Number : 80-001096  
08/29/2008

10:45:20

Standard Lot#: 706502I  
Last Changed By: S. MARTINKO #25

QAS: R. GROVER #22  
CAMP VERDE MO

Operator: R. GROVER #22  
CAMP VERDE MO

Subject: SQAP, TEST

DOB: 08/28/2008

Sex: M

Weight: 225

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	RFI*	10:45:59
Air Blank	0.000	10:46:27

\*RFI Detect

Not a Successfully Completed Test Sequence

*Ronald Brown #22*

*8/29/08*

INTOXILYZER 8000

Location: CVMO

Serial Number : 80-001096

08/29/2008

10:47:51

Standard Lot#: 706502I

Last Changed By: S. MARTINKO #25

QAS: R. GROVER #22

CAMP VERDE MO

Operator: R. GROVER #22

CAMP VERDE MO

Subject: SQAP,TEST

DOB: 08/28/2008

Sex: M

Weight: 225

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	0.000	10:48:45
Diagnostic Test	Pass	10:49:13
Air Blank	0.000	10:49:40
0.100 Cal Check	0.096	10:49:59
Air Blank	0.000	10:50:31
Subject Test	0.064*	10:53:39
Air Blank	0.000	10:54:13
Air Blank	0.000	10:54:41
Subject Test	XXX**	10:54:58
Air Blank	0.000	10:55:31

\*Deficient Sample

\*\*Invalid Sample

Begin new deprivation period

Not a Successfully Completed Test Sequence

*Ronald Grover #22*

*8/29/08*



EXHIBIT G-2  
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY  
STANDARD QUALITY ASSURANCE PROCEDURES  
INTOXILYZER MODEL 8000  
STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST Russell Grover AGENCY C.V.M.O.  
DATE 8-29-08 TIME 11:11  
INTOXILYZER SERIAL # 80-001736 05 LOCATION C.V.M.O.

- (✓) 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution .100 AC.  
OR  
Pour a standard alcohol concentration solution \_\_\_\_\_ AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach  $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$
- (✓) 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
- (✓) 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
- (✓) 4. Air blank completed.
- (✓) 5. Calibration check completed. Test results 0. 096 AC.
- (✓) 6. Air blank completed.
- (✓) 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE Russell Grover #22 8/29/08

DPS Form Exh G-2 (Rev 05-01)

Lot # 706502 I

EXP DATE 3-9-09

EXHIBIT G-4  
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY

STANDARD QUALITY ASSURANCE PROCEDURES  
INTOXILYZER MODEL 8000

STANDARD QUALITY ASSURANCE PROCEDURE

QA SPECIALIST Russell Grover AGENCY C.V.M.O.  
DATE 8-29-08 TIME 11:18  
INTOXILYZER SERIAL # 80-001736 DPS LOCATION C.V.M.O.

() 1. Display Reads "PUSH BUTTON TO START"

DIAGNOSTIC TESTS

- () 1. Clock time check.  
() 2. Date check.

OPERATIONAL TESTS

- () 1. Alcohol-free subject test result 0. 000 AC.  
() 2. Error recognition logic system functioning.  
Not a Successfully Completed Test Sequence printed  
() 3. Proper sample recognition system.  
Not a Successfully Completed Test Sequence printed  
Deficient sample printed.  
() 4. Standard Calibration Check standard 0. 100 AC. Result 0. 094 AC.

Instrument is operating properly and accurately. YES  NO

COMMENTS

NONE

SIGNATURE Russell Grover #22

DPS Form Exh G-4 (Rev 05-01)

Lot # 706502 I

EXP DATE 3-9-09

INTOXILYZER 8000  
Location: DPS  
Serial Number: 80-001736  
08/29/2008 11:10:54

-----  
Standard Lot#: 706502I  
Last Changed By: S. MARTINKO #25  
-----

QAS: R. GROVER #22  
CAMP VERDE MO  
-----

Test	g/210L	Time
Air Blank	0.000	11:11:21
Cal Check	0.096	11:11:40
Air Blank	0.000	11:12:09

*R. Grover #22 8/29/08*  
-----  
Operator's Signature

Lot # 706502 I

EXP DATE 3-9-09

INTOXILYZER 8000  
Location: DPS  
Serial Number : 80-001736  
08/29/2008 11:15:21

---

Standard Lot#: 706502I  
Last Changed By: S. MARTINKO #25

---

QAS: R. GROVER #22  
CAMP VERDE MO

---

Operator: R. GROVER #22  
CAMP VERDE MO

---

Subject: SQAP, TEST  
DOB: 08/29/2008  
Sex: M Weight: 225

---

15 Minute Deprivation Period? Yes

---

Test	g/210L	Time
Air Blank	0.000	11:16:19
Diagnostic Test	Pass	11:16:46
Air Blank	0.000	11:17:14
0.100 Cal Check	0.096	11:17:32
Air Blank	0.000	11:18:02
Subject Test	0.000	11:18:20
Air Blank	0.000	11:18:48
Five Minute Wait	ABT*	
Air Blank	0.000	11:19:46

\*Sequence Aborted

**Not a Successfully Completed Test Sequence**

*Russell Shaw #22*

*8/29/08*

INTOXILYZER 8000  
Location: DPS  
Serial Number : 80-001736  
08/29/2008 11:22:04

---

Standard Lot#: 706502I  
Last Changed By: S. MARTINKO #25

---

QAS: R. GROVER #22  
CAMP VERDE MO

---

Operator: R. GROVER #22  
CAMP VERDE MO

---

Subject: SQAP, TEST  
DOB: 08/29/2008  
Sex: M Weight: 225

---

15 Minute Deprivation Period? Yes

---

Test	g/210L	Time
Air Blank	IPS*	11:22:54
Air Blank	0.000	11:23:22

\*Improper Sample

Not a Successfully Completed Test Sequence

*Ramon M #22*

*8/29/08*

INTOXILYZER 8000  
Location: DPS  
Serial Number : 80-001736  
08/29/2008 11:24:55

---

Standard Lot#: 706502I  
Last Changed By: S. MARTINKO #25

---

QAS: R. GROVER #22  
CAMP VERDE MO

---

Operator: R. GROVER #22  
CAMP VERDE MO

---

Subject: SOAP,TEST  
DOB: 08/29/2008  
Sex: M Weight: 225

---

15 Minute Deprivation Period? Yes

---

Test	g/210L	Time
Air Blank	RFI*	11:25:36
Air Blank	0.000	11:26:05

\*RFI Detect

**Not a Successfully Completed Test Sequence**

*Small the #22*

*8/29/08*

INTOXILYZER 8000  
Location: DPS  
Serial Number : 80-001736  
08/29/2008

11:30:15

---

Standard Lot#: 706502I  
Last Changed By: S. MARTINKO #25

---

QAS: R. GROVER #22  
CAMP VERDE MO

---

Operator: R. GROVER #22  
CAMP VERDE MO

---

Subject: SQAP,TEST  
DOB: 08/29/2008  
Sex: M Weight: 225

---

15 Minute Deprivation Period? Yes

---

Test	g/210L	Time
Air Blank	0.000	11:31:04
Diagnostic Test	Pass	11:31:31
Air Blank	0.000	11:31:58
0.100 Cal Check	0.097	11:32:17
Air Blank	0.000	11:32:46
Subject Test	0.248*	11:35:55
Air Blank	0.000	11:36:30
Air Blank	0.000	11:36:57
Subject Test	XXX**	11:37:14
Air Blank	0.000	11:37:48

\*Deficient Sample

\*\*Invalid Sample

Begin new deprivation period

**Not a Successfully Completed Test Sequence**

*Ramon M... #22*


*8/29/08*

INTOXILYZER 8000  
Location: DPS  
Serial Number: 80-001736  
08/29/2008 11:40:10

-----  
Standard Lot#: 706502I  
Last Changed By: S. MARTINKO #25  
-----

QAS: R. GROVER #22  
CAMP VERDE MO  
-----

Test	g/210L	Time
Air Blank	0.000	11:40:37
Cal Check	0.094	11:40:56
Air Blank	0.000	11:41:25

 8/29/08  
-----  
Operator's Signature

Lot # 706502 I

EXP DATE 3-9-09