

EXHIBIT G-2  
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY  
STANDARD QUALITY ASSURANCE PROCEDURES  
INTOXILYZER MODEL 8000  
STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST Russell Crower AGENCY C.V.M.O.  
DATE 9-29-08 TIME 10:22  
INTOXILYZER SERIAL # 80-001096 DPS LOCATION C.V.M.O.

- (✓) 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution 0.100 AC.  
OR  
Pour a standard alcohol concentration solution \_\_\_\_\_ AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach  $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$
- (✓) 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
- (✓) 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
- (✓) 4. Air blank completed.
- (✓) 5. Calibration check completed. Test results 0. 097 AC.
- (✓) 6. Air blank completed.
- (✓) 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE Russell Crower #22

DPS Form Exh G-2 (Rev 05-01)

Lot # 706502-I  
EXP DATE: 3-9-09

INTOXILYZER 8000

Location: CVMO

Serial Number: 80-001096

09/29/2008

10:22:03


-----  
Standard Lot#: 706502I

Last Changed By: S. MARTINKO #25

-----  
QAS: R. GROVER #22

CAMP VERDE MO

-----  
Test                      g/210L                      Time  
-----  
Air Blank                  0.000                      10:22:30  
Cal Check                  0.097                      10:22:49  
Air Blank                  0.000                      10:23:21



-----  
Operator's Signature

EXHIBIT G-2  
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY  
STANDARD QUALITY ASSURANCE PROCEDURES  
INTOXILYZER MODEL 8000  
STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST Russell Wover AGENCY C.V.M.O.  
DATE 9-29-08 TIME 10:29  
INTOXILYZER SERIAL # 80-001736 DS LOCATION C.V.M.O.

- (✓) 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution 1100 AC.  
OR  
Pour a standard alcohol concentration solution \_\_\_\_\_ AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach  $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$
- (✓) 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
- (✓) 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
- (✓) 4. Air blank completed.
- (✓) 5. Calibration check completed. Test results 0. 097 AC.
- (✓) 6. Air blank completed.
- (✓) 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE Russell Wover #22

DPS Form Exh G-2 (Rev 05-01)

Lot # 706502 I  
EXP DATE: 3-9-09

INTOXILYZER 8000

Location: DPS

Serial Number: 80-001736

09/29/2008

10:28:53


-----  
Standard Lot#: 706502I

Last Changed By: S. MARTINKO #25

-----  
QAS: R. GROVER #22

CAMP VERDE MO

-----  
Test                      g/210L                      Time  
-----  
Air Blank                  0.000                      10:29:21  
Cal Check                  0.097                      10:29:39  
Air Blank                  0.000                      10:30:09



Operator's Signature