

EXHIBIT G-2  
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY  
STANDARD QUALITY ASSURANCE PROCEDURES  
INTOXILYZER MODEL 8000  
STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST Russell Greer AGENCY C.V.M.O.  
DATE 2-24-09 TIME 1321  
INTOXILYZER SERIAL # 90-001096 DPS LOCATION C.V.M.O.

- (✓) 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution 100 AC.  
OR  
Pour a standard alcohol concentration solution \_\_\_\_\_ AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach  $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$
- (✓) 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
- (✓) 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
- (✓) 4. Air blank completed.
- (✓) 5. Calibration check completed. Test results 0. 096 AC.
- (✓) 6. Air blank completed.
- (✓) 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE [Signature] #22

DPS Form Exh G-2 (Rev 05-01)

Lot # 706502 I  
Exp Date 3-9-09

*[Faint, illegible text at the bottom of the page]*

EXHIBIT G-4  
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY  
STANDARD QUALITY ASSURANCE PROCEDURES  
INTOXILYZER MODEL 8000

STANDARD QUALITY ASSURANCE PROCEDURE

QA SPECIALIST Roman Rivera AGENCY C.V.M.O.  
DATE 2-24-09 TIME 1325  
INTOXILYZER SERIAL # 86-001096 DPS LOCATION C.V.M.O.

() 1. Display Reads "PUSH BUTTON TO START"

DIAGNOSTIC TESTS

- () 1. Clock time check.  
() 2. Date check.

OPERATIONAL TESTS

- () 1. Alcohol-free subject test result 0. 000 AC.  
() 2. Error recognition logic system functioning.  
Not a Successfully Completed Test Sequence printed  
() 3. Proper sample recognition system.  
Not a Successfully Completed Test Sequence printed  
Deficient sample printed.  
() 4. Standard Calibration Check standard 0. 100 AC. Result 0. 096 AC.

Instrument is operating properly and accurately. YES  NO

COMMENTS

NONE

SIGNATURE Roman Rivera #22

INTOXILYZER 8000

Location: CVMO

Serial Number: 89-001096

02/24/2009

13:20:47

-----  
Standard Lot#: 7065021

Last Changed By: S. MARTINEZ #26

-----  
OAS: K. GROVER #22

CAMP VERDE MO

-----  
Test                      g/210L                      Time  
-----  
Air Blank                  0.000                      13:21:10  
Cal Check                  0.096                      13:21:29  
Air Blank                  0.000                      13:22:01



-----  
Operator's Signature

INTOXILYZER 8000  
Location: CVMO  
Serial Number : 80-001096  
02/24/2009

13:23:46

Standard Lot#: 706502I  
Last Changed By: S. MARTINKO #25

QAS: R. GROVER #22  
CAMP VERDE MO

Operator: R. GROVER #22  
CAMP VERDE MO

Subject: SOAP,TEST

DOB: 08/19/1971

Sex: M Weight: 225

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	0.000	13:24:45
Diagnostic Test	Pass	13:25:13
Air Blank	0.000	13:25:40
0.100 Cal Check	0.096	13:25:59
Air Blank	0.000	13:26:30
Subject Test	0.000	13:26:49
Air Blank	0.000	13:27:16
Five Minute Wait	ABT*	
Air Blank	0.000	13:27:59

\*Sequence Aborted

Not a Successfully Completed Test Sequence

*Russell M. #22*

INTOXILYZER 8000

Location: CVMO

Serial Number : 80-001096

02/24/2009

13:28:54

Standard Lot#: 7065021

Last Changed By: S. MARTINKO

#25

OAS: R. GROVER

#22

CAMP VERDE MO

Operator: R. GROVER

#22

CAMP VERDE MO

Subject: SQAP, TEST

DOB: 08/19/1971

Sex: M

Weight: 225

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	0.000	13:29:44
Diagnostic Test	Pass	13:30:12
Air Blank	0.000	13:30:39
0.100 Cal Check	0.096	13:30:58
Air Blank	IPS*	13:31:05
Air Blank	0.000	13:31:34

\*Improper Sample

Not a Successfully Completed Test Sequence

*Ronald M*

INTOXILYZER 8000

Location: CVMO

Serial Number : 80-001096

02/24/2009

13:33:33

Standard Lot#: 706502I

Last Changed By: S. MARTINKO #25

QAS: R. GROVER #22

CAMP VERDE MO

Operator: R. GROVER #22

CAMP VERDE MO

Subject: SOAP,TEST

DOB: 08/19/1971

Sex: M

Weight: 225

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	RFI*	13:34:38
Air Blank	0.000	13:34:46

\*RFI Detect

*R. Grover #22*

Not a Successfully Completed Test Sequence

INTOXILYZER 8000

Location: CVM0

Serial Number : 80-00109a

02/24/2009

13:38:38

Standard Lot#: 7065021

Last Changed By: R. MARJANO

431

OAS: R. GROVER

427

CAMP VERDE MO

Operator: R. GROVER

422

CAMP VERDE MO

Subject: SOBE, TRIST

DOB: 08/19/1971

Sex: M

Weight: 225

15 Minute Deprivation Period? Yes

Test	g/210l	Time
Air Blank	0.000	13:36:37
Diagnostic Test	Pass	13:36:38
Air Blank	0.000	13:37:02
0.100 Cal Check	0.057	13:37:21
Air Blank	0.000	13:37:47
Subject Test	0.244*	13:38:01
Air Blank	0.000	13:38:16
Air Blank	0.000	13:38:30
Subject Test	0.244*	13:38:45
Air Blank	0.000	13:38:59

\*Deficient Sample

\*\*Invalid Sample

Begin new deprivation period

Not a Successfully Completed Test Sequence

*R. Grover*



EXHIBIT G-2  
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY  
STANDARD QUALITY ASSURANCE PROCEDURES  
INTOXILYZER MODEL 8000  
STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST Russell Groves AGENCY C.V.M.O.  
DATE 2-24-09 TIME 1356  
INTOXILYZER SERIAL # 80-001736 D15 LOCATION C.V.M.O.

- (✓) 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution .100 AC.  
OR  
Pour a standard alcohol concentration solution \_\_\_\_\_ AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach  $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$
- (✓) 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
- (✓) 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
- (✓) 4. Air blank completed.
- (✓) 5. Calibration check completed. Test results 0.096 AC.
- (✓) 6. Air blank completed.
- (✓) 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE Russell Groves #22

DPS Form Exh G-2 (Rev 05-01)

Lot # 706502 I

EH DATE 3-9-09

EXHIBIT G-4  
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY  
STANDARD QUALITY ASSURANCE PROCEDURES  
INTOXILYZER MODEL 8000  
STANDARD QUALITY ASSURANCE PROCEDURE

QA SPECIALIST Russell Grever AGENCY C.V.M.O.  
DATE 2-24-09 TIME 1402  
INTOXILYZER SERIAL # 80-001736 DPS LOCATION C.V.M.O.

() 1. Display Reads "PUSH BUTTON TO START"

DIAGNOSTIC TESTS

- () 1. Clock time check.  
() 2. Date check.

OPERATIONAL TESTS

- () 1. Alcohol-free subject test result 0. 000 AC.  
() 2. Error recognition logic system functioning.  
Not a Successfully Completed Test Sequence printed  
() 3. Proper sample recognition system.  
Not a Successfully Completed Test Sequence printed  
Deficient sample printed.  
() 4. Standard Calibration Check standard 0. 100 AC. Result 0. 097 AC.

Instrument is operating properly and accurately. YES  NO

COMMENTS

none

SIGNATURE Russell Grever

INTOXILYZER 8000

Location: DPS

Serial Number: 80-001736

02/24/2009

13:56:10

-----  
Standard Lot#: 706502I

Last Changed By: S. MARTINKO #25

-----  
QAS: R. GROVER #22

CAMP VERDE MO

-----  
Test                      g/210L                      Time  
-----  
Air Blank                  0.000                      13:56:38  
Cal Check                  0.096                      13:56:56  
Air Blank                  0.000                      13:57:26



-----  
Operator's Signature

Lot # 706502 I

EXP DATE 3-9-09

INTOXILYZER 8000  
Location: DPS  
Serial Number : 80-001736  
02/24/2009 13:59:36

---

Standard Lot#: 706502I  
Last Changed By: S. MARTINKO #25

---

QAS: R. GROVER #22  
CAMP VERDE MO

---

Operator: R. GROVER #22  
CAMP VERDE MO

---

Subject: SQAP,TEST  
DOB: 08/19/1971  
Sex: M Weight: 225

---

15 Minute Deprivation Period? Yes

---

Test	g/210L	Time
Air Blank	0.000	14:00:22
Diagnostic Test	Pass	14:00:49
Air Blank	0.000	14:01:17
0.100 Cal Check	0.097	14:01:36
Air Blank	0.000	14:02:05
Subject Test	0.000	14:02:20
Air Blank	0.000	14:02:47
Five Minute Wait	ABT*	
Air Blank	0.000	14:03:28

\*Sequence Aborted

**Not a Successfully Completed Test Sequence**

*Ronald M. #22*

INTOXILYZER 8000  
Location: DPS  
Serial Number : 80-001736  
02/24/2009 14:05:10

---

Standard Lot#: 706502I  
Last Changed By: S. MARTINKO #25

---

QAS: R. GROVER #22  
CAMP VERDE MO

---

Operator: R. GROVER #22  
CAMP VERDE MO

---

Subject: SQAP,TEST  
DOB: 08/19/1971  
Sex: M Weight: 225

---

15 Minute Deprivation Period? Yes

---

Test	g/210L	Time
Air Blank	IPS*	14:05:42
Air Blank	0.000	14:06:11

\*Improper Sample

**Not a Successfully Completed Test Sequence**

*Amelie L. H...*

INTOXILYZER 8000  
Location: DPS  
Serial Number : 80-001736  
02/24/2009 14:06:46

---

Standard Lot#: 706502I  
Last Changed By: S. MARTINKO #25

---

QAS: R. GROVER #22  
CAMP VERDE MO

---

Operator: R. GROVER #22  
CAMP VERDE MO

---

Subject: SQAP,TEST  
DOB: 08/19/1971  
Sex: M Weight: 225

---

15 Minute Deprivation Period? Yes

---

Test	g/210L	Time
Air Blank	RFI*	14:07:16
Air Blank	0.000	14:07:44

\*RFI Detect

**Not a Successfully Completed Test Sequence**

*Randy Martin*

INTOXILYZER 8000  
Location: DPS  
Serial Number : 80-001736  
02/24/2009 14:08:15

---

Standard Lot#: 706502I  
Last Changed By: S. MARTINKO #25

---

QAS: R. GROVER #22  
CAMP VERDE MO

---

Operator: R. GROVER #22  
CAMP VERDE MO

---

Subject: SQAP,TEST  
DOB: 08/19/1971  
Sex: M Weight: 225

---

15 Minute Deprivation Period? Yes

---

Test	g/210L	Time
Air Blank	0.000	14:09:02
Diagnostic Test	Pass	14:09:29
Air Blank	0.000	14:09:56
0.100 Cal Check	0.096	14:10:15
Air Blank	0.000	14:10:44
Subject Test	0.076*	14:13:53
Air Blank	0.000	14:14:24
Air Blank	0.000	14:14:51
Subject Test	XXX**	14:15:11
Air Blank	0.000	14:15:45

\*Deficient Sample

\*\*Invalid Sample

Begin new deprivation period

**Not a Successfully Completed Test Sequence**

*Russell M. #22*

