

90 DAY

EXHIBIT G-2  
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY  
STANDARD QUALITY ASSURANCE PROCEDURES  
INTOXILYZER MODEL 8000  
STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST Russell A. Gower AGENCY C.V.M.O.  
DATE 8-25-10 TIME 1439 "  
INTOXILYZER SERIAL # 80-001096 DPS LOCATION CVMO

- (✓) 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution .100 AC.  
OR  
Pour a standard alcohol concentration solution \_\_\_\_\_ AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach  $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$
- (✓) 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
- (✓) 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
- (✓) 4. Air blank completed.
- (✓) 5. Calibration check completed. Test results 0. 098 AC.
- (✓) 6. Air blank completed.
- (✓) 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE R. Gower

DPS Form Exh G-2 (Rev 05-01)

Lot # F329016  
Exp DATE 12-4-10

EXHIBIT G-4  
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY  
STANDARD QUALITY ASSURANCE PROCEDURES  
INTOXILYZER MODEL 8000  
STANDARD QUALITY ASSURANCE PROCEDURE

QA SPECIALIST RUSSELL GROVER AGENCY CVMO  
DATE 8-25-10 TIME 1504  
INTOXILYZER SERIAL # 80-001096 OPS LOCATION CVMO

() 1. Display Reads "PUSH BUTTON TO START"

DIAGNOSTIC TESTS

- () 1. Clock time check.
- () 2. Date check.

OPERATIONAL TESTS

- () 1. Alcohol-free subject test result 0. 000 AC.
- () 2. Error recognition logic system functioning.  
Not a Successfully Completed Test Sequence printed
- () 3. Proper sample recognition system.  
Not a Successfully Completed Test Sequence printed  
Deficient sample printed.
- () 4. Standard Calibration Check standard 0. 100 AC. Result 0. 098 AC.

Instrument is operating properly and accurately. YES  NO

COMMENTS BOTH INTOX HAVE 90 DAY TEST DONE

SIGNATURE R. W #22

INTOXILYZER 8000

Location: CVMO

Serial Number: 80-001096

Core Version: 8105.48

08/25/2010

14:39:05

-----  
Standard Lot#: 832901G

Last Changed By: R. GROVER #22  
-----

QAS: R. GROVER #22

CAMP VERDE MO  
-----

Operator: R. GROVER #22

CAMP VERDE MO  
-----

Test	g/210L	Time
Air Blank	0.000	14:39:38
0.100 Cal Check	0.098	14:39:59
Air Blank	0.000	14:40:30

Preventative Maintenance  
Performed -- Timer Reset

*R. Grover #22*

*8/25/10*

INTOXILYZER 8000  
Location: CVMO  
Serial Number: 80-001096  
Core Version: 8105.48  
08/25/2010

14:44:00

---

Standard Lot#: 832901G  
Last Changed By: R. GROVER #22

---

QAS: R. GROVER #22  
CAMP VERDE MO

---

Operator: R. GROVER #22  
CAMP VERDE MO

---

Subject: SQAP,TEST  
DOB: 08/19/1971  
Sex: M Weight: 225

---

15 Minute Deprivation Period? Yes

---

Test	g/210L	Time
Air Blank	0.000	14:45:04
Diagnostic Test	Pass	14:45:39
Air Blank	0.000	14:46:06
0.100 Cal Check	0.098	14:46:27
Air Blank	0.000	14:46:59
Subject Test	0.000	14:47:28
Air Blank	0.000	14:47:55
Wait	ABT*	
Air Blank	0.000	14:48:33

\*Sequence Aborted

**Not a Successfully Completed Test Sequence**

*R. W. # 8/25/10*

INTOXILYZER 8000  
Location: CVMO  
Serial Number: 80-001096  
Core Version: 8105.48  
08/25/2010

14:49:35

---

Standard Lot#: 832901G  
Last Changed By: R. GROVER #22

---

QAS: R. GROVER #22  
CAMP VERDE MO

---

Operator: R. GROVER #22  
CAMP VERDE MO

---

Subject: SQAP,TEST  
DOB: 08/19/1971  
Sex: M Weight: 225

---

15 Minute Deprivation Period? Yes

---

Test	g/210L	Time
Air Blank	IPS*	14:50:12
Air Blank	0.000	14:50:41

\*Improper Sample

**Not a Successfully Completed Test Sequence**

*R. Grover #22*

INTOXILYZER 8000

Location: CVMO

Serial Number: 80-001096

Core Version: 8105.48

08/25/2010

14:51:51

---

Standard Lot#: 832901G

Last Changed By: R. GROVER #22

---

QAS: R. GROVER #22  
CAMP VERDE MO

---

Operator: R. GROVER #22  
CAMP VERDE MO

---

Subject: SQAP,TEST

DOB: 08/19/1971

Sex: M Weight: 225

---

15 Minute Deprivation Period? Yes

---

Test	g/210L	Time
Air Blank	RFI*	14:52:25
Air Blank	0.000	14:52:53

\*RFI Detect

**Not a Successfully Completed Test Sequence**

*R. Gro #22*

*8/25/10*

INTOXILYZER 8000  
Location: CVMO  
Serial Number: 80-001096  
Core Version: 8105.48  
08/25/2010

14:53:54

---

Standard Lot#: 832901G  
Last Changed By: R. GROVER #22

---

QAS: R. GROVER #22  
CAMP VERDE MO

---

Operator: R. GROVER #22  
CAMP VERDE MO

---

Subject: SQAP,TEST  
DOB: 08/19/1971  
Sex: M Weight: 225

---

15 Minute Deprivation Period? Yes

---

Test	g/210L	Time
Air Blank	0.000	14:54:41
Diagnostic Test	Pass	14:55:16
Air Blank	0.000	14:55:43
0.100 Cal Check	0.098	14:56:04
Air Blank	0.000	14:56:36
Subject Test	0.090*	14:59:43
Air Blank	0.000	15:00:14
Air Blank	0.000	15:00:42
Subject Test	XXX**	15:01:01
Air Blank	0.000	15:01:33

\*Deficient Sample

\*\*Invalid Sample

Begin new deprivation period

**Not a Successfully Completed Test Sequence**

*R. Grover #22*

*8/25/10*

INTOXILYZER 8000

Location: CVMO

Serial Number: 80-001096

Core Version: 8105.48

08/25/2010

15:03:37

-----  
Standard Lot#: 832901G

Last Changed By: R. GROVER #22  
-----

QAS: R. GROVER #22

CAMP VERDE MO  
-----

Operator: R. GROVER #22

CAMP VERDE MO  
-----

Test	g/210L	Time
Air Blank	0.000	15:04:07
0.100 Cal Check	0.098	15:04:28
Air Blank	0.000	15:05:00

Preventative Maintenance  
Performed -- Timer Reset

*R. Grover #22*

*8/25/10*