


EXHIBIT G-2
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY
STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000
STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST RUSSELL GROVER AGENCY CVMD
DATE 9-27-10 TIME 0942
INTOXILYZER SERIAL # 80-0010960PS LOCATION CVMD

- (✓) 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution .100 AC.
OR
Pour a standard alcohol concentration solution _____ AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$
- (✓) 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
- (✓) 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
- (✓) 4. Air blank completed.
- (✓) 5. Calibration check completed. Test results 0. 098 AC.
- (✓) 6. Air blank completed.
- (✓) 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE 

DPS Form Exh G-2 (Rev 05-01)

EXP DATE 12-4-10
LOT # 8329016

31 DAY

INTOXILYZER 8000

Location: CVMO

Serial Number: 80-001096

Core Version: 8105.48

09/27/2010

09:41:07

Standard Lot#: 767153

Last Changed By: R. GROVER #22

QAS: R. GROVER #22

CAMP VERDE MO

Operator: R. GROVER #22

CAMP VERDE MO

Test	g/210L	Time
Air Blank	0.000	09:41:56
0.100 Cal Check	0.098	09:42:17
Air Blank	0.000	09:42:49

Preventative Maintenance
Performed -- Timer Reset