

31 DAY Both

EXHIBIT G-2

THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY

STANDARD QUALITY ASSURANCE PROCEDURES  
INTOXILYZER MODEL 8000.

STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST RUSSELL GROVER AGENCY CVMO  
DATE 12-21-10 TIME 2335  
INTOXILYZER SERIAL # 80-001096 LOCATION CVMO

- () 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution .100 AC.  
OR  
Pour a standard alcohol concentration solution \_\_\_\_\_ AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach  $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$
- () 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
- () 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
- () 4. Air blank completed.
- () 5. Calibration check completed. Test results 0. 098 AC.
- () 6. Air blank completed.
- () 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE R. Grover #22

INTOXILYZER 8000  
Location: CVMO  
Serial Number: 80-001096  
Core Version: 8105.48  
12/21/2010 23:35:02

-----  
Standard Lot#: 767153  
Last Changed By: J. SLOWINSKI #6101  
-----

QAS: R. GROVER #22  
CAMP VERDE MO  
-----

Operator: R. GROVER #22  
CAMP VERDE MO  
-----

Test	g/210L	Time
Air Blank	0.000	23:35:33
0.100 Cal Check	0.098	23:35:54
Air Blank	0.000	23:36:25

Preventative Maintenance  
Performed -- Timer Reset

EXHIBIT G-2  
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY  
STANDARD QUALITY ASSURANCE PROCEDURES  
INTOXILYZER MODEL 8000  
STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST RUSSELL GROVER AGENCY CVMO

DATE 12-21-10 TIME 2326

INTOXILYZER SERIAL # 80-001736 DPS LOCATION CVMO

- () 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution 100 AC.  
OR  
Pour a standard alcohol concentration solution \_\_\_\_\_ AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach  $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$
- () 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
- () 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
- () 4. Air blank completed.
- () 5. Calibration check completed. Test results 0. 098 AC.
- () 6. Air blank completed.
- () 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE R. Grover #22

INTOXILYZER 8000  
Location: DPS  
Serial Number: 80-001736  
Core Version: 8105.48  
12/21/2010 23:25:19

-----  
Standard Lot#: 767153  
Last Changed By: R. GROVER #22

-----  
QAS: R. GROVER #22  
CAMP VERDE MO

-----  
Operator: R. GROVER #22  
CAMP VERDE MO

-----  
Test g/210L Time  
-----  
Air Blank 0.000 23:25:50  
0.100 Cal Check 0.098 23:26:11  
Air Blank 0.000 23:26:40

Preventative Maintenance  
Performed -- Timer Reset