

EXHIBIT G-2  
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY  
STANDARD QUALITY ASSURANCE PROCEDURES  
INTOXILYZER MODEL 8000

STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST RUSSELL WOVER AGENCY CUMM  
DATE 7-20-11 TIME 0848  
INTOXILYZER SERIAL # 80-001096 DPS LOCATION CUMM

- (✓) 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution 100 AC  
OR  
Pour a standard alcohol concentration solution \_\_\_\_\_ AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach  $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$
- (✓) 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
- (✓) 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
- (✓) 4. Air blank completed
- (✓) 5. Calibration check completed. Test results 0.096 AC
- (✓) 6. Air blank completed.
- ( ) 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE R. Wover D-12

DPS Form Exh G-2 (Rev 05-01)

Lot # 767153

EXP DATE 5-21-12

INTOXILYZER 8000  
Location: CVMO  
Serial Number: 80-001096  
Core Version: 8105.48  
07/20/2011 08:47:05

-----  
Standard Lot#: 767153  
Last Changed By: J. SLOWINSKI #6101  
-----

QAS: R. GROVER #22  
CAMP VERDE MO  
-----

Operator: R. GROVER #22  
CAMP VERDE MO  
-----

Test	g/210L	Time
Air Blank	0.000	08:47:50
0.100 Cal Check	0.098	08:48:11
Air Blank	0.000	08:48:42

Preventative Maintenance  
Performed -- Timer Reset

Lot # 767153  
EXP DATE 5-21-12