

31 DAY TEST

EXHIBIT G-2
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY

STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000.

STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST RUSSELL COOPER AGENCY CVMO

DATE 1-17-12 TIME _____

INTOXILYZER SERIAL # 80-001096DPS LOCATION CVMO

- () 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution .100 AC
OR
Pour a standard alcohol concentration solution _____ AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach 34°C ± 0.2°C
- () 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
- () 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
- () 4. Air blank completed.
- () 5. Calibration check completed. Test results 0. 099 AC.
- () 6. Air blank completed.
- () 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE R. Cooper D-12

DPS Form Exh G-2 (Rev 05-01)

Lot # 767153
EXN 5/21/12
DATE

INTOXILYZER 8000
Location: CVMO
Serial Number: 80-001096
Core Version: 8105.48
01/17/2012 08:12:32

Standard Lot#: 767153
Last Changed By: J. SLOWINSKI #6101

QAS: R. GROVER #22
CAMP VERDE MO

Operator: R. GROVER #22
CAMP VERDE MO

Test	g/210L	Time
Air Blank	0.000	08:13:31
0.100 Cal Check	0.099	08:13:52
Air Blank	0.000	08:14:23

Preventative Maintenance
Performed -- Timer Reset

R. Grover 1-17-12