

31 DAY BIRTH

EXHIBIT G-2

THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY

STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000

STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST RUSSELL GROVER AGENCY CVMO

DATE 3-20-12 TIME 1353

INTOXILYZER SERIAL # 80-001096DR LOCATION CVMO

- (✓) 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution 5100 AC.
OR
Pour a standard alcohol concentration solution _____ AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$
- (✓) 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
- (✓) 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
- (✓) 4. Air blank completed.
- (✓) 5. Calibration check completed. Test results 0. 0.98 AC.
- (✓) 6. Air blank completed.
- (✓) 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE R. Grover 3-20-12

INTOXILYZER 8000
Location: CVMO
Serial Number: 80-001096
Core Version: 8105.48
03/20/2012 13:52:21

Standard Lot#: 767153
Last Changed By: J. SLOWINSKI #6101

QAS: R. GROVER #22
CAMP VERDE MO

Operator: R. GROVER #22
CAMP VERDE MO

Test	g/210L	Time
Air Blank	0.000	13:53:34
0.100 Cal Check	0.098	13:53:55
Air Blank	0.000	13:54:26

Preventative Maintenance
Performed -- Timer Reset

R. Grover D-12

3/20/12