

EXHIBIT G-2
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY
STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000
STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST Deborah Hughes AGENCY Camp Verde Marshals Ofc
DATE 2/11/15 TIME 0849
INTOXILYZER SERIAL # 80-601616 LOCATION CWMO - DUI Trailer

- (✓) 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution 0.100 AC.
OR
Pour a standard alcohol concentration solution _____ AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$
- (✓) 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
- (✓) 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
- (✓) 4. Air blank completed.
- (✓) 5. Calibration check completed. Test results 0.099 AC.
- (✓) 6. Air blank completed.
- (✓) 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE [Signature] #53

INTOXILYZER 8000
Location: CAMP VERDE MO
Serial Number: 80-001096
Core Version: 8105.48
02/11/2015 08:49:53

Standard Lot#: 14014100A1
Last Changed By: D. JACOBS #CV23

QAS: D. HUGHES #53
CAMP VERDE MO

Operator: D. HUGHES #53
CAMP VERDE MO

Test	g/210L	Time
Air Blank	0.000	08:50:24
0.100 Cal Check	0.099	08:50:45
Air Blank	0.000	08:51:15

Preventative Maintenance
Performed -- Timer Reset

EXHIBIT G-4
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY
STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000

STANDARD QUALITY ASSURANCE PROCEDURE

QA SPECIALIST Don Jacobs AGENCY CVMO
DATE 2/11/15 TIME 0853
INTOXILYZER SERIAL # 80-001096 LOCATION DUI Trailer

1. Display Reads "PUSH BUTTON TO START"

DIAGNOSTIC TESTS

1. Clock time check.
 2. Date check.

OPERATIONAL TESTS

1. Alcohol-free subject test result 0. 000 AC.
 2. Error recognition logic system functioning.
Not a Successfully Completed Test Sequence printed
 3. Proper sample recognition system.
Not a Successfully Completed Test Sequence printed
Deficient sample printed.
 4. Standard Calibration Check standard 0. 100 AC. Result 0. 099 AC.

Instrument is operating properly and accurately. YES NO

COMMENTS

SIGNATURE  0123

INTOXILYZER 8000

Location: CAMP VERDE MO

Serial Number: 80-001096

Core Version: 8105.48

02/11/2015

09:01:32

Standard Lot#: 14014100A1

Last Changed By: D. JACOBS

#CV23

QAS: D. HUGHES

#53

CAMP VERDE MO

Operator: D. JACOBS

#CV23

CAMP VERDE MO

Subject: TEST,TEST

DOB: 02/02/1986

Sex: M

Weight: 180

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	0.000	09:02:37
Diagnostic Test	Pass	09:03:12
Air Blank	0.000	09:03:39
0.100 Cal Check	0.099	09:04:00
Air Blank	0.000	09:04:30
Subject Test	RGE*	09:04:52
Air Blank	0.000	09:05:34

*Range Exceeded

Not a Successfully Completed Test Sequence

INTOXILYZER 8000
Location: CAMP VERDE MO
Serial Number: 80-001096
Core Version: 8105.48
02/11/2015

08:55:19

Standard Lot#: 14014100A1
Last Changed By: D. JACOBS #CV23

QAS: D. HUGHES #53
CAMP VERDE MO

Operator: D. JACOBS #CV23
CAMP VERDE MO

Subject: TEST,TEST
DOB: 02/02/1986
Sex: M Weight: 180

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	0.000	08:56:32
Diagnostic Test	Pass	08:57:07
Air Blank	0.000	08:57:35
0.100 Cal Check	0.099	08:57:56
Air Blank	0.000	08:58:25
Subject Test	0.000	08:58:45
Air Blank	0.000	08:59:12
Wait	IPS*	
Air Blank	0.000	09:00:14

*Improper Sample

Not a Successfully Completed Test Sequence

INTOXILYZER 8000
Location: CAMP VERDE MO
Serial Number: 80-001096
Core Version: 8105.48
02/11/2015

09:11:50

Standard Lot#: 14014100A1
Last Changed By: D. JACOBS #CV23

QAS: D. HUGHES #53
CAMP VERDE MO

Operator: D. JACOBS #CV23
CAMP VERDE MO

Subject: TEST,TEST
DOB: 02/02/1986
Sex: M Weight: 180

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	0.000	09:12:52
Diagnostic Test	Pass	09:13:27
Air Blank	0.000	09:13:54
0.100 Cal Check	0.099	09:14:15
Air Blank	0.000	09:14:45
Subject Test	0.129*	09:17:52
Air Blank	0.000	09:18:24
Air Blank	0.000	09:18:52
Subject Test	ABT**	09:19:09
Air Blank	RFI***	09:19:17

*Deficient Sample
**Sequence Aborted
***RFI Detect

Not a Successfully Completed Test Sequence

EXHIBIT G-2
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY
STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000

STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST Dan Sachs AGENCY CVMO
DATE 2/11/15 TIME 0921
INTOXILYZER SERIAL # 80-001096 LOCATION DUI Trailer

- () 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution 0.100 AC.
OR
Pour a standard alcohol concentration solution _____ AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$
- (X) 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
(X) 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
(X) 4. Air blank completed.
(X) 5. Calibration check completed. Test results 0. 099 AC.
(X) 6. Air blank completed.
(X) 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE 

INTOXILYZER 8000

Location: CAMP VERDE MO

Serial Number: 80-001096

Core Version: 8105.48

02/11/2015

09:21:19

Standard Lot#: 14014100A1

Last Changed By: D. JACOBS #CV23

QAS: D. JACOBS #CV23

CAMP VERDE MO

Operator: D. JACOBS #CV23

CAMP VERDE MO

Test	g/210L	Time
Air Blank	0.000	09:21:50
0.100 Cal Check	0.099	09:22:11
Air Blank	0.000	09:22:40

Preventative Maintenance
Performed -- Timer Reset