

EXHIBIT G-2  
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY  
STANDARD QUALITY ASSURANCE PROCEDURES  
INTOXILYZER MODEL 8000

STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST Deborah Hughes AGENCY Camp Verde Marshal's Ofc

DATE 12/21/18 TIME 1143

INTOXILYZER SERIAL # 80-001096 LOCATION CVMO - Booking Ofc

(X) 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution 0.100 AC.  
OR

Pour a standard alcohol concentration solution \_\_\_\_\_ AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach  $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$

- (X) 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
- (X) 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
- (X) 4. Air blank completed.
- (X) 5. Calibration check completed. Test results 0.100 AC.
- (X) 6. Air blank completed.
- (X) 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE [Signature] #63

DPS Form Exh G-2 (Rev 05-01)

Before Standard Change

INTOXILYZER 8000

Location: CAMP VERDE MO

Serial Number: 80-001096

Core Version: 8105.48

12/21/2018

11:45:22

-----  
Standard Lot#: 32316100A1

Last Changed By: D. JACOBS #CV08

-----  
QAS: D. HUGHES #63

CAMP VERDE MO

-----  
Operator: D. HUGHES #63

CAMP VERDE MO

-----  
Test                      g/210L                      Time  
-----  
Air Blank                      0.000                      11:46:09  
0.100 Cal Check                0.100                      11:46:30  
Air Blank                      0.000                      11:47:00

Preventative Maintenance  
Performed -- Timer Reset

EXHIBIT G-2  
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY  
STANDARD QUALITY ASSURANCE PROCEDURES  
INTOXILYZER MODEL 8000

STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST Deborah Hughes AGENCY Camp Verde Marshal's ofc  
DATE 12-21-18 TIME 1150  
INTOXILYZER SERIAL # 80-001094 LOCATION CUMO - Booking

1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution 0.100 AC.  
OR  
Pour a standard alcohol concentration solution \_\_\_\_\_ AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach  $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$
2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
4. Air blank completed.
5. Calibration check completed. Test results 0. 099 AC.
6. Air blank completed.
7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE [Signature] #63

DPS Form Exh G-2 (Rev 05-01)

After Standard change

INTOXILYZER 8000  
Location: CAMP VERDE MO  
Serial Number: 80-001096  
Core Version: 8105.48  
12/21/2018 11:50:01

-----  
Standard Lot#: 13518100A1  
Last Changed By: D. HUGHES #63  
-----

QAS: D. HUGHES #63  
CAMP VERDE MO  
-----

Operator: D. HUGHES #63  
CAMP VERDE MO  
-----

Test	g/210L	Time
Air Blank	0.000	11:50:37
0.100 Cal Check	0.099	11:50:57
Air Blank	0.000	11:51:27

Preventative Maintenance  
Performed -- Timer Reset

EXHIBIT G-4  
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY  
STANDARD QUALITY ASSURANCE PROCEDURES  
INTOXILYZER MODEL 8000

STANDARD QUALITY ASSURANCE PROCEDURE

QA SPECIALIST Deborah Hughes AGENCY Camp Verde Marshal's Ofc  
DATE 12-21-18 TIME 1154  
INTOXILYZER SERIAL # 80-001094 LOCATION CWMO-Booking Ofc

1. Display Reads "PUSH BUTTON TO START"

DIAGNOSTIC TESTS

- 1. Clock time check.
- 2. Date check.


OPERATIONAL TESTS

- 1. Alcohol-free subject test result 0. 000 AC.
- 2. Error recognition logic system functioning.  
Not a Successfully Completed Test Sequence printed
- 3. Proper sample recognition system.  
Not a Successfully Completed Test Sequence printed  
Deficient sample printed.
- 4. Standard Calibration Check standard 0. 100 AC. Result 0. 098 AC.

Instrument is operating properly and accurately. YES  NO

COMMENTS

\_\_\_\_\_  
\_\_\_\_\_  
\_\_\_\_\_

SIGNATURE  #63

INTOXILYZER 8000  
Location: CAMP VERDE MO  
Serial Number: 80-001096  
Core Version: 8105.48  
12/21/2018 11:56:45

---

Standard Lot#: 13518100A1  
Last Changed By: D. HUGHES #63

---

QAS: D. HUGHES #63  
CAMP VERDE MO

---

Operator: D. HUGHES #63  
CAMP VERDE MO

---

Subject: TEST,TEST  
DOB: 12/21/1963  
Sex: F Weight: 190

---

15 Minute Deprivation Period? Yes

---

Test	g/210L	Time
Air Blank	0.000	11:58:00
Diagnostic Test	Pass	11:58:35
Air Blank	0.000	11:59:02
0.100 Cal Check	0.099	11:59:23
Air Blank	0.000	11:59:53
Subject Test	0.000	12:00:14
Air Blank	ABT*	12:00:22
Air Blank	IPS**	12:00:31

\*Sequence Aborted

\*\*Improper Sample

**Not a Successfully Completed Test Sequence**

INTOXILYZER 8000  
Location: CAMP VERDE MO  
Serial Number: 80-001096  
Core Version: 8105.48  
12/21/2018

12:00:57

---

Standard Lot#: 13518100A1  
Last Changed By: D. HUGHES #63

---

QAS: D. HUGHES #63  
CAMP VERDE MO

---

Operator: D. HUGHES #63  
CAMP VERDE MO

---

Subject: TEST,TEST  
DOB: 12/21/1963  
Sex: F Weight: 190

---

15 Minute Deprivation Period? Yes

---

Test	g/210L	Time
Air Blank	0.000	12:02:03
Diagnostic Test	Pass	12:02:38
Air Blank	0.000	12:03:06
0.100 Cal Check	0.099	12:03:26
Air Blank	0.000	12:03:56
Subject Test	0.107*	12:07:04
Air Blank	0.000	12:07:35
Air Blank	0.000	12:08:02
Subject Test	XXX**	12:08:42
Air Blank	RFI***	12:08:50

\*Deficient Sample

\*\*Invalid Sample

Begin new deprivation period

\*\*\*RFI Detect

**Not a Successfully Completed Test Sequence**

EXHIBIT G-2  
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY  
STANDARD QUALITY ASSURANCE PROCEDURES  
INTOXILYZER MODEL 8000

STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST Deborah Hughes AGENCY Camp Verde Marshals ofc  
DATE 12-21-18 TIME 1210  
INTOXILYZER SERIAL # 80-00196 LOCATION CV MO Booking ofc

- (X) 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution 0.100 AC.  
OR  
Pour a standard alcohol concentration solution \_\_\_\_\_ AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach  $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$
- (X) 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
- (X) 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
- (X) 4. Air blank completed.
- (X) 5. Calibration check completed. Test results 0. 098 AC.
- (X) 6. Air blank completed.
- (X) 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE [Signature] #63

DPS Form Exh G-2 (Rev 05-01)

Ambient Fail 1st Attempt  
2nd Attempt Good

INTOXILYZER 8000

Location: CAMP VERDE MO

Serial Number: 80-001096

Core Version: 8105.48

12/21/2018

12:10:12

-----  
Standard Lot#: 13518100A1

Last Changed By: D. HUGHES #63  
-----

QAS: D. HUGHES #63

CAMP VERDE MO  
-----

Operator: D. HUGHES #63

CAMP VERDE MO  
-----

Test	g/210L	Time
Air Blank	AMB*	12:10:47
Air Blank	0.000	12:11:15

\*Ambient Fail

