

EXHIBIT G-2  
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY  
STANDARD QUALITY ASSURANCE PROCEDURES  
INTOXILYZER MODEL 8000

STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST Daniel McGinnis AGENCY YCSO

DATE 08-22-09 TIME 2346

INTOXILYZER SERIAL # 80-001222 LOCATION Verde Jail

1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution .10 AC.  
OR

~~Pour a standard alcohol concentration solution \_\_\_\_\_ AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach  $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$~~

2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"  
3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.  
4. Air blank completed.  
5. Calibration check completed. Test results 0. .099 AC.  
6. Air blank completed.  
7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE Daniel McGinnis 543/QAS# 33426

INTOXILYZER 8000

Location: YCSO-VERDE JAIL

Serial Number: 80-001222

08/22/2009

23:46:54


-----  
Standard Lot#: 901401I

Last Changed By: D. MCGINNIS #8086  
-----

QAS: D. MCGINNIS #8086

YAVAPAI CSO  
-----

| Test      | g/210L | Time     |
|-----------|--------|----------|
| Air Blank | 0.000  | 23:47:22 |
| Cal Check | 0.099  | 23:47:40 |
| Air Blank | 0.000  | 23:48:10 |

  
Operator's Signature