

EXHIBIT G-2
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY
STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000

STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST Daniel McGinnis AGENCY YCSO

DATE 10-17-09 TIME 0030

INTOXILYZER SERIAL # 80-001222 LOCATION Verde Jail

DM 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution .10 AC.
OR

~~Pour a standard alcohol concentration solution _____ AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach 34°C ± 0.2°C~~

- DM 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
DM 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
DM 4. Air blank completed.
DM 5. Calibration check completed. Test results 0. 10 AC.
DM 6. Air blank completed.
DM 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE

Daniel McGinnis S43/QAS# 33426

INTOXILYZER 8000
Location: YCSO-VERDE JAIL
Serial Number: 80-001222
10/17/2009 00:30:07

Standard Lot#: 901401I
Last Changed By: D. MCGINNIS #8086

QAS: D. MCGINNIS #8086
YAVAPAI CSO

Test	g/210L	Time
Air Blank	0.000	00:30:34
Cal Check	0.100	00:30:53
Air Blank	0.000	00:31:23


Operator's Signature