

EXHIBIT G-2  
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY  
STANDARD QUALITY ASSURANCE PROCEDURES  
INTOXILYZER MODEL 8000  
STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST D. McGILLIS AGENCY YCSO

DATE 9-17-10 TIME 2130

INTOXILYZER SERIAL # 80-001222 LOCATION Verde

- CM 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution .16 AC.  
OR  
Pour a standard alcohol concentration solution \_\_\_\_\_ AC into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach  $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$
- CM 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
- CM 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
- CM 4. Air blank completed.
- CM 5. Calibration check completed. Test results 0. 10 AC.
- CM 6. Air blank completed.
- CM 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE [Signature] S43/QAS # 33426

INTOXILYZER 8000  
Location: YCSO-VERDE JAIL  
Serial Number: 80-001222  
Core Version: 8105.48  
09/17/2010

21:30:35

-----  
Standard Lot#: 712322  
Last Changed By: D. MCGINNIS #8086  
-----

QAS: D. MCGINNIS #8086  
YAVAPAI CSO  
-----

Operator: D. MCGINNIS #8086  
YAVAPAI CSO  
-----

Test	g/210L	Time
Air Blank	0.000	21:31:08
0.100 Cal Check	0.100	21:31:29
Air Blank	0.000	21:31:59

Preventative Maintenance  
Performed -- Timer Reset