

EXHIBIT G-2  
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY  
STANDARD QUALITY ASSURANCE PROCEDURES  
INTOXILYZER MODEL 8000  
STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST LOGAN MOODY AGENCY YAVAPAI CO. S.O.

DATE 01-13-2012 TIME 2123

INTOXILYZER SERIAL # 80-001222 LOCATION CAMP VERDE JAIL

- (  ) 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution .100 AC.  
OR  
Pour a standard alcohol concentration solution \_\_\_\_\_ AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach  $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$
- (  ) 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
- (  ) 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
- (  ) 4. Air blank completed.
- (  ) 5. Calibration check completed. Test results 0. 099 AC.
- (  ) 6. Air blank completed.
- (  ) 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE Logan Moody

INTOXILYZER 8000  
Location: YCSO-VERDE JAIL  
Serial Number: 80-001222  
Core Version: 8105.48  
01/13/2012 21:23:28

-----  
Standard Lot#: 15811100A1  
Last Changed By: L. MOODY #7486  
-----

QAS: L. MOODY #7486  
YAVAPAI CSO

-----  
Operator: L. MOODY #7486  
YAVAPAI CSO  
-----

Test	g/210L	Time
Air Blank	0.000	21:23:57
0.100 Cal Check	0.099	21:24:18
Air Blank	0.000	21:24:48

Preventative Maintenance  
Performed -- Timer Reset

EXHIBIT G-4  
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY

STANDARD QUALITY ASSURANCE PROCEDURES  
INTOXILYZER MODEL 8000

STANDARD QUALITY ASSURANCE PROCEDURE

QA SPECIALIST LOGAN MOODY AGENCY YAVAPAI CO. S.O.  
DATE 01-13-2012 TIME 2125  
INTOXILYZER SERIAL # 80-001222 LOCATION CAMP VERDE JAIL

() 1. Display Reads "PUSH BUTTON TO START"

DIAGNOSTIC TESTS

- () 1. Clock time check.  
() 2. Date check.

OPERATIONAL TESTS

- () 1. Alcohol-free subject test result 0. 000 AC.  
() 2. Error recognition logic system functioning.  
Not a Successfully Completed Test Sequence printed  
() 3. Proper sample recognition system.  
Not a Successfully Completed Test Sequence printed  
Deficient sample printed.  
() 4. Standard Calibration Check standard 0. 100 AC. Result 0. 098 AC.

Instrument is operating properly and accurately. YES  NO

COMMENTS

TAUK @ 800 psi

SIGNATURE Logan Moody

INTOXILYZER 8000  
Location: YCSO-VERDE JAIL  
Serial Number: 80-001222  
Core Version: 8105.48  
01/13/2012

21:25:11

---

Standard Lot#: 15811100A1  
Last Changed By: L. MOODY #7486

---

QAS: L. MOODY #7486  
YAVAPAI CSO

---

Operator: L. MOODY #7486  
YAVAPAI CSO

---

Subject: SQAP,TEST  
DOB: 01/01/2000  
Sex: M Weight: 190

---

15 Minute Deprivation Period? Yes

---

Test	g/210L	Time
Air Blank	0.000	21:25:59
Diagnostic Test	Pass	21:26:33
Air Blank	0.000	21:27:01
0.100 Cal Check	0.100	21:27:21
Air Blank	0.000	21:27:51
Subject Test	0.000	21:28:10
Air Blank	0.000	21:28:37
Wait	ABT*	
Air Blank	RFI**	21:31:01

\*Sequence Aborted  
\*\*RFI Detect

**Not a Successfully Completed Test Sequence**

INTOXILYZER 8000  
Location: YCSO-VERDE JAIL  
Serial Number: 80-001222  
Core Version: 8105.48  
01/13/2012

21:31:07

---

Standard Lot#: 15811100A1  
Last Changed By: L. MOODY #7486

---

QAS: L. MOODY #7486  
YAVAPAI CSO

---

Operator: L. MOODY #7486  
YAVAPAI CSO

---

Subject: SQAP,TEST  
DOB: 01/01/2000  
Sex: M Weight: 190

---

15 Minute Deprivation Period? Yes

---

Test	g/210L	Time
Air Blank	0.000	21:31:58
Diagnostic Test	Pass	21:32:32
Air Blank	0.000	21:32:59
0.100 Cal Check	0.099	21:33:20
Air Blank	0.000	21:33:50
Subject Test	0.062*	21:36:57
Air Blank	0.000	21:37:29
Air Blank	0.000	21:37:56
Subject Test	XXX**	21:39:00
Air Blank	IPS***	21:39:08

\*Deficient Sample

\*\*Invalid Sample

Begin new deprivation period

\*\*\*Improper Sample

**Not a Successfully Completed Test Sequence**

EXHIBIT G-2  
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY  
STANDARD QUALITY ASSURANCE PROCEDURES  
INTOXILYZER MODEL 8000  
STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST LOGAN MOODY AGENCY YAVAPAI CO. S.O.

DATE 01-13-2012 TIME 2139

INTOXILYZER SERIAL # 80-001222 LOCATION CAMP VERDE JAIL

- (✓) 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution .100 AC.  
OR  
Pour a standard alcohol concentration solution \_\_\_\_\_ AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach  $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$
- (✓) 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
- (✓) 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
- (✓) 4. Air blank completed.
- (✓) 5. Calibration check completed. Test results 0.098 AC.
- (✓) 6. Air blank completed.
- (✓) 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE Logan Moody

INTOXILYZER 8000

Location: YCSO-VERDE JAIL

Serial Number: 80-001222

Core Version: 8105.48

01/13/2012

21:39:22

-----  
Standard Lot#: 15811100A1

Last Changed By: L. MOODY #7486  
-----

QAS: L. MOODY #7486

YAVAPAI CSO  
-----

Operator: L. MOODY #7486

YAVAPAI CSO  
-----

Test	g/210L	Time
Air Blank	0.000	21:39:51
0.100 Cal Check	0.098	21:40:12
Air Blank	0.000	21:40:41

Preventative Maintenance  
Performed -- Timer Reset