

EXHIBIT G-2
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)
ARIZONA DEPARTMENT OF PUBLIC SAFETY
STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000
STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST LOGAN MOODY 7486 AGENCY YAVAPAI CO. S.O.
DATE 4-2-2014 TIME 2115
INTOXILYZER SERIAL # 80-001222 LOCATION YCSO - CAMP VERDE JAIL

- (05) 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution, 100 AC.
OR
Pour a standard alcohol concentration solution AC into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach 34°C ± 0.2°C
- (31) 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
- (01) 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
- (04) 4. Air blank completed.
- (04) 5. Calibration check completed. Test results 0. 096 AC.
- (04) 6. Air blank completed.
- (01) 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE Logan Moody 7486

INTOXILYZER 8000

Location: YCSO-VERDE JAIL

Serial Number: 80-001222

Core Version: 8105.48

04/02/2014

21:15:02

Standard Lot#: 15013100A2

Last Changed By: R. POOLE #9014

QAS: L. MOODY #7486

YAVAPAI CSO

Operator: L. MOODY #7486

YAVAPAI CSO

Test	g/210L	Time
Air Blank	0.000	21:15:31
0.100 Cal Check	0.098	21:15:52
Air Blank	0.000	21:16:22

Preventative Maintenance
Performed -- Timer Reset

EXHIBIT G-4
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY
STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000
STANDARD QUALITY ASSURANCE PROCEDURE

QA SPECIALIST LOGAN MOORE 7486 AGENCY YAVAPAI CO. S.O.

DATE 4-2-2014 TIME 2116

INTOXILYZER SERIAL # 80-001222 LOCATION YCSO - CAMP VERDE STATE

(M) 1. Display Reads "PUSH BUTTON TO START"

DIAGNOSTIC TESTS

- (M) 1. Clock time check.
- (M) 2. Date check.

OPERATIONAL TESTS

- (M) 1. Alcohol-free subject test result 0. 000 AC.
- (M) 2. Error recognition logic system functioning.
Not a Successfully Completed Test Sequence printed
- (M) 3. Proper sample recognition system.
Not a Successfully Completed Test Sequence printed
Deficient sample printed.
- (M) 4. Standard Calibration Check standard 0. 100 AC. Result 0. 098 AC.

Instrument is operating properly and accurately. YES NO

COMMENTS

450 psi

SIGNATURE Logan Moore 7486

INTOXILYZER 8000

Location: YCSO-VERDE JAIL

Serial Number: 80-001222

Core Version: 8105.48

04/02/2014

21:16:33

Standard Lot#: 15013100A2

Last Changed By: R. POOLE

#9014

QAS: L. MOODY

#7486

YAVAPAI CSO

Operator: L. MOODY

#7486

YAVAPAI CSO

Subject: SQAP, TEST

DOB: 01/01/2000

Sex: M

Weight: 190

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	0.000	21:17:22
Diagnostic Test	Pass	21:17:56
Air Blank	0.000	21:18:24
0.100 Cal Check	0.099	21:18:45
Air Blank	0.000	21:19:14
Subject Test	0.000	21:19:30
Air Blank	ABT*	21:19:44
Air Blank	RFI**	21:19:55

*Sequence Aborted

**RFI Detect

Not a Successfully Completed Test Sequence

INTOXILYZER 8000

Location: YCSO-VERDE JAIL

Serial Number: 80-001222

Core Version: 8105.48

04/02/2014

21:20:08

Standard Lot#: 15013100A2

Last Changed By: R. POOLE

#9014

QAS: L. MOODY

#7486

YAVAPAI CSO

Operator: L. MOODY

#7486

YAVAPAI CSO

Subject: SQAP,TEST

DOB: 01/01/2000

Sex: M

Weight: 190

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	0.000	21:20:56
Diagnostic Test	Pass	21:21:29
Air Blank	0.000	21:21:57
0.100 Cal Check	0.099	21:22:18
Air Blank	0.000	21:22:47
Subject Test	0.068*	21:25:55
Air Blank	0.000	21:26:27
Air Blank	0.000	21:26:54
Subject Test	XXX**	21:27:13
Air Blank	IPS***	21:27:20

*Deficient Sample

**Invalid Sample

Begin new deprivation period

***Improper Sample

Not a Successfully Completed Test Sequence

EXHIBIT G-2
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)
ARIZONA DEPARTMENT OF PUBLIC SAFETY
STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000
STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST LOGAN MOODY AGENCY YAVAPAI CO. S.O.

DATE 4-2-14 TIME 2128

INTOXILYZER SERIAL # 80-601222 LOCATION YCSO - CAMP VERDE TRIL

- (67) 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution 100 AC.
OR
Pour a standard alcohol concentration solution _____ AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$
- (67) 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
- (67) 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
- (67) 4. Air blank completed.
- (67) 5. Calibration check completed. Test results 0.098 AC.
- (67) 6. Air blank completed.
- (67) 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE Logan Moody

DPS Form Exh G-2 (Rev 05-01)

AMBIENT FAIL IMMEDIATELY FOLLOWING SWAP.

INTOXILYZER 8000
Location: YCSO-VERDE JAIL
Serial Number: 80-001222
Core Version: 8105.48
04/02/2014

21:27:41

Standard Lot#: 15013100A2
Last Changed By: R. POOLE #9014

QAS: L. MOODY #7486
YAVAPAI CSO

Operator: L. MOODY #7486
YAVAPAI CSO

Test	g/210L	Time
Air Blank	AMB*	21:28:10
Air Blank	0.000	21:28:38

*Ambient Fail

INTOXILYZER 8000

Location: YCSQ-VERDE JAIL

Serial Number: 80-001222

Core Version: 8105.48

04/02/2014

21:28:44

Standard Lot#: 15013100A2

Last Changed By: R. POOLE #9014

QAS: L. MOODY #7486

YAVAPAI CSO

Operator: L. MOODY #7486

YAVAPAI CSO

Test	g/210L	Time
Air Blank	0.000	21:29:12
0.100 Cal Check	0.098	21:29:33
Air Blank	0.000	21:30:03

Preventative Maintenance
Performed -- Timer Reset