

EXHIBIT G-2  
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY  
STANDARD QUALITY ASSURANCE PROCEDURES  
INTOXILYZER MODEL 8000

STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST Chris Krogedal AGENCY YCSO  
DATE 03/08/08 TIME 2057  
INTOXILYZER SERIAL # 80-001455 LOCATION Bagdad

- (Ch) 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution .100 AC.  
OR  
Pour a standard alcohol concentration solution \_\_\_\_\_ AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach  $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$
- (Ch) 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
- (Ch) 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
- (Ch) 4. Air blank completed.
- (Ch) 5. Calibration check completed. Test results 0. 101 AC.
- (Ch) 6. Air blank completed.
- (Ch) 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE

Ch Krogedal #570/6535

INTOXILYZER 8000

Location: YAVAPAI SO, BAGDAD

Serial Number: 80-001455

03/08/2008

20:50:51


-----  
Standard Lot#: 730901I

Last Changed By: C. KROGEDAL #6535

-----  
QAS: C. KROGEDAL #6535

YAVAPAI CSO  
-----

Test	g/210L	Time
Air Blank	0.000	20:51:18
Cal Check	0.101	20:51:37
Air Blank	0.000	20:52:07

 #570/6535

Operator's Signature

QAS# 33200