

EXHIBIT G-2

THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY

STANDARD QUALITY ASSURANCE PROCEDURES  
INTOXILYZER MODEL 8000

STANDARD CALIBRATION CHECK PROCEDURE

QA SPEC # 562 AGENCY YCSO

DATE 11/10/05 TIME 0742

INTOXILYZER SERIAL # 80-001455 LOCATION Bagdad Substation

(✓) 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution 100 AC.  
OR

Pour a standard alcohol concentration solution \_\_\_\_\_ AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach  $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$

- (✓) 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
- (✓) 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
- (✓) 4. Air blank completed.
- (✓) 5. Calibration check completed. Test results 0. 100 AC.
- (✓) 6. Air blank completed.
- (✓) 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE Deann J. Ashby #562/QAS#33400



EXHIBIT G-4

THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY

STANDARD QUALITY ASSURANCE PROCEDURES  
INTOXILYZER MODEL 8000

STANDARD QUALITY ASSURANCE PROCEDURE

QA SPECIALIST D. Ashby #S162 AGENCY YCSO

DATE 6-24-09 TIME 0745

INTOXILYZER SERIAL # 80-001455 LOCATION Bagdad Substation

(  ) 1. Display Reads "PUSH BUTTON TO START"

DIAGNOSTIC TESTS

(  ) 1. Clock time check.

(  ) 2. Date check.

OPERATIONAL TESTS

(  ) 1. Alcohol-free subject test result 0. 000 AC.

(  ) 2. Error recognition logic system functioning.  
Not a Successfully Completed Test Sequence printed

(  ) 3. Proper sample recognition system.  
Not a Successfully Completed Test Sequence printed  
Deficient sample printed.

(  ) 4. Standard Calibration Check standard 0. 100 AC. Result 0. 099 AC.

Instrument is operating properly and accurately. YES  NO

COMMENTS

SIGNATURE Deann J Ashby #S162/QAS# 33400

1971-1972  
1973-1974  
1975-1976  
1977-1978

1979-1980  
1981-1982

1983-1984  
1985-1986

1987-1988  
1989-1990

1991-1992  
1993-1994  
1995-1996

1997-1998

1999-2000  
2001-2002  
2003-2004  
2005-2006  
2007-2008  
2009-2010  
2011-2012  
2013-2014  
2015-2016  
2017-2018

2019-2020

EMPOKILYER: 0003

Location: YAVAPAI CO, STODAF

Serial Number: 80-001455

06/24/2000

07/01/00

Standard Lot#: 20000

Cost Charged to: 00000

QAS: D. ASHBY

YAVAPAI CO

Operator: D. ASHBY

YAVAPAI CO

Subject: SOAP TEST

LOP: 01/01/1999

Sex: M

Weight: 150

15 Minute Deprivation Period: Yes

Test	Y/2101	Y/N
Air Blank	1904	07/01/00
Air Blank	0000	07/01/00

\*Improper Sample

Not a Successfully Completed Test Report

INSTRUMENT: 8000

Location: YAVAPAI SO, SACDAD

Serial Number : 86-001455

06/24/2009

07:45:11

Standard Lot#: 200011

Last Changed by: L. ASHBY

0.00

QAC: D. ASHEY

0.00

YAVAPAI CSO

Operator: D. ASHEY

0.00

YAVAPAI CSO

Subject: ROAD, ROAD

DOB: 01/01/1999

Sex: M

Weight: 150

15 Minute Deprivation Period: Yes

Test	g/210l	Time
Air Blank	0.00	0.00
Air Blank	0.000	0.00

\*RFI Detect

Vol e Susceptibility (unpublished) 0.00

CRUCIALSERV 0009

LOCAL TEL: YAVAPAI 001 846000

Serial Number: 000001175

06/27/2009

00000000

Sampled in: 00000000

Anal Channel: 00000000

Serial: 00000000

YAVAPAI 000

Operator: 00000000

YAVAPAI 000

Subject: 00000000

POB: 01/01/1998

Sex: M

Weight: 130

15 Minute Deprivation Period: 000

-----

Dist	Weight	Unit
Air Blank	0.000	0.000000
Composite Test	0.000	0.000000
Air Blank	0.000	0.000000
0.100 Chl Check	0.100	0.000000
Air Blank	0.000	0.000000
Subject Test	0.100	0.000000
Air Blank	0.000	0.000000
Air Blank	0.000	0.000000
Subject Test	0.100	0.000000
Air Blank	0.000	0.000000

\*Deficient Sample

\*Invalid Sample

Serial: 0000000000000000

Not a Successful or Deficient Sample

EXHIBIT G-2  
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY

STANDARD QUALITY ASSURANCE PROCEDURES  
INTOXILYZER MODEL 8000

STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST D. Ashby #562 AGENCY YCSO  
DATE 10-24-09 TIME 0800  
INTOXILYZER SERIAL # 80-001455 LOCATION Bagdad Substation

- (✓) 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution 100 AC.  
OR

Pour a standard alcohol concentration solution \_\_\_\_\_ AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach  $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$

- (✓) 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"  
(✓) 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.  
(✓) 4. Air blank completed.  
(✓) 5. Calibration check completed. Test results 0.099 AC.  
(✓) 6. Air blank completed.  
(✓) 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE Deann Ashby #562 / QAS# 33400

