

EXHIBIT G-2  
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY  
STANDARD QUALITY ASSURANCE PROCEDURES  
INTOXILYZER MODEL 8000  
STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST R. RIVERA AGENCY YAVAPAI CO. SO  
DATE 5-28-17 TIME 1753  
INTOXILYZER SERIAL # 80-001455 LOCATION BAGDAD

- (✓) 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution .100 AC.  
OR  
Pour a standard alcohol concentration solution \_\_\_\_\_ AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach 34°C ± 0.2°C
- (✓) 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
- (✓) 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
- (✓) 4. Air blank completed.
- (✓) 5. Calibration check completed. Test results 0.002 AC.
- (✓) 6. Air blank completed.
- (✓) 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE 

## INTOXILYZER 8000

Location: YAVAPAI SO, BAGDAD

Serial Number: 80-001455

Core Version: 8105.48

05/28/2017

17:53:26

-----  
Standard Lot#: 15615100A1Last Changed By: R. RIVERA #6356  
-----

QAS: R. RIVERA #6356

YAVAPAI CSO  
-----

Operator: R. RIVERA #6356

YAVAPAI CSO  
-----

Test	g./210L	Time
Air Blank	0.000	17:54:00
0.100 Cal Check	0.097	17:54:21
Air Blank	0.000	17:54:50

Preventative Maintenance  
Performed -- Timer Reset

EXHIBIT G-4  
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY

STANDARD QUALITY ASSURANCE PROCEDURES  
INTOXILYZER MODEL 8000

STANDARD QUALITY ASSURANCE PROCEDURE

QA SPECIALIST R. ROJERA AGENCY YAVAPAI CO. SO  
DATE 5-28-17 TIME 1759  
INTOXILYZER SERIAL # 80-001455 LOCATION BAGDAD

() 1. Display Reads "PUSH BUTTON TO START"

DIAGNOSTIC TESTS

- () 1. Clock time check.
- () 2. Date check.

OPERATIONAL TESTS

- () 1. Alcohol-free subject test result 0. 000 AC.
- () 2. Error recognition logic system functioning.  
Not a Successfully Completed Test Sequence printed
- () 3. Proper sample recognition system.  
Not a Successfully Completed Test Sequence printed  
Deficient sample printed.
- () 4. Standard Calibration Check standard 0. 100 AC. Result 0. 094 AC.

Instrument is operating properly and accurately. YES  NO

COMMENTS Ambien Fail

\_\_\_\_\_  
\_\_\_\_\_  
\_\_\_\_\_

SIGNATURE [Signature]

INTOXILYZER 8000  
 Location: YAVAPAI SO, BAGDAD  
 Serial Number: 80-001455  
 Core Version: 8105.48  
 05/28/2017 17:59:33

Standard Lot#: 15615100A1  
 Last Changed By: R. RIVERA #6356

QAS: R. RIVERA #6356  
 YAVAPAI CSO

Operator: R. RIVERA #6356  
 YAVAPAI CSO

Subject: SOAP,TEST  
 DOB: 07/18/1980  
 Sex: M Weight: 123

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	0.000	18:00:44
Diagnostic Test	Pass	18:01:19
Air Blank	0.000	18:01:47
0.100 Cal Check	0.098	18:02:08
Air Blank	0.000	18:02:37
Subject Test	0.000	18:02:52
Air Blank	ABT*	18:02:59
Air Blank	IPS**	18:03:07

\*Sequence Aborted

\*\*Improper Sample

**Not a Successfully Completed Test Sequence**

INTOXILYZER 8000  
 Location: YAVAPAI SO, BAGDAD  
 Serial Number: 80-001455  
 Core Version: 8105.48  
 05/28/2017 18:03:32

Standard Lot#: 15615100A1  
 Last Changed By: R. RIVERA #6356

OAS: R. RIVERA #6356  
 YAVAPAI CSO

Operator: R. RIVERA #6356  
 YAVAPAI CSO

Subject: SQAP,TEST  
 DOB: 07/18/1980  
 Sex: M Weight: 123

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	0.000	18:04:39
Diagnostic Test	Pass	18:05:14
Air Blank	0.000	18:05:41
0.100 Cal Check	0.098	18:06:02
Air Blank	0.000	18:06:32
Subject Test	0.030*	18:09:39
Air Blank	0.000	18:10:14
Air Blank	0.000	18:10:42
Subject Test	XXX**	18:11:16
Air Blank	RFI***	18:11:23

\*Deficient Sample

\*\*Invalid Sample

Begin new deprivation period

\*\*\*RFI Detect

**Not a Successfully Completed Test Sequence**

EXHIBIT G-2  
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY  
STANDARD QUALITY ASSURANCE PROCEDURES  
INTOXILYZER MODEL 8000  
STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST R. RIVERA AGENCY YAVAPAI CO. SO  
DATE 5-28-17 TIME 1813  
INTOXILYZER SERIAL # 80-001455 LOCATION BAGDAD

- (✓) 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution .100 AC.  
OR  
Pour a standard alcohol concentration solution \_\_\_\_\_ AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach 34°C ± 0.2°C
- (✓) 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
- (✓) 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
- (✓) 4. Air blank completed.
- (✓) 5. Calibration check completed. Test results 0.094 AC.
- (✓) 6. Air blank completed.
- (✓) 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE [Signature]

INTOXILYZER 8000  
 Location: YAVAPAI SO, BAGDAD  
 Serial Number: 80-001455  
 Core Version: 8105.48  
 05/28/2017 18:13:27

-----  
 Standard Lot#: 15615100A1  
 Last Changed By: R. RIVERA #6356  
 -----

QAS: R. RIVERA #6356  
 YAVAPAI CSO  
 -----

Operator: R. RIVERA #6356  
 YAVAPAI CSO  
 -----

Test	g/210L	Time
Air Blank	ANG*	18:13:58
Air Blank	0.000	18:14:26

\*Ambient Fail

INTOXILYZER 8000

Location: YAVAPAI SO, BAGDAD

Serial Number: 80-001455

Core Version: 8105.48

05/28/2017

18:14:29

-----  
Standard Lot#: 15615100A1

Last Changed By: R. RIVERA #6356

-----  
QAS: R. RIVERA #6356

YAVAPAI CSO

-----  
Operator: R. RIVERA #6356

YAVAPAI CSO

-----

Test	g/210L	Time
Air Blank	0.000	18:15:00
0.100 Cal Check	0.094	18:15:21
Air Blank	0.000	18:15:50

-----

Preventative Maintenance

Performed -- Timer Reset