


EXHIBIT G-2
 THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)
 ARIZONA DEPARTMENT OF PUBLIC SAFETY
 STANDARD QUALITY ASSURANCE PROCEDURES
 INTOXILYZER MODEL 8000
 STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST R. RIVERA AGENCY YAVAPAI CO. SO
 DATE 7-30-17 TIME 2350
 INTOXILYZER SERIAL # 80-001455 LOCATION Bagdad

- 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution .100 AC.
 OR
 Pour a standard alcohol concentration solution _____ AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach 34°C ± 0.2°C
- 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
- 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
- 4. Air blank completed.
- 5. Calibration check completed. Test results 0. 098 AC.
- 6. Air blank completed.
- 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE 

INTOXILYZER 8000
Location: YAVAPAI SO, BAGDAD
Serial Number: 80-001455
Core Version: 8105.48
07/30/2017 23:50:06

Standard Lot#: 27016100A1
Last Changed By: R. RIVERA #6356

QAS: R. RIVERA #6356
YAVAPAI CSO

Operator: R. RIVERA #6356
YAVAPAI CSO

Test	g/210L	Time
Air Blank	0.000	23:50:49
0.100 Cal Check	0.098	23:51:10
Air Blank	0.000	23:51:40

Preventative Maintenance
Performed -- Timer Reset