

EXHIBIT G-2
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY
STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000
STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST R. RIVERA AGENCY YAVAPAI CO. SO
DATE 8-28-17 TIME 1955
INTOXILYZER SERIAL # 80-001455 LOCATION BAGDAD

- () 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution .100 AC.
OR
Pour a standard alcohol concentration solution _____ AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$
- () 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
- () 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
- () 4. Air blank completed.
- () 5. Calibration check completed. Test results 0.098 AC.
- () 6. Air blank completed.
- () 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE 

INTOXILYZER 8000
Location: YAUAPAI SO, BAGDAD
Serial Number: 80-001455
Core Version: 8105.48
08/28/2017

19:55:58

Standard Lot#: 27016100A1
Last Changed By: R. RIVERA #6356

QAS: R. RIVERA #6356
YAUAPAI CSO

Operator: R. RIVERA #6356
YAUAPAI CSO

Test	g/210L	Time
Air Blank	0.000	19:56:27
0.100 Cal Check	0.098	19:56:48
Air Blank	0.000	19:57:18

Preventative Maintenance
Performed -- Timer Reset