

EXHIBIT G-2
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104(A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY
STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000
STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST R. Rivera

AGENCY YAVAPAI CO SO

DATE 10-3-17

TIME 2321

INTOXILYZER SERIAL # 80-001455

LOCATION Bagdad

- 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution, 100 AC.
OR
Pour a standard alcohol concentration solution _____ AC into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach 34 C ± 0.2°C
- 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
- 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
- 4. Air blank completed.
- 5. Calibration check completed. Test results 0.098 AC
- 6. Air blank completed
- 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE 

17-035934

INTOX:LYZER 8000
Location: YAUAPAI SO, BAGDAD
Serial Number: 80-001455
Core Version: 8105.48
10/03/2017

23:21:17

Standard Lot#: 27016100A1
Last Changed By: R. RIVERA #6356

QAS: R. RIVERA #6356
YAUAPAI CSO

Operator: R. RIVERA #6356
YAUAPAI CSO

Test	g/210L	Time
Air Blank	0.000	23:21:48
0.100 Cal Check	0.098	23:22:09
Air Blank	0.000	23:22:39

Preventative Maintenance
Performed -- Timer Reset