

EXHIBIT G-2  
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY  
STANDARD QUALITY ASSURANCE PROCEDURES  
INTOXILYZER MODEL 8000  
STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST R. Rivera AGENCY YAVAPAI S.O  
DATE 12-11-18 TIME 2141  
INTOXILYZER SERIAL # 80-001455 LOCATION Bagdad

1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution .100 AC.  
OR  
Pour a standard alcohol concentration solution \_\_\_\_\_ AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach  $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$
2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
4. Air blank completed.
5. Calibration check completed. Test results 0. 098 AC.
6. Air blank completed.
7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE [Signature]

INTOXILYZER 8000

Location: YAVAPAI SO, BAGDAD

Serial Number: 80-001455

Core Version: 8105.48

12/11/2018

21:41:29

-----  
Standard Lot#: 10818100A3

Last Changed By: R. RIVERA #6356  
-----

QAS: R. RIVERA #6356

YAVAPAI CSO  
-----

Operator: R. RIVERA #6356

YAVAPAI CSO  
-----

Test	g/210L	Time
Air Blank	0.000	21:42:02
0.100 Cal Check	0.098	21:42:23
Air Blank	0.000	21:42:52

Preventative Maintenance  
Performed -- Timer Reset