

EXHIBIT G-2
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY
STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000
STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST D. Ashby #562 AGENCY YCSO
DATE 10-22-10 TIME 0816
INTOXILYZER SERIAL # 80-001456 LOCATION Yarnell Substation

- (✓) 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution 0.100 AC.
OR
Pour a standard alcohol concentration solution _____ AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$
- (✓) 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
(✓) 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
(✓) 4. Air blank completed.
(✓) 5. Calibration check completed. Test results 0. 097 AC.
(✓) 6. Air blank completed.
(✓) 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE Dennis P. Ashby #562 / GAS # 33400

INTOXILYZER 8000

Location: YAVAPAI SO YARNELL

Serial Number: 80-001456

Core Version: 8105.48

06/22/2010

08:16:37

Standard Lot#: 585220

Last Changed By: D. ASHBY #S62

QAS: D. ASHBY #S62

YAVAPAI CSO

Operator: D. ASHBY #S62

YAVAPAI CSO

Test g/210L Time

Air Blank 0.000 08:17:27
0.100 Cal Check 0.097 08:17:48
Air Blank 0.000 08:18:18

Preventative Maintenance
Performed -- Timer Reset

EXHIBIT G-4

THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY

STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000

STANDARD QUALITY ASSURANCE PROCEDURE

QA SPECIALIST D. Ashby #562 AGENCY MOSE

DATE 10-22-10 TIME 0819

INTOXILYZER SERIAL # 30-001456 LOCATION Yarnell Substation

1. Display Reads "PUSH BUTTON TO START"

DIAGNOSTIC TESTS

- 1. Clock time check.
- 2. Date check.

OPERATIONAL TESTS

- 1. Alcohol-free subject test result 0. ACC AC.
- 2. Error recognition logic system functioning.
Not a Successfully Completed Test Sequence printed
- 3. Proper sample recognition system.
Not a Successfully Completed Test Sequence printed
Deficient sample printed.
- 4. Standard Calibration Check standard 0. 100 AC. Result 0. 097 AC.

Instrument is operating properly and accurately. YES NO

COMMENTS

SIGNATURE Diana L. Ashby #562 / QAS #33400

INTOXILYZER 8000

Location: YAVAPAI SO YARNELL

Serial Number: 80-001456

Core Version: 8105.48

06/22/2010

08:19:14

Standard Lot#: 585220

Last Changed By: D. ASHBY #862

QAS: D. ASHBY #862

YAVAPAI CSO

Operator: D. ASHBY #862

YAVAPAI CSO

Subject: SQAP,TEST

DOR: 01/01/1999

Sex: M

Weight: 150

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	0.000	08:20:03
Diagnostic Test	Pass	08:20:37
Air Blank	0.000	08:21:05
0.100 Cal Check	0.097	08:21:35
Air Blank	0.000	08:21:55
Subject Test	0.000	08:22:17
Air Blank	0.000	08:22:40
Wait	ABT*	
Air Blank	0.000	08:23:11

*Sequence Aborted

Not a Successfully Completed Test Sequence

INTOXILYZER 8000

Location: YAVAPAI CO YARNELL

Serial Number: 80-001456

Core Version: 8105.48

06/22/2010

08:23:23

Standard Lot#: 585220

Last Changed By: D. ASHBY

1362

OAS: D. ASHBY

1362

YAVAPAI CSO

Operator: D. ASHBY

1362

YAVAPAI CSO

Subject: SQAP,TEST

DOB: 01/01/1999

Sex: M

Weight: 150

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	IPS*	08:23:57
Air Blank	0.000	08:24:25

*Improper Sample

Not a Successfully Completed Test Sequence

INTOXILYZER 8000

Location: YAVAPAI SO YARNELL

Serial Number: 80-001456

Core Version: 8105.4B

06/22/2010

08:24:36

Standard Lot#: 585220

Last Changed By: D. ASHBY

#862

QAS: D. ASHBY

#862

YAVAPAI CSO

Operator: D. ASHBY

#862

YAVAPAI CSO

Subject: SOAP,TEST

DOB: 01/01/1999

Sex: M

Weight: 150

15 Minute Deprivation Period? Yes

Test g/210L Time

Air Blank RFI* 08:25:19
Air Blank 0.000 08:25:19

*RFI Detect

Not a Successfully Completed Test Report

INFOILYZER 8000

Location: YAVAPAI SO YARHPID

Serial Number: 80-001456

Core Version: 8105.148

08/22/2010

08:25:50

Standard Lot#: 585220

Last Changed By: D. ASHBY

4562

QAS: D. ASHBY

8102

YAVAPAI CSO

Operator: D. ASHBY

1063

YAVAPAI CSO

Subject: SOAP,TEST

DOB: 01/01/1999

Sex: M

Weight: 150

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	0.000	08:26:05
Diagnostic Test	Pass	08:27:17
Air Blank	0.000	08:27:46
0.100 Cal Check	0.096	08:28:07
Air Blank	0.000	08:28:29
Subject Test	0.443*	08:31:03
Air Blank	0.000	08:32:22
Air Blank	0.000	08:32:56
Subject Test	XXX**	08:33:15
Air Blank	0.000	08:33:40

*Deficient Sample

**Invalid Sample

Begin new deprivation period

Not a Successfully Completed Test Log

EXHIBIT G-2
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY
STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000
STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST D. Ashby #5162 AGENCY YCSO
DATE 10-22-10 TIME 0834
INTOXILYZER SERIAL # 92-003451 LOCATION Yarnell Substation

- () 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution 100 AC.
OR
Pour a standard alcohol concentration solution _____ AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$
- (✓) 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
(✓) 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
(✓) 4. Air blank completed.
(✓) 5. Calibration check completed. Test results 0. 097 AC.
(✓) 6. Air blank completed.
(✓) 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE D. Ashby #5162 / JAS #2345

INPOXLYERR 8000

Location: YAVAPAI SO YARBRELL

Serial Number: 89-001976

Core Version: 8105.42

06/22/2010

08:34:18

Standard Lot# 575220

Last Changed By: D. ASHBY #862

OAS: D. ASHBY #862

YAVAPAI CSO

Operator: D. ASHBY #862

YAVAPAI CSO

Test g/2UGL Time

Air Blank 0.000 08:34:50
0.100 Cal Check 0.097 08:35:10
Air Blank 0.000 08:35:41

Preventative Maintenance
Performed -- Timer Reset