

EXHIBIT G-2
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY
STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000
STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST R. BURKE #6350 AGENCY AZ DPS
DATE 07/26/2013 TIME 1501
INTOXILYZER SERIAL # 80-101156 ^{YES} ~~NO~~ LOCATION YARNELL PATROL

- (14) 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution 1.60 AC.
OR
Pour a standard alcohol concentration solution N/A AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach 34°C ± 0.2°C
- (15) 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
- (16) 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
- (17) 4. Air blank completed.
- (18) 5. Calibration check completed. Test results 0.077 AC.
- (19) 6. Air blank completed.
- (20) 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE *[Signature]* # 6350

INTOXILYZER 8000

Location: YAVAPAI SO YARNELL

Serial Number: 80-001456

Core Version: 8105.48

07/26/2013

14:58:27

Standard Lot#: 13712100A2

Last Changed By: D. ASHBY #862

QAS: R. BYRNE III #6386

AZ DPS

Operator: R. BYRNE III #6386

AZ DPS

Test g/210L Time

Air Blank 0.000 14:59:00
0.100 Cal Check 0.096 14:59:00
Air Blank RFI* 14:59:00
Air Blank 0.000 14:59:00

*RFI Detect

AS

INTOXILYZER 8000

Location: YAVAPAI SO YARNELL

Serial Number: 80-001456

Core Version: 8105.48

07/26/2013

15:01:55

Standard Lot#: 13712100A2

Last Changed By: D. ASHBY #S62

QAS: R. BYRNE III #6386

AZ DPS

Operator: R. BYRNE III #6386

AZ DPS

Test g/210L Time

Air Blank 0.000 15:02:28
0.100 Cal Check 0.097 15:02:37
Air Blank 0.000 15:03:17

Preventative Maintenance
Performed -- Timer Reset

RECEIVED
AUG 05 2013
TC