

EXHIBIT G-2  
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY  
STANDARD QUALITY ASSURANCE PROCEDURES  
INTOXILYZER MODEL 8000

STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST D. Ashby #2271 AGENCY YCSO

DATE 1-15-14 TIME 1406

INTOXILYZER SERIAL # 80-001456 LOCATION Yarnell Substation

- (✓) 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution 100 AC.  
OR  
Pour a standard alcohol concentration solution \_\_\_\_\_ AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach  $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$
- (✓) 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
- (✓) 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
- (✓) 4. Air blank completed.
- (✓) 5. Calibration check completed. Test results 0. 097 AC.
- (✓) 6. Air blank completed.
- (✓) 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE

D. Ashby #2271 / QAS #33400

INTOXILYZER 8000

Location: YAVAPAI SO YARNELL

Serial Number: 30-001456

Core Version: 8105.48

01/15/2014

14:06:01

-----  
Standard Lot#: 13712100A2

Last Changed By: D. ASHBY #S62

-----  
QAS: D. ASHBY #2271

YAVAPAI CSO

-----  
Operator: D. ASHBY #2271

YAVAPAI CSO

-----  
Test                      g/210L                      Time  
-----  
Air Blank                      0.000                      14:07:01  
0.100 Cal Check                0.097                      14:07:01  
Air Blank                      0.000                      14:07:01

Preventative Maintenance  
Performed -- Timer Reset

EXHIBIT G-4  
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY  
STANDARD QUALITY ASSURANCE PROCEDURES  
INTOXILYZER MODEL 8000

STANDARD QUALITY ASSURANCE PROCEDURE

QA SPECIALIST D. Ashby #2271 AGENCY YCSO  
DATE 1.15.14 TIME 1408  
INTOXILYZER SERIAL # 80-001456 LOCATION Yarnell Substation

(✓) 1. Display Reads "PUSH BUTTON TO START"

DIAGNOSTIC TESTS

- (✓) 1. Clock time check.  
(✓) 2. Date check.

OPERATIONAL TESTS

- (✓) 1. Alcohol-free subject test result 0. 000 AC.  
(✓) 2. Error recognition logic system functioning.  
Not a Successfully Completed Test Sequence printed  
(✓) 3. Proper sample recognition system.  
Not a Successfully Completed Test Sequence printed  
Deficient sample printed.  
(✓) 4. Standard Calibration Check standard 0. 100 AC. Result 0. 096 AC.

Instrument is operating properly and accurately. YES  NO

COMMENTS

\*\* Range Exceeded - I used too much mouthwash  
\* RFI Detect - My cell phone set it off when I received  
a text message \* RFI Detect - DPS Patrol  
Car Drove By  
SIGNATURE D. Ashby #2271 / QAS #33400

INTOXILYZER 8000  
Location: YAVAPAI SO YARNELL  
Serial Number: 80-001456  
Core Version: 8105.48  
01/15/2014

14:08:53

---

Standard Lot#: 13712100A2  
Last Changed By: D. ASHBY #S62

---

QAS: D. ASHBY #2271  
YAVAPAI CSO

---

Operator: D. ASHBY #2271  
YAVAPAI CSO

---

Subject: SOAP,TEST  
DOB: 01/01/1999  
Sex: M Weight: 150

---

15 Minute Deprivation Period? Yes

---

Test	g/210L	Time
Air Blank	0.000	14:10:14
Diagnostic Test	Pass	14:10:49
Air Blank	0.000	14:11:16
0.100 Cal Check	0.097	14:11:37
Air Blank	0.000	14:12:07
Subject Test	0.000	14:13:04
Air Blank	0.000	14:13:32
Wait	ABT*	
Air Blank	0.000	14:14:15

\*Sequence Aborted

**Not a Successfully Completed Test Sequence**

INTOXILYZER 8000

Location: YAVAPAI SO YARNELL

Serial Number: 80-001456

Core Version: 8105.48

01/15/2014

14:14:19

---

Standard Lot#: 13712100A2

Last Changed By: D. ASHBY

#302

---

QAS: D. ASHBY

#2771

YAVAPAI CSO

---

Operator: D. ASHBY

#2771

YAVAPAI CSO

---

Subject: SQAP,TEST

DOB: 01/01/1999

Sex: M

Weight: 150

---

15 Minute Deprivation Period? Yes

---

Test	g/210L	Time
Air Blank	IPS*	14:14:55
Air Blank	0.000	14:15:17

\*Improper Sample

**Not a Successfully Completed Test Sequence**

INTOXILYZER 8000  
Location: YAVAPAI SO YARNELL  
Serial Number: 80-001456  
Core Version: 8105.48  
01/15/2014

14:16:27

---

Standard Lot#: 13712100A2  
Last Changed By: D. ASHBY #862

---

QAS: D. ASHBY #2271  
YAVAPAI CSO

---

Operator: D. ASHBY #2271  
YAVAPAI CSO

---

Subject: SQAP,TEST  
DOB: 01/01/1999  
Sex: M Weight: 150

---

15 Minute Deprivation Period? Yes

---

Test	g/210L	Time
Air Blank	RFI*	14:16:59
Air Blank	0.000	14:17:27

\*RFI Detect

**Not a Successfully Completed Test Sequence**

INTOXILYZER 8000  
Location: YAVAPAI SO YARNELL  
Serial Number: 80-001456  
Core Version: 8105.48  
01/15/2014

14:17:50

---

Standard Lot#: 13712100A2  
Last Changed By: D. ASHBY #862

---

QAS: D. ASHBY #2271  
YAVAPAI CSO

---

Operator: D. ASHBY #2271  
YAVAPAI CSO

---

Subject: SQAP,TEST  
DOB: 01/01/1999  
Sex: M Weight: 150

---

15 Minute Deprivation Period? Yes

---

Test	g/210L	Time
Air Blank	0.000	14:18:40
Diagnostic Test	Pass	14:19:15
Air Blank	0.000	14:19:42
0.100 Cal Check	0.097	14:20:05
Air Blank	0.000	14:20:12
Subject Test	0.215*	14:23:00
Air Blank	0.000	14:24:12
Air Blank	0.000	14:24:42
Subject Test	RGE**	14:25:21
Air Blank	0.000	14:26:01

\*Deficient Sample

\*\*Range Exceeded

Not a Successfully Completed Test Sequence

INTOXILYZER 8000  
Location: YAVAPAI SO YARNELL  
Serial Number: 80-001456  
Core Version: 8105.48  
01/15/2014

14:26:19

---

Standard Lot#: 13712100A2  
Last Changed By: D. ASHBY #862

---

QAS: D. ASHBY #2271  
YAVAPAI CSO

---

Operator: D. ASHBY #2271  
YAVAPAI CSO

---

Subject: SQAP,TEST  
DOB: 01/01/1999  
Sex: M Weight: 150

---

15 Minute Deprivation Period? Yes

---

Test	g/210L	Time
Air Blank	0.000	14:27:09
Diagnostic Test	Pass	14:27:44
Air Blank	RFI*	14:28:03
Air Blank	0.000	14:28:32

\*RFI Detect

**Not a Successfully Completed Test Sequence**

INTOXILYZER 8000

Location: YAVAPAI SO YARNELL

Serial Number: 80-001456

Core Version: 8105.48

01/15/2014

14:29:22

---

Standard Lot#: 13712100A2

Last Changed By: D. ASHBY

#562

---

QAS: D. ASHBY

#2271

YAVAPAI CSO

---

Operator: D. ASHBY

#2271

YAVAPAI CSO

---

Subject: SOAP,TEST

DOB: 01/01/1999

Sex: M

Weight: 150

---

15 Minute Deprivation Period? Yes

---

Test	g/210L	Time
Air Blank	0.000	14:30:07
Diagnostic Test	Pass	14:30:42
Air Blank	0.000	14:31:10
0.100 Cal Check	0.097	14:31:31
Air Blank	0.000	14:32:00
Subject Test	RFI*	14:33:00
Air Blank	0.000	14:33:48

\*RFI Detect

**Not a Successfully Completed Test Sequence**

INTOXILYZER 8000  
Location: YAVAPAI SO YARNELL  
Serial Number: 80-001456  
Core Version: 8105.48  
01/15/2014

14:35:30

---

Standard Lot#: 13712100A2  
Last Changed By: D. ASHBY #862

---

QAS: D. ASHBY #2271  
YAVAPAI CSO

---

Operator: D. ASHBY #2271  
YAVAPAI CSO

---

Subject: SQAP,TEST  
DOB: 01/01/1999  
Sex: M Weight: 150

---

15 Minute Deprivation Period? Yes

---

Test	g/210L	Time
Air Blank	0.000	14:36:38
Diagnostic Test	Pass	14:36:53
Air Blank	0.000	14:37:26
0.100 Cal Check	0.096	14:37:43
Air Blank	0.000	14:38:10
Subject Test	0.091*	14:41:18
Air Blank	0.000	14:41:50
Air Blank	0.000	14:42:17
Subject Test	XXX**	14:43:02
Air Blank	0.000	14:43:30

\*Deficient Sample

\*\*Invalid Sample

Begin new deprivation period

**Not a Successfully Completed Test Sequence**

EXHIBIT G-2  
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY  
STANDARD QUALITY ASSURANCE PROCEDURES  
INTOXILYZER MODEL 8000

STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST D. Ashby #2271 AGENCY YCSO  
DATE 1.15.14 TIME 1444  
INTOXILYZER SERIAL # 80.001456 LOCATION Yarnell Substation

- (✓) 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution .100 AC.  
OR  
Pour a standard alcohol concentration solution \_\_\_\_\_ AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach  $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$
- (✓) 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
- (✓) 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
- (✓) 4. Air blank completed.
- (✓) 5. Calibration check completed. Test results 0.096 AC.
- (✓) 6. Air blank completed.
- (✓) 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE D. Ashby #2271 / QAS#33400

INTOXILYZER 8000

Location: YAVAPAI SO YARNELL

Serial Number: 90-001456

Core Version: 3105.48

01/15/2014

14:44:02

-----  
Standard Lot#: 13712100A2

Last Changed By: D. ASHBY #362

-----  
QAS: D. ASHBY #2271

YAVAPAI CSO

-----  
Operator: D. ASHBY #2271

YAVAPAI CSO

-----  
Test                                   g/210L                                   Time  
-----  
Air Blank                               0.000                                   14:44:13  
0.100 Cal Check                        0.096                                   14:44:13  
Air Blank                               0.000                                   14:45:23

Preventative Maintenance  
Performed -- Timer Reset