

EXHIBIT G-2
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY
STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000

STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST Murray NSZ AGENCY YCSO
DATE 4-27-08 TIME 9 50 10
INTOXILYZER SERIAL # 80-01734 LOCATION Seligman

1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution 100 AC.
OR
Pour a standard alcohol concentration solution _____ AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$
2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
4. Air blank completed.
5. Calibration check completed. Test results 0.099 AC.
6. Air blank completed.
7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE J. Murray

TEST NUMBER: 0000

Location: ACID ANALYSIS

Serial Number: 00001750

DATE/TIME

08/10/00

Standard Lot#: 731104

Part analyzed by: J. M. W. [unclear]

DATE OF RECEIPT: 08/10/00

ANALYST: J. M. W.

Time	g (100)	Min
Air Blank	0.000	09:51:28
Cal Check	0.000	09:51:28
Air Blank	0.000	09:51:28

J. M. W. NSZ

Operator's Signature

EXHIBIT G-4
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY
STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000
STANDARD QUALITY ASSURANCE PROCEDURE

QA SPECIALIST Murray NSB AGENCY YCSO
DATE 4-27-08 TIME 9 57 51
INTOXILYZER SERIAL # 80-01734 LOCATION Seligman

1. Display Reads "PUSH BUTTON TO START"

DIAGNOSTIC TESTS

- 1. Clock time check.
- 2. Date check.

OPERATIONAL TESTS

- 1. Alcohol-free subject test result 0. 000 AC.
- 2. Error recognition logic system functioning.
Not a Successfully Completed Test Sequence printed
- 3. Proper sample recognition system.
Not a Successfully Completed Test Sequence printed
Deficient sample printed.
- 4. Standard Calibration Check standard 0. 100 AC. Result 0. 099 AC.

Instrument is operating properly and accurately. YES NO

COMMENTS

Gas Tank:
735401G exp 12/09

SIGNATURE J. Murray NSB

INTOXILYZER 8000

Location: YCSO SELIGMAN

Serial Number : 80-001734

04/27/2008

09:57:51

Standard Lot#: 735401G

Last Changed By: J. MURRAY

#N-56

OAS: J. MURRAY

#N-56

YAVAPAI CSO

Operator: J. MURRAY

#N-56

YAVAPAI CSO

Subject: SOAP,TEST

DOB: 10/20/1960

Sex: F

Weight: 115

15 Minute Deprivation Period? Yes

Test	g/210I.	Time
Air Blank	0.000	09:59:28
Diagnostic Test	Pass	09:59:56
Air Blank	0.000	10:00:23
0.100 Cal Check	0.101	10:00:42
Air Blank	0.000	10:01:11
Subject Test	0.000	10:01:51
Air Blank	IPS*	10:02:03
Air Blank	RFI**	10:02:11

*Improper Sample

**RFI Detect

Not a Successfully Completed Test Sequence

INTOXILYZER 8000

Location: YCSO SELIGMAN

Serial Number : 80-001734

04/27/2008

10:03:40

Standard Lot#: 735401G

Last Changed By: J. MURRAY

#N-56

QAS: J. MURRAY

#N-56

YAVAPAI CSO

Operator: J. MURRAY

#N-56

YAVAPAI CSO

Subject: AQAP,TEST

DOB: 10/20/1960

Sex: F

Weight: 115

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	0.000	10:04:44
Diagnostic Test	Pass	10:05:12
Air Blank	0.000	10:05:39
0.100 Cal Check	0.100	10:05:58
Air Blank	0.000	10:06:27
Subject Test	0.000	10:07:03
Air Blank	0.000	10:07:30
Five Minute Wait	ABT*	
Air Blank	IPS**	10:07:47

*Sequence Aborted

**Improper Sample

Not a Successfully Completed Test Sequence

INTOXILYZER 8000

Location: YCSO SELIGMAN

Serial Number : 80-001734

04/27/2008

10:08:26

Standard Lot#: 735401G

Last Changed By: J. MURRAY

#N-56

QAS: J. MURRAY

#N-56

YAVAPAI CSO

Operator: J. MURRAY

#N-56

YAVAPAI CSO

Subject: SQAP,TEST

DOB: 10/20/1960

Sex: F

Weight: 115

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	0.000	10:09:22
Diagnostic Test	Pass	10:09:50
Air Blank	0.000	10:10:17
0.100 Cal Check	0.100	10:10:36
Air Blank	0.000	10:11:05
Subject Test	0.000*	10:14:14
Air Blank	0.000	10:14:46
Air Blank	0.000	10:15:13
Subject Test	XXX**	10:15:40
Air Blank	0.000	10:16:14

*Deficient Sample

**Invalid Sample

Begin new deprivation period

Not a Successfully Completed Test Sequence

EXHIBIT G-2
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY
STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000
STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST MURRAY NSL AGENCY YCSO
DATE 4-27-08 TIME 10 17 29
INTOXILYZER SERIAL 80-01734 LOCATION Seligman

- (X) 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution. 0.00 AC.
OR
Pour a standard alcohol concentration solution _____ AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$
- (X) 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
- (X) 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
- (+) 4. Air blank completed.
- (X) 5. Calibration check completed. Test results 0. 099 AC.
- (X) 6. Air blank completed.
- (X) 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE J. Murray

INTOXILYZER 8000

Serial Ion: Y630 41010004

Recd. Number: 100-60112

04/27/008

04/17/09


Standard Lot#: 2,54612

Last Checked by: J. MURRAY 10-56

QAS: J. MURRAY 10-56

YAVAPAI 030

Test	g/210L	Time
Air Blank	0.000	19:17:56
Cal. Check	0.000	19:18:15
Air Blank	0.000	19:18:41



Operator's Signature