

EXHIBIT G-2
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY
STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000
STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST Murray AGENCY YCSO
DATE 7-13-08 TIME 7:31:35
INTOXILYZER SERIAL # 80-01734 LOCATION Seligman

- (X) 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution 160 AC.
OR
Pour a standard alcohol concentration solution _____ AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$
- (X) 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
- (X) 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
- (X) 4. Air blank completed.
- (X) 5. Calibration check completed. Test results 0. 098 AC.
- (X) 6. Air blank completed.
- (X) 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE J. Murray

INTOXILYZER 8000

Location: YCSO SELIGMAN

Serial Number: 80-001734

07/13/2008

07:31:35

Standard Lot#: 735401G

Last Changed By: J. MURRAY #N-56

QAS: J. MURRAY #N-56

YAVAPAI CSO

Test	g/210L	Time
Air Blank	0.000	07:32:03
Cal Check	0.098	07:32:22
Air Blank	0.000	07:32:51

J. Murray

Operator's Signature

EXHIBIT G-4
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY
STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000

STANDARD QUALITY ASSURANCE PROCEDURE

QA SPECIALIST Murray AGENCY YCSO
DATE 7-13-08 TIME 7:45:44
INTOXILYZER SERIAL # 86-01734 LOCATION Seligman

1. Display Reads "PUSH BUTTON TO START"

DIAGNOSTIC TESTS

- 1. Clock time check.
- 2. Date check.

OPERATIONAL TESTS

- 1. Alcohol-free subject test result 0. 000 AC.
- 2. Error recognition logic system functioning.
Not a Successfully Completed Test Sequence printed
- 3. Proper sample recognition system.
Not a Successfully Completed Test Sequence printed
Deficient sample printed.
- 4. Standard Calibration Check standard 0. 100 AC. Result 0. 098 AC.

Instrument is operating properly and accurately. YES NO

COMMENTS

Bottle SN# 735401G

SIGNATURE

J. Murray NS6

INTOXILYZER 8000

Location: YCSO SELIGMAN

Serial Number : 80-001734

07/13/2008

07:45:44

Standard Lot#: 735401G

Last Changed By: J. MURRAY

#N-56

QAS: J. MURRAY

#N-56

YAVAPAI CSO

Operator: J. MURRAY

#N-56

YAVAPAI CSO

Subject: SQAP,TEST

DOB: 11/20/1960

Sex: F

Weight: 115

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	0.000	07:47:23
Diagnostic Test	Pass	07:47:50
Air Blank	0.000	07:48:18
0.100 Cal Check	0.098	07:48:36
Air Blank	0.000	07:49:06
Subject Test	0.000	07:49:23
Air Blank	IPS*	07:49:31
Air Blank	RFI**	07:49:43

*Improper Sample

**RFI Detect

Not a Successfully Completed Test Sequence

INTOXILYZER 8000

Location: YCSO SELIGMAN

Serial Number : 80-001734

07/13/2008

07:50:05

Standard Lot#: 735401G

Last Changed By: J. MURRAY

#N-56

QAS: J. MURRAY

#N-56

YAVAPAI CSO

Operator: J. MURRAY

#N-56

YAVAPAI CSO

Subject: SQAP,TEST

DOB: 10/20/1960

Sex: F

Weight: 115

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	ABT*	07:51:07
Air Blank	0.000	07:51:35

*Sequence Aborted

Not a Successfully Completed Test Sequence

INTOXILYZER 8000
Location: YCSO SELIGMAN
Serial Number : 80-001734
07/13/2008 07:51:45

Standard Lot#: 735401G
Last Changed By: J. MURRAY #N-56

QAS: J. MURRAY #N-56
YAVAPAI CSO

Operator: J. MURRAY #N-56
YAVAPAI CSO

Subject: SQAP,TEST
DOB: 10/20/1960
Sex: F Weight: 115

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	0.000	07:52:41
Diagnostic Test	Pass	07:53:08
Air Blank	0.000	07:53:36
0.100 Cal Check	0.099	07:53:54
Air Blank	0.000	07:54:24
Subject Test	0.168*	07:57:32
Air Blank	0.000	07:58:06
Air Blank	0.000	07:58:33
Subject Test	XXX**	07:59:00
Air Blank	0.000	07:59:36

*Deficient Sample
**Invalid Sample
Begin new deprivation period

Not a Successfully Completed Test Sequence

EXHIBIT G-2
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY
STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000

STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST Murray AGENCY YCSO
DATE 7-13-08 TIME 8:00:39
INTOXILYZER SERIAL # 80-01234 LOCATION Seligman

- (X) 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution 100 AC.
OR
Pour a standard alcohol concentration solution _____ AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$
- (X) 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
- (+) 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
- (+) 4. Air blank completed.
- (+) 5. Calibration check completed. Test results 0. 098 AC.
- (+) 6. Air blank completed.
- (+) 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE J. Murray NSB

INTOXILYZER 8000

Location: YCSO SELIGMAN

Serial Number: 80-001734

07/13/2008

08:00:39

Standard Lot#: 735401G

Last Changed By: J. MURRAY #N-56

QAS: J. MURRAY #N-56

YAVAPAI CSO

Test	g/210L	Time
Air Blank	0.000	08:01:07
Cal Check	0.098	08:01:26
Air Blank	0.000	08:01:55

J. Murray N56

Operator's Signature