

EXHIBIT G-2
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY
STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000

STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST J. Murray NSB AGENCY YCSO
DATE 10-13-08 TIME 1231 17
INTOXILYZER SERIAL # 80-01734 LOCATION Sellman

- (~~7~~) 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution .100 AC.
OR
Pour a standard alcohol concentration solution _____ AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$
- (~~X~~) 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
(~~X~~) 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
(~~X~~) 4. Air blank completed.
(~~X~~) 5. Calibration check completed. Test results 0.099 AC.
(~~7~~) 6. Air blank completed.
(~~X~~) 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE J. Murray NSB

INTOXILYZER 8000

Location: YCSO SELIGMAN

Serial Number: 80-001721

10/12/2008

12:31:17

Standard Lot#: 705401G

Last Changed By: J. MURRAY #N-56

QAS: J. MURRAY #N-56

YAVAPAI CSO

Test	g/210L	Time
Air Blank	0.000	12:31:44
Cal Check	0.099	12:32:03
Air Blank	0.000	12:32:33

J. Murray N56

Operator's Signature

EXHIBIT G-4
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY
STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000

STANDARD QUALITY ASSURANCE PROCEDURE

QA SPECIALIST Murray NSG AGENCY YCSO
DATE 10-13-08 TIME 1243 43
INTOXILYZER SERIAL # 80-01734 LOCATION Seligman

() 1. Display Reads "PUSH BUTTON TO START"

DIAGNOSTIC TESTS

- (~~X~~) 1. Clock time check.
(~~X~~) 2. Date check.

OPERATIONAL TESTS

- (~~+~~) 1. Alcohol-free subject test result 0. 000 AC.
(~~+~~) 2. Error recognition logic system functioning.
Not a Successfully Completed Test Sequence printed
(~~+~~) 3. Proper sample recognition system.
Not a Successfully Completed Test Sequence printed
Deficient sample printed.
(~~+~~) 4. Standard Calibration Check standard 0. 100 AC. Result 0. 098 AC.

Instrument is operating properly and accurately. YES X NO

COMMENTS

Dry Gas Serial # 7354016 12-28-09

SIGNATURE J. Murray NSG

INTOXILYZER 8000

Location: YCSO SELIGMAN

Serial Number : 80-001734

10/13/2008

12:48:33

Standard Lot#: 735401G

Last Changed By: J. MURRAY

#N-56

QAS: J. MURRAY

#N-56

YAVAPAI CSO

Operator: J. MURRAY

#N-56

YAVAPAI CSO

Subject: SQAP,TEST

DOB: 10/20/1960

Sex: F

Weight: 115

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	IPS*	12:49:36
Air Blank	0.000	12:50:04

*Improper Sample

Not a Successfully Completed Test Sequence

INTOXILYZER 8000

Location: YCSO SELIGMAN

Serial Number : 80-001734

10/13/2008

12:43:43

Standard Lot#: 735401G

Last Changed By: J. MURRAY

#N-56

OAS: J. MURRAY

#N-56

YAVAPAI CSO

Operator: J. MURRAY

#N-56

YAVAPAI CSO

Subject: SQAP,TEST

DOB: 10/20/1960

Sex: F

Weight: 115

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	0.000	12:45:47
Diagnostic Test	Pass	12:46:15
Air Blank	0.000	12:46:42
0.100 Cal Check	0.099	12:47:01
Air Blank	0.000	12:47:30
Subject Test	0.000	12:47:45
Air Blank	RFI*	12:47:59
Air Blank	0.000	12:48:27

*RFI Detect

Not a Successfully Completed Test Sequence

INTOXILYZER 8000
Location: YCSO SELIGMAN
Serial Number : 80-001734
10/13/2008 12:50:18

Standard Lot#: 735401G
Last Changed By: J. MURRAY #N-56

QAS: J. MURRAY #N-56
YAVAPAI CSO

Operator: J. MURRAY #N-56
YAVAPAI CSO

Subject: SQAP,TEST
DOB: 10/20/1960
Sex: F Weight: 115
15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	0.000	12:51:09
Diagnostic Test	Pass	12:51:37
Air Blank	0.000	12:52:04
0.100 Cal Check	0.099	12:52:23
Air Blank	0.000	12:52:52
Subject Test	0.000	12:53:10
Air Blank	0.000	12:53:37
Five Minute Wait	ABT*	
Air Blank	IPS**	12:54:00

*Sequence Aborted
**Improper Sample

Not a Successfully Completed Test Sequence

INTOXILYZER 8000

Location: YCSO SELIGMAN

Serial Number : 80-001734

10/13/2008

12:54:47

Standard Lot#: 735401G

Last Changed By: J. MURRAY

#N-56

QAS: J. MURRAY

#N-56

YAVAPAI CSO

Operator: J. MURRAY

#N-56

YAVAPAI CSO

Subject: SOAP,TEST

DOB: 10/20/1960

Sex: F

Weight: 115

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	0.000	12:55:43
Diagnostic Test	Pass	12:56:10
Air Blank	0.000	12:56:38
0.100 Cal Check	0.100	12:56:57
Air Blank	0.000	12:57:26
Subject Test	0.091*	13:00:34
Air Blank	0.000	13:01:09
Air Blank	0.000	13:01:37
Subject Test	XXX**	13:04:05
Air Blank	0.000	13:04:43

*Deficient Sample

**Invalid Sample

Begin new deprivation period

Not a Successfully Completed Test Sequence

EXHIBIT G-2
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY
STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000

STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST J. Murray AGENCY YCSO
DATE 10-13-08 TIME 130532
INTOXILYZER SERIAL # 80-01734 LOCATION Saligman

- (7) 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution .100 AC.
OR
Pour a standard alcohol concentration solution _____ AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$
- (+) 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
(+) 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
(+) 4. Air blank completed.
(+) 5. Calibration check completed. Test results 0. 098 AC.
(+) 6. Air blank completed.
(+) 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE J. Murray NSB

INTOXILYZER 8000

Location: YCSO BELMONT

Serial Number: 80-001734

10/13/2005

13:05:32

Standard Lot#: 7354010

Last Changed By: J. MURRAY #N-56

QAS: J. MURRAY #N-56

YAVAPAI CSO

Test	g/210L	Time
Air Blank	0.000	13:05:59
Cal Check	0.038	13:06:18
Air Blank	0.000	13:06:48

J. MURRAY #N-56

Operator's Signature