

EXHIBIT G-2
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY
STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000
STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST MURPHY AGENCY YCSO
DATE 2/17/09 TIME 12:19
INTOXILYZER SERIAL # 80-01734 LOCATION Seligman

- (X) 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution 100 AC.
OR
Pour a standard alcohol concentration solution _____ AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$
- (X) 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
- (X) 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
- (X) 4. Air blank completed.
- (X) 5. Calibration check completed. Test results 0. 099 AC.
- (X) 6. Air blank completed.
- (X) 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE J. MURPHY

INTOXILYZER 8000

Location: YCSO SELLINGEN

Serial Number: 80-001734

02/17/2009

12:19:20

Standard Lot#: 754010

Last Changed By: J. MURRAY #8597

OAC: J. MURRAY #8597

YAVAPAI CNO

Test	g/210L	Time
Air Blank	0.000	12:19:48
Cal Check	0.000	12:20:07
Air Blank	0.000	12:20:36

J. Murray

Operator's Signature

EXHIBIT G-4
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY
STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000

STANDARD QUALITY ASSURANCE PROCEDURE

QA SPECIALIST M. V. ASH AGENCY YESO
DATE 2-17-09 TIME 1229
INTOXILYZER SERIAL # 80-01734 LOCATION Seligman

1. Display Reads "PUSH BUTTON TO START"

DIAGNOSTIC TESTS

1. Clock time check.
 2. Date check.

OPERATIONAL TESTS

1. Alcohol-free subject test result 0. 050 AC.
 2. Error recognition logic system functioning.
Not a Successfully Completed Test Sequence printed
 3. Proper sample recognition system.
Not a Successfully Completed Test Sequence printed
Deficient sample printed.
 4. Standard Calibration Check standard 0. 100 AC. Result 0. 097 AC.

Instrument is operating properly and accurately. YES NO

COMMENTS 735401G 12-28-09 800+

SIGNATURE J. Murphy

INTOXILYZER 8000
Location: YCSO SELIGMAN
Serial Number : 80-001734
02/17/2009 12:29:29

Standard Lot#: 735401G
Last Changed By: J. MURRAY #N-56

QAS: J. MURRAY #5597
YAVAPAI CSO

Operator: J. MURRAY #5597
YAVAPAI CSO

Subject: SQAP,TEST
DOB: 10/20/1960
Sex: F Weight: 115

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	0.000	12:30:53
Diagnostic Test	Pass	12:31:20
Air Blank	0.000	12:31:47
0.100 Cal Check	0.098	12:32:06
Air Blank	0.000	12:32:36
Subject Test	0.000	12:33:00
Air Blank	RFI*	12:33:14
Air Blank	IPS**	12:33:41

*RFI Detect
**Improper Sample

Not a Successfully Completed Test Sequence

INTOXILYZER 8000

Location: YCSO SELIGMAN

Serial Number : 80-001734

02/17/2009

12:36:04

Standard Lot#: 735401G

Last Changed By: J. MURRAY

#N-56

QAS: J. MURRAY

#5597

YAVAPAI CSO

Operator: J. MURRAY

#5597

YAVAPAI CSO

Subject: SQAP, TEST

DOB: 10/20/1960

Sex: F

Weight: 115

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	0.000	12:37:03
Diagnostic Test	Pass	12:37:30
Air Blank	0.000	12:37:58
0.100 Cal Check	0.099	12:38:16
Air Blank	0.000	12:38:46
Subject Test	0.000	12:39:06
Air Blank	0.000	12:39:33
Five Minute Wait	ABT*	
Air Blank	0.000	12:40:24

*Sequence Aborted

Not a Successfully Completed Test Sequence

INTOXILYZER 8000

Location: YCSO SELIGMAN

Serial Number : 80-001734

02/17/2009

12:41:33

Standard Lot#: 735401G

Last Changed By: J. MURRAY

#N-56

OAS: J. MURRAY

#5597

YAVAPAI CSO

Operator: J. MURRAY

#5597

YAVAPAI CSO

Subject: SQAP,TEST

DOB: 10/20/1960

Sex: F

Weight: 115

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	0.000	12:43:40
Diagnostic Test	Pass	12:44:08
Air Blank	0.000	12:44:35
0.100 Cal Check	0.099	12:44:54
Air Blank	0.000	12:45:23
Subject Test	0.027*	12:48:32
Air Blank	0.000	12:49:06
Air Blank	0.000	12:49:33
Subject Test	XXX**	12:50:06
Air Blank	0.000	12:50:42

*Deficient Sample

**Invalid Sample

Begin new deprivation period

Not a Successfully Completed Test Sequence

EXHIBIT G-2
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY
STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000
STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST J. Murray AGENCY YCSO

DATE 2-17-09 TIME 1252

INTOXILYZER SERIAL # 80-01734 LOCATION Seligman

- () 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution _____ AC.
OR
Pour a standard alcohol concentration solution _____ AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$
- () 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
- () 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
- () 4. Air blank completed.
- () 5. Calibration check completed. Test results 0. 097 AC.
- () 6. Air blank completed.
- () 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE J. Murray

INTOXILYZER 6600

Location: YCSO SELLIGMAN

Serial Number: 00-001334

02/17/2009

12:52:21

Standard Lot#: 058010

Last Changed By: J. MURRAY #856

QAS: J. MURRAY #5597

YAVAPAI CSO

Test	g/210L	Time
Air Blank	0.000	12:52:58
Cal Check	0.097	12:53:17
Air Blank	0.000	12:53:46

J. Murray

Operator's Signature