

COPY

EXHIBIT G-2  
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY

STANDARD QUALITY ASSURANCE PROCEDURES  
INTOXILYZER MODEL 8000

STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST Murray NSB AGENCY YCSO  
DATE 12/24/09 TIME 10 03 01  
INTOXILYZER SERIAL # 80-01734 LOCATION Seligman

- (X) 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution .100 AC.  
OR  
Pour a standard alcohol concentration solution        AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach  $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$
- (X) 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
- (X) 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
- (X) 4. Air blank completed.
- (X) 5. Calibration check completed. Test results 0. 099 AC.
- (X) 6. Air blank completed.
- (X) 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE Jim Murray

INTOXILYZER 8000

Location: YCSO SELIGMAN

Serial Number: 80-001734

12/26/2009

10:03:01

-----  
Standard Lot#: 9014011

Last Changed By: J. MURRAY #5597

-----  
QAS: J. MURRAY #5597

YAVAPAI CSO  
-----

Test	g/210L	Time
Air Blank	0.000	10:03:29
Cal Check	0.099	10:03:47
Air Blank	0.000	10:04:17

*J. Murray*

\_\_\_\_\_  
Operator's Signature

EXHIBIT G-4  
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY

STANDARD QUALITY ASSURANCE PROCEDURES  
INTOXILYZER MODEL 8000

STANDARD QUALITY ASSURANCE PROCEDURE

QA SPECIALIST J. Murray AGENCY YESO  
DATE 12/26/09 TIME 10 06 28  
INTOXILYZER SERIAL # 80-01234 LOCATION Seligman

1. Display Reads "PUSH BUTTON TO START"

DIAGNOSTIC TESTS

1. Clock time check.  
 2. Date check.

OPERATIONAL TESTS

1. Alcohol-free subject test result 0. 000 AC.  
 2. Error recognition logic system functioning.  
Not a Successfully Completed Test Sequence printed  
 3. Proper sample recognition system.  
Not a Successfully Completed Test Sequence printed  
Deficient sample printed.  
 4. Standard Calibration Check standard 0. 100 AC. Result 0. 099 AC.

Instrument is operating properly and accurately. YES  NO

COMMENTS Bottle # 9014011 500PSI +-

NAME J. Murray

INTOXILYZER 8000  
Location: YCSO SELIGMAN  
Serial Number : 80-001734  
12/26/2009 10:12:25

Standard Lot#: 9014011  
Last Changed By: J. MURRAY #5597

QAS: J. MURRAY #5597  
YAVAPAI CSO

Operator: J. MURRAY #5597  
YAVAPAI CSO

Subject: TEST,SQAP  
DOB: 10/20/1960  
Sex: F Weight: 115

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	0.000	10:13:36
Diagnostic Test	Pass	10:14:03
Air Blank	0.000	10:14:30
0.100 Cal Check	0.099	10:14:49
Air Blank	0.000	10:15:18
Subject Test	0.074*	10:18:27
Air Blank	0.000	10:19:01
Air Blank	0.000	10:19:28
Subject Test	XXX**	10:20:01
Air Blank	RFI***	10:20:14

\*Deficient Sample

\*\*Invalid Sample

Begin new deprivation period

\*\*\*RFI Detect

**Not a Successfully Completed Test Sequence**

INTOXILYZER 8000  
Location: YCSO SELIGMAN  
Serial Number : 80-001734  
12/26/2009 10:06:28

---

Standard Lot#: 9014011  
Last Changed By: J. MURRAY #5597

---

QAS: J. MURRAY #5597  
YAVAPAI CSO

---

Operator: J. MURRAY #5597  
YAVAPAI CSO

---

Subject: SQAP,TEST  
DOB: 10/20/1960  
Sex: F Weight: 115

---

15 Minute Deprivation Period? Yes

---

Test	g/210L	Time
Air Blank	0.000	10:07:44
Diagnostic Test	Pass	10:08:11
Air Blank	0.000	10:08:38
0.100 Cal Check	0.099	10:08:57
Air Blank	0.000	10:09:26
Subject Test	0.000	10:09:50
Air Blank	0.000	10:10:18
Five Minute Wait	ABT*	
Air Blank	IPS**	10:11:06

\*Sequence Aborted

\*\*Improper Sample

**Not a Successfully Completed Test Sequence**

INTOXILYZER 8000  
Location: YCSO SELIGMAN  
Serial Number : 80-001734  
12/26/2009 10:23:38

---

Standard Lot#: 9014011  
Last Changed By: J. MURRAY #5597

---

QAS: J. MURRAY #5597  
YAVAPAI CSO

---

Operator: J. MURRAY #5597  
YAVAPAI CSO

---

Subject: TEST,SQAP  
DOB: 10/20/1960  
Sex: F Weight: 115

---

15 Minute Deprivation Period? Yes

---

Test	g/210L	Time
Air Blank	0.000	10:24:33
Diagnostic Test	Pass	10:25:00
Air Blank	0.000	10:25:27
0.100 Cal Check	0.100	10:25:46
Air Blank	0.000	10:26:16
Subject Test	RGE*	10:26:39
Air Blank	0.000	10:27:18

\*Range Exceeded

**Not a Successfully Completed Test Sequence**

EXHIBIT G-2  
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY  
STANDARD QUALITY ASSURANCE PROCEDURES  
INTOXILYZER MODEL 8000  
STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST J. Murray AGENCY Y250  
DATE 12/26/09 TIME 10 28 23  
INTOXILYZER SERIAL # 80-01734 LOCATION Sally MIA

- (X) 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution .100 AC.  
OR  
Pour a standard alcohol concentration solution ✓ AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach  $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$
- (X) 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
- (X) 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
- (X) 4. Air blank completed.
- (X) 5. Calibration check completed. Test results 0. 099 AC.
- (X) 6. Air blank completed.
- (X) 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE J. Murray


INTOXILYZER 8000  
Location: YCSO SELIGMAN  
Serial Number: 80-001734  
12/26/2009

10:28:23

-----  
Standard Lot#: 9014011  
Last Changed By: J. MURRAY #5597  
-----

QAS: J. MURRAY #5597  
YAVAPAI CSO  
-----

Test	g/210L	Time
Air Blank	0.000	10:28:51
Cal Check	0.099	10:29:09
Air Blank	0.000	10:29:39



\_\_\_\_\_  
Operator's Signature