

EXHIBIT G-2

THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY

STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000

STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST MURPHY AGENCY YCSO
DATE 7-24-10 TIME 8 18 46
INTOXILYZER SERIAL # 80-01734 LOCATION Selma

1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution 100 AC.
OR
Pour a standard alcohol concentration solution AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$
2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
4. Air blank completed.
5. Calibration check completed. Test results 0.098 AC.
6. Air blank completed.
7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE J. Murphy NSB

INTOXILYZER 8000

Location: YCSO SELIGMAN

Serial Number: 80-001734

Core Version: 8105.48

07/24/2010

08:18:46

Standard Lot#: 901401I

Last Changed By: J. MURRAY #5597

QAS: J. MURRAY #N56

YAVAPAI YCSO

Operator: J. MURRAY #N56

YAVAPAI YCSO

Test	g/210L	Time
Air Blank	0.000	08:19:26
0.100 Cal Check	0.098	08:19:47
Air Blank	0.000	08:20:16

Preventative Maintenance
Performed -- Timer Reset

EXHIBIT G-4
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY
STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000

STANDARD QUALITY ASSURANCE PROCEDURE

QA SPECIALIST MURPHY AGENCY YCSO
DATE 7-24-10 TIME 8 21 57
INTOXILYZER SERIAL # 80-01734 LOCATION Seligman

1. Display Reads "PUSH BUTTON TO START"

DIAGNOSTIC TESTS

1. Clock time check.
 2. Date check.

OPERATIONAL TESTS

1. Alcohol-free subject test result 0. 000 AC.
 2. Error recognition logic system functioning.
Not a Successfully Completed Test Sequence printed
 3. Proper sample recognition system.
Not a Successfully Completed Test Sequence printed
Deficient sample printed.
 4. Standard Calibration Check standard 0. 100 AC. Result 0. 096 AC.

Instrument is operating properly and accurately. YES NO

COMMENTS Bottle 901401I +-300PST

SIGNATURE J. Murray NSB

INTOXILYZER 8000
Location: YCSO SELIGMAN
Serial Number: 80-001734
Core Version: 8105.48
07/24/2010

08:21:57

Standard Lot#: 901401I
Last Changed By: J. MURRAY #5597

QAS: J. MURRAY #N56
YAVAPAI YCSO

Operator: J. MURRAY #N56
YAVAPAI YCSO

Subject: SQAP,TEST
DOB: 10/20/1960
Sex: F Weight: 115

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	0.000	08:23:01
Diagnostic Test	Pass	08:23:36
Air Blank	0.000	08:24:04
0.100 Cal Check	0.098	08:24:24
Air Blank	0.000	08:24:54
Subject Test	DEF*	08:28:01
Air Blank	0.000	08:28:29
Air Blank	0.000	08:28:56
Subject Test	0.000	08:29:32
Air Blank	0.000	08:30:00
Wait	ABT**	
Air Blank	IPS***	08:31:06

*Deficient Sample
**Sequence Aborted
***Improper Sample

Not a Successfully Completed Test Sequence

INTOXILYZER 8000
Location: YCSO SELIGMAN
Serial Number: 80-001734
Core Version: 8105.48
07/24/2010

08:31:19

Standard Lot#: 901401I
Last Changed By: J. MURRAY #5597

QAS: J. MURRAY #N56
YAVAPAI YCSO

Operator: J. MURRAY #N56
YAVAPAI YCSO

Subject: SQAP, TWST
DOB: 10/20/1960
Sex: F Weight: 115

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	0.000	08:32:20
Diagnostic Test	Pass	08:32:55
Air Blank	0.000	08:33:22
0.100 Cal Check	0.096	08:33:43
Air Blank	0.000	08:34:13
Subject Test	0.337*	08:37:20
Air Blank	0.000	08:37:55
Air Blank	0.000	08:38:22
Subject Test	XXX**	08:39:19
Air Blank	RFI***	08:39:27

*Deficient Sample

**Invalid Sample

Begin new deprivation period

***RFI Detect

Not a Successfully Completed Test Sequence

EXHIBIT G-2
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY

STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000

STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST Murray AGENCY YCSO
DATE 7-24-10 TIME 8:45:19
INTOXILYZER SERIAL # 80-01734 LOCATION Seligman

- (x) 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution 1100 AC.
OR
Pour a standard alcohol concentration solution _____ AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$
- (x) 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
- (x) 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
- (x) 4. Air blank completed.
- (x) 5. Calibration check completed. Test results 0.096 AC.
- (x) 6. Air blank completed.
- (x) 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE J. Murray NSG

INTOXILYZER 8000

Location: YCSO SELIGMAN

Serial Number: 80-001734

Core Version: 8105.48

07/24/2010

08:44:11

Standard Lot#: 901401I

Last Changed By: J. MURRAY #5597

QAS: J. MURRAY #N56

YAVAPAI YCSO

Operator: J. MURRAY #N56

YAVAPAI YCSO

Test	g/210L	Time
Air Blank	AMB*	08:44:48
Air Blank	0.000	08:45:16

*Ambient Fail

INTOXILYZER 8000

Location: YCSO SELIGMAN

Serial Number: 80-001734

Core Version: 8105.48

07/24/2010

08:45:19

Standard Lot#: 901401I

Last Changed By: J. MURRAY #5597

QAS: J. MURRAY #N56

YAVAPAI YCSO

Operator: J. MURRAY #N56

YAVAPAI YCSO

Test	g/210L	Time
Air Blank	0.000	08:45:49
0.100 Cal Check	0.096	08:46:10
Air Blank	0.000	08:46:40

Preventative Maintenance
Performed -- Timer Reset