

COPY

EXHIBIT G-2
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY
STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000

STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST J. MURPHY NSB AGENCY Yavapai County Sheriff's Office
DATE 5/4/11 TIME 131802
INTOXILYZER SERIAL # 80-01734 LOCATION Seligman 103

- (~~X~~) 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution 100 AC.
OR
Pour a standard alcohol concentration solution AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$
- (~~X~~) 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
- (~~X~~) 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
- (~~X~~) 4. Air blank completed.
- (~~X~~) 5. Calibration check completed. Test results 0. 099 AC.
- (~~X~~) 6. Air blank completed.
- (~~X~~) 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE J. Murphy NSB

INTOXILYZER 8000

Location: YCSO SELIGMAN

Serial Number: 80-001734

Core Version: 8105.48

05/04/2011

13:18:02

Standard Lot#: 779386

Last Changed By: J. MURRAY #N56

QAS: J. MURRAY #N56

YAVAPAI YCSO

Operator: J. MURRAY #N56

YAVAPAI YCSO

Test g/210L Time

Air Blank 0.000 13:18:39
0.100 Cal Check 0.099 13:19:00
Air Blank 0.000 13:19:30

Preventative Maintenance
Performed -- Timer Reset

EXHIBIT G-4
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY
STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000

STANDARD QUALITY ASSURANCE PROCEDURE

QA SPECIALIST J. Murray NSB AGENCY Yavapai County Sheriff's Office
DATE 5/4/2011 TIME 13:21:50
INTOXILYZER SERIAL # 80-01734 LOCATION Seligman 103

1. Display Reads "PUSH BUTTON TO START"

DIAGNOSTIC TESTS

1. Clock time check.
 2. Date check.

OPERATIONAL TESTS

1. Alcohol-free subject test result 0. 000 AC.
 2. Error recognition logic system functioning.
Not a Successfully Completed Test Sequence printed
 3. Proper sample recognition system.
Not a Successfully Completed Test Sequence printed
Deficient sample printed.
 4. Standard Calibration Check standard 0. 100 AC. Result 0. 097 AC.

Instrument is operating properly and accurately. YES NO

COMMENTS

779386 +900 PSI

SIGNATURE J. Murray NSB

INTOXILYZER 8000
Location: YCSO SELIGMAN
Serial Number: 80-001734
Core Version: 8105.48
05/04/2011

13:21:50

Standard Lot#: 779386
Last Changed By: J. MURRAY #N56

QAS: J. MURRAY #N56
YAVAPAI YCSO

Operator: J. MURRAY #N56
YAVAPAI YCSO

Subject: SQAP,TEST
DOB: 10/20/1960
Sex: F Weight: 115

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	0.000	13:22:50
Diagnostic Test	Pass	13:23:25
Air Blank	0.000	13:23:52
0.100 Cal Check	0.099	13:24:13
Air Blank	0.000	13:24:43
Subject Test	0.000	13:25:26
Air Blank	0.000	13:25:53
Wait	ABT*	
Air Blank	IPS**	13:26:32

*Sequence Aborted

**Improper Sample

Not a Successfully Completed Test Sequence

INTOXILYZER 8000
Location: YCSO SELIGMAN
Serial Number: 80-001734
Core Version: 8105.48
05/04/2011

13:26:40

Standard Lot#: 779386
Last Changed By: J. MURRAY #N56

QAS: J. MURRAY #N56
YAVAPAI YCSO

Operator: J. MURRAY #N56
YAVAPAI YCSO

Subject: SOAP,TEST
DOB: 10/20/1960
Sex: F Weight: 115

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	RFI*	13:27:25
Air Blank	IPS**	13:27:34

*RFI Detect
**Improper Sample

Not a Successfully Completed Test Sequence

INTOXILYZER 8000
Location: YCSO SELIGMAN
Serial Number: 80-001734
Core Version: 8105.48
05/04/2011 13:28:26

Standard Lot#: 779386
Last Changed By: J. MURRAY #N56

QAS: J. MURRAY #N56
YAVAPAI YCSO

Operator: J. MURRAY #N56
YAVAPAI YCSO

Subject: SQAP,EST
DOB: 10/20/1960
Sex: F Weight: 115

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	0.000	13:29:27
Diagnostic Test	Pass	13:30:02
Air Blank	0.000	13:30:29
0.100 Cal Check	0.098	13:30:50
Air Blank	0.000	13:31:19
Subject Test	0.130*	13:34:27
Air Blank	0.000	13:35:00
Air Blank	0.000	13:35:27
Subject Test	XXX**	13:35:57
Air Blank	0.000	13:36:34

*Deficient Sample
**Invalid Sample
Begin new deprivation period

Not a Successfully Completed Test Sequence

EXHIBIT G-2
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY
STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000
STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST J. MURRAY NSL AGENCY Yavapai County Sheriff's Office
DATE 5/4/2011 TIME 13:36:48
INTOXILYZER SERIAL # 80-01734 LOCATION Seligman 103

- (x) 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution 100 AC.
OR
Pour a standard alcohol concentration solution AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$
- (x) 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
- (x) 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
- (x) 4. Air blank completed.
- (x) 5. Calibration check completed. Test results 0.097 AC.
- (x) 6. Air blank completed.
- (x) 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE J. MURRAY NSL

DPS Form Exh G-2 (Rev 05-01)

SPECIALIST
EXHIBIT
1. Ensure the
2. Intoxilyzer
3. Calibration
4. Air blank
5. Air blank
6. Remove

SIGNATURE

INTOXILYZER 8000

Location: YCSO SELIGMAN

Serial Number: 80-001734

Core Version: 8105.48

05/04/2011

13:36:48

Standard Lot#: 779386

Last Changed By: J. MURRAY #N56

QAS: J. MURRAY #N56

YAVAPAI YCSO

Operator: J. MURRAY #N56

YAVAPAI YCSO

Test	g/210L	Time
Air Blank	0.000	13:37:31
0.100 Cal Check	0.097	13:37:51
Air Blank	0.000	13:38:21

Preventative Maintenance
Performed -- Timer Reset