

EXHIBIT G-2  
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY  
STANDARD QUALITY ASSURANCE PROCEDURES  
INTOXILYZER MODEL 8000

STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST J. MURPHY NSB/5597 AGENCY Yavapai County Sheriff's Office  
DATE 7/15/12 TIME 14:53:14  
INTOXILYZER SERIAL # 80-01734 LOCATION Selby Man 103

- (X) 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution .100 AC.  
OR  
Pour a standard alcohol concentration solution        AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach  $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$
- (X) 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
- (X) 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
- (X) 4. Air blank completed.
- (X) 5. Calibration check completed. Test results 0.098 AC.
- (X) 6. Air blank completed.
- (X) 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE J. Murphy NSB/5597

INTOXILYZER 8000

Location: YCSO SELIGMAN

Serial Number: 80-001734

Core Version: 8105.48

07/15/2012

14:53:14

-----  
Standard Lot#: 09312100A1

Last Changed By: J. MURRAY #N56  
-----

QAS: J. MURRAY #N56

YAVAPAI 5597  
-----

Operator: J. MURRAY #N56

YAVAPAI 5597  
-----

Test	g/210L	Time
Air Blank	0.000	14:53:51
0.100 Cal Check	0.098	14:54:12
Air Blank	0.000	14:54:41

Preventative Maintenance

Performed -- Timer Reset

EXHIBIT G-4  
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY  
STANDARD QUALITY ASSURANCE PROCEDURES  
INTOXILYZER MODEL 8000

STANDARD QUALITY ASSURANCE PROCEDURE

QA SPECIALIST J. MURRAY NSB/5597 AGENCY Yavapai County Sheriff's Office  
DATE 7/15/12 TIME 14:57:07  
INTOXILYZER SERIAL # 80-01734 LOCATION Sherman 103

1. Display Reads "PUSH BUTTON TO START"

DIAGNOSTIC TESTS

1. Clock time check.  
 2. Date check.

OPERATIONAL TESTS

1. Alcohol-free subject test result 0. 000 AC.  
 2. Error recognition logic system functioning.  
Not a Successfully Completed Test Sequence printed  
 3. Proper sample recognition system.  
Not a Successfully Completed Test Sequence printed  
Deficient sample printed.  
 4. Standard Calibration Check standard 0. 100 AC. Result 0. 097 AC.

Instrument is operating properly and accurately. YES  NO

COMMENTS

09312150A1 exp. 6/1/2014  
+1000 PSI

SIGNATURE J. MURRAY NSB/5597

INTOXILYZER 8000

Location: YCSO SELIGMAN

Serial Number: 80-001734

Core Version: 8105.48

07/15/2012

14:57:07

---

Standard Lot#: 09312100A1

Last Changed By: J. MURRAY

#N56

---

OAS: J. MURRAY

#N56

YAVAPAI 5597

---

Operator: J. MURRAY

#N56

YAVAPAI YCSO

---

Subject: TEST, SQAP

DOB: 10/20/1960

Sex: F

Weight: 115

---

15 Minute Deprivation Period? Yes

---

Test	g/210L	Time
Air Blank	0.000	14:58:01
Diagnostic Test	Pass	14:58:39
Air Blank	0.000	14:59:06
0.100 Cal Check	0.098	14:59:27
Air Blank	0.000	14:59:56
Subject Test	0.000	15:00:19
Air Blank	0.000	15:00:46
Wait	ABT*	
Air Blank	RFI**	15:01:03

\*Sequence Aborted

\*\*RFI Detect

**Not a Successfully Completed Test Sequence**

INTOXILYZER 8000

Location: YCSO SELIGMAN

Serial Number: 80-001734

Core Version: 8105.48

07/15/2012

15:01:20

---

Standard Lot#: 09312100A1

Last Changed By: J. MURRAY

#N56

---

QAS: J. MURRAY

#N56

YAVAPAI 5597

---

Operator: J. MURRAY

#N56

YAVAPAI YCSO

---

Subject: SQAP,TEST

DOB: 10/20/1960

Sex: F

Weight: 115

---

15 Minute Deprivation Period? Yes

---

Test	g/210L	Time
Air Blank	IPS*	15:02:03
Air Blank	RFI**	15:02:10

\*Improper Sample

\*\*RFI Detect

**Not a Successfully Completed Test Sequence**

INTOXILYZER 8000

Location: YCSO SELIGMAN

Serial Number: 80-001734

Core Version: 8105.48

07/15/2012

15:02:26

---

Standard Lot#: 09312100A1

Last Changed By: J. MURRAY

#N56

---

OAS: J. MURRAY

#N56

YAVAPAI 5597

---

Operator: J. MURRAY

#N56

YAVAPAI YCSO

---

Subject: SQAP, TEST

DOB: 10/20/1960

Sex: F

Weight: 115

---

15 Minute Deprivation Period? Yes

---

Test	g/210L	Time
Air Blank	0.000	15:03:23
Diagnostic Test	Pass	15:03:58
Air Blank	0.000	15:04:25
0.100 Cal Check	0.097	15:04:46
Air Blank	0.000	15:05:16
Subject Test	0.276*	15:08:23
Air Blank	0.000	15:08:57
Air Blank	0.000	15:09:24
Subject Test	XXX**	15:09:48
Air Blank	0.000	15:10:24

\*Deficient Sample

\*\*Invalid Sample

Begin new deprivation period

**Not a Successfully Completed Test Sequence**

EXHIBIT G-2  
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY

STANDARD QUALITY ASSURANCE PROCEDURES  
INTOXILYZER MODEL 8000

STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST J. Murray NSB/5597 AGENCY Yavapai County Sheriff's Office  
DATE 7/15/12 TIME 15:11:58  
INTOXILYZER SERIAL # 80-01734 LOCATION Selma 103

1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution (5) 60 AC.  
OR  
Pour a standard alcohol concentration solution          AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach  $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$
2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
4. Air blank completed.
5. Calibration check completed. Test results 0.097 AC.
6. Air blank completed.
7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE J. Murray NSB/5597

INTOXILYZER 8000

Location: YCSO SELIGMAN

Serial Number: 80-001734

Core Version: 8105.48

07/15/2012

15:11:58

-----  
Standard Lot#: 09312100A1

Last Changed By: J. MURRAY #N56  
-----

QAS: J. MURRAY #N56

YAVAPAI YCSO  
-----

Operator: J. MURRAY #N56

YAVAPAI YCSO  
-----

Test	g/210L	Time
Air Blank	0.000	15:12:34
0.100 Cal Check	0.097	15:12:54
Air Blank	0.000	15:13:24

Preventative Maintenance  
Performed -- Timer Reset